

INDEX

- ABRADABILITY INDEX 20
 ABS RESIN 28
 ACID BASE ACCOUNTING 2
 AEROSOLS 43
 AIR PARTICULATE 2
 ALUMINA REFRACTORIES 29
 ASH 9, 10, 13, 31
 ATTRITION INDEX 2
 AUTOMOBILE CATALYST 14
- BASIC SLAG 36
 BLAINE 5
 BLAST FURNACE SLAG 37
 BORON CARBIDE 42
 BORON NITRIDE 34
 BRINELL 22
 BURNT REFRACTORIES 30
- CALCIUM ALUMINATE 2
 CALCIUM CARBONATE 2
 CARBIDE 34, 42
 CASTING POWDER 2, 11
 CEMENT 3, 4, 5
 CEMENT SIEVING 5
 CENOSPHERE 6
 CERAMIC POWDER 32
 CHARPY IMPACT 22
 CHROMIUM CARBIDE 42
 CLINKER 5
 COAL 6, 7, 8, 9
 COAL ASH 9, 10
 COAL FLY ASH 10
 COAL WASTE ROCK 10
 COAL-TAR PITCH 11
 COKE 12, 13
 COKE ASH 13
 CONTINUOUS CASTING POWDER 11
 CONVERTER SLAG 38
 COVER POWDER 11
 Cr-Mg REFRACTOIREIS 30
- DRY ANALYSIS 9
 DUST 13, 14
- ELECTRIC FURNACE SLAG 38
 ELECTRODE CARBON 14
 ELECTROLYTIC MANGANESE 23
 ELECTRONIC SCRAP 15
 EXHAUST CATALYST 14
- FERROALLOY SLAG 36
 FERROBORON 15
 FERROCHROMIUM 15
 FERROMANGANESE 16
 FERROMOLYBDENUM 16
 FERRONICKEL 17
 FERRONIUBIUM 17
 FERROPHOSPHORUS 17
 FERROSILICOALUMINUM 18
 FERROSILICOCHROMIUM 18
 FERROSILICON 18, 19
 FERROSILICOTITANIUM 18
 FERROTITANIUM 17
 FERROTUNGSTEN 18
 FERROVANADIUM 18
 FILTER MEDIA 2
 FINENESS 5
 FIRECLAY REFRACTORIES 30
- FLUE DUST 14
 FLUORINE SLAG 38
 FLUX 20
 FLY ASH 10, 11
 FOUNDRY SAND 32
 FURNACE DUST 14
 FURNACE SLAG 38
 FUSIBILITY OF COAL 8
 FUSIBILITY OF COAL ASH 9
- GLASS 20, 21
 GLASS SAND 21
 GRAVEL 21
- HARDGROVE GRINDABILITY INDEX 22
 HARDNESS 22, 42
 HEXAVALENT CHROMIUM 20
- IMPACT 22
 INCINERATED WASTE 22
 INDOOR DUST 14
 INDUSTRIAL FLY ASH 11
 IRON MAKING SLAG 36, 40
 IRON OXIDE 25
 IZOD 22
- LAYER THICKNESS 23
 LEAD PAINT 23
 LEEB 22
- MAGNESIA REFRACTORIES 31
 MAGNETIC CENOSPHERE 6
 MANGANESE 23
 MANGANESE SLAG 38
 MASS VOLUME 26
 MELTING POINT 23
- NANOSCALE 23
 NICKEL OXIDE 25
 NITRIDE 34, 43
- OXIDE 24, 25, 26
- PAPER 26
 PARTICLE DENSITY 27
 PARTICLE SIZE 26, 27
 pH 27
 PHOSPHATE SLAG 39
 PLASTER 27
 PLASTIC 27
 PLATINUM GROUP 14
 POLYETHYLENE 28
 POLYVINYL CHLORIDE 28
 POROUS MATERIALS 27, 28
- RED SLURRY 41
 REFRACTORIES 29, 30, 31
 REFRACTORY 31
 RICE STRAW ASH 31
 ROAD DUST 14
 ROCKWELL 22
 RoHS/WEEE 28
- SAND 21, 32
 SAND FOR SLIDING GATES 32
 SHORE 22
 SIEVING 5
 SILICA BRICK 32
 SILICA POWDER 32
 SILICATE 33
 SILICEOUS MATERIAL 32
 SILICOALUMINUM 34
 SILICOBARIUM 34
 SILICOCALCIUM 35
 SILICOCROMIUM 35
 SILICOMANGANESE 35
 SILICON CARBIDE 33, 34, 42
 SILICON CARBIDE REFRACTORY 31
 SILICON METAL 33
 SILICON NITRIDE 34
 SILICOZIRCONIUM 36
 SIZE DISTRIBUTION 27
 SLAG 36, 37, 38, 39, 40
 SLUDGE 40
 SLURRY 41
 SODA ASH 41
 STEEL MAKING SLAG 36
 SURFACE AREA 27, 28, 42
 SYNTHETIC SILICATE 33
- TENSILE CREEP 42
 TENSILE STRENGTH 42
 THERMOSTIL 31
 TIN SLAG 39
 TITANIUM DIBORIDE 26
 TITANIUM DIOXIDE 26
 TITANIUM SLAG 39
 TUNDISH SLAG 40
 TUNGSTEN CARBIDE 42
- URANIUM 21
 URBAN AEROSOLS 43
 URBAN PARTICULATE MATTER 43
- VACUUM SLAG 40
 VANADIUM NITROGEN ALLOY 43
 VANADIUM PENTOXIDE 26
 VANADIUM SLAG 40
 VICKERS 22
- WASTE 22
 WELDING FLUX 20
- YTTRIUM 43
- ZINC OXIDE 26
 ZIRCON CONCENTRATE 43
 ZIRCON REFRACTORIES 31
 ZIRCON SAND 32
 ZIRCONIA 43
 ZIRCONIUM OXIDE 43

CRM ACID BASE ACCOUNTING

| certified values | | informational values listed in mass % | | | | | | | | | | | | | | | | 100 g units | |
|------------------|----------|---------------------------------------|-------|------|-----------------|-----------------|------|------|------|------|-------|------|------|----------------------|-------|-------|------|-------------|-------|
| Number | Total S% | Al | Ba | C | CO ₂ | CO ₃ | Ca | Fe | K | Mg | Mn | Na | P | S as SO ₄ | Si | Ti | LOI | LOM | Total |
| CAN NBM-1 | 0.28 | 7.86 | 0.117 | 0.79 | . | 0.50 | 2.30 | 4.09 | 2.36 | 1.39 | 0.046 | 2.70 | 0.10 | 0.02 | 28.47 | 0.335 | 3.45 | 0.32 | 98.38 |
| CAN KZK-1 | 0.80 | 7.37 | 0.27 | 0.95 | 3.37 | 4.22 | 1.80 | 3.30 | 3.55 | 0.95 | 0.07 | 1.18 | 0.08 | 0.01 | 29.38 | 0.35 | . | . | . |

values listed in kgCaCO₃/t

| Number | Paste PH | Acid Producing Potential | | Neutralization Potential | | Fizz Rating | |
|-----------|----------|--------------------------|----------------|--------------------------|----------|-----------------------|----------|
| | | Sobek | Modified Sobek | Sobek Slight | Moderate | Modified Sobek Slight | Moderate |
| CAN NBM-1 | 8.45 | 8.73 | 8.46 | (49.6) | (70.9) | (46.6) | (52.3) |
| CAN KZK-1 | (8.8) | 24.9 | (24.6) | 59.0 | 64.8 | 58.9 | (61.6) |

CRM AIR PARTICULATE ON FILTER MEDIA

SRM 278e is supplied as 2 loaded + 2 blank filters, analysis in ng, good for nondestructive analysis

| Number | Al | As | Ba | Ca | Ce | Co | Cr | Cu | Fe | K | Mg | Mn | Na | Ni | Pb |
|-----------------|-------|------|-------|-------|--------|--------|------|-----|-------|------|------|-----|------|-----|-------|
| SRM 2783 blank | (30) | . | (0.4) | . | . | (0.04) | (70) | . | . | . | . | . | (15) | (8) | (0.4) |
| SRM 2783 loaded | 23210 | 11.8 | 335 | 13200 | (23.4) | 7.7 | 135 | 404 | 26500 | 5280 | 8620 | 320 | 1860 | 68 | 317 |

| Number | Rb | S | Sb | Sc | Si | Sm | Th | Ti | U | V | W | Zn |
|-----------------|--------|--------|------|--------|---------|--------|--------|------|---------|------|-------|------|
| SRM 2783 blank | . | (100) | . | . | . | . | . | . | . | . | . | (50) |
| SRM 2783 loaded | (24.0) | (1050) | 71.8 | (3.54) | (58600) | (2.04) | (3.23) | 1490 | (1.234) | 48.5 | (5.0) | 1790 |

CRM ATTRITION INDEX

| Number | Attrition Index (AI units) | Standard Deviation | Uncertainty @ 95% CL | Units |
|-----------|----------------------------|--------------------|----------------------|-------|
| ASCRM 025 | 18.8 | ± 1.3 | ± 2.6 | 750 g |

RM CALCIUM ALUMINATE

typical analysis

100 g

| Number | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | MoO ₃ | S | SiO ₂ | SrO | TiO ₂ | V ₂ O ₅ |
|----------|--------------------------------|-------|--------------------------------|--------------------------------|------------------|-------|-------|------------------|-------|------------------|-------|------------------|-------------------------------|
| DH X0101 | 72.2 | 26.74 | 0.006 | 0.118 | . | 0.191 | 0.008 | . | 0.011 | 0.17 | . | . | <0.005 |
| DH X0103 | 68.8 | 23.38 | 0.028 | 0.289 | 0.296 | 3.53 | 0.024 | 0.014 | . | 0.450 | 0.009 | 0.067 | 2.36 |
| DH X0102 | 64.30 | 18.34 | 0.054 | 0.708 | . | 12.54 | 0.114 | . | 0.020 | 2.02 | 0.024 | 0.165 | 1.48 |

CRM CALCIUM CARBONATE

certified analysis in mass % and mg/kg

analysis in mg/kg

BAM: 100 g

SRM: 20 g units

| Number | certified analysis in mass % and mg/kg | | | | | | | | | | | | analysis in mg/kg | | | | | | | | | | | | |
|----------|--|-----------------|---------|------|----|----|----|------|-----|------|-------|----|------------------------|----|------|----|------|-----|------|----|------|-----|----|------|------|
| | CaCO ₃ | CO ₃ | Ca | Ba | Cr | Cu | Fe | Mg | Mn | Na | Sr | Zn | Al | B | Cd | Co | Ga | K | La | Ni | Pb | Si | Sn | Ti | Zr |
| SRM 915b | 99.907 | 59.923 | 40.0104 | (2) | . | . | . | (40) | . | (17) | (150) | . | Cl: (8) P: (3) S: (30) | | | | | | | | | | | | |
| BAM RS 3 | 99.79 | . | . | 45.3 | <1 | <1 | <5 | 183 | 3.0 | 47.5 | 173 | <2 | <5 | <1 | <0.5 | <1 | <1.5 | <20 | <0.5 | <3 | <0.1 | <20 | <1 | <0.5 | <0.2 |

CRM CASTING POWDER

analysis listed in mass %

powder 30 g

| Number | Al ₂ O ₃ | BaO | C | CaO | CO ₂ | Cr ₂ O ₃ | F | Fe ₂ O ₃ | K ₂ O | MgO | MnO ₂ | Na ₂ O | P ₂ O ₅ | SiO ₂ | TiO ₂ | ZnO | ZrO ₂ |
|------------|--------------------------------|-------|---------|-------|-----------------|--------------------------------|-----|--------------------------------|------------------|------|------------------|-------------------|-------------------------------|------------------|------------------|---------|------------------|
| FLX CRM127 | 7.82 | 0.301 | (0.12) | 34.85 | (0.4) | 0.021 | 8.7 | 0.57 | 0.09 | 2.59 | 0.032 | 10.45 | 0.037 | 37.27 | 0.241 | 0.079 | 0.016 |
| FLX CRM124 | 7.36 | 0.287 | (7.10) | 32.83 | (10.4) | 0.009 | 5.2 | 1.73 | 0.36 | 0.90 | 3.845 | 5.84 | 0.111 | 28.26 | 0.337 | 0.010 | 0.020 |
| FLX CRM125 | 7.12 | 0.207 | (9.14) | 32.07 | (12.0) | 0.011 | 4.6 | 0.77 | 0.21 | 0.95 | 0.259 | 3.90 | 0.065 | 33.29 | 0.216 | (0.005) | 0.018 |
| FLX CRM126 | 5.49 | 0.061 | (15.83) | 23.72 | (6.5) | 0.008 | 4.5 | (1.19) | 0.36 | 2.47 | 0.082 | 7.84 | 0.066 | 33.45 | 0.330 | (0.007) | 0.020 |
| FLX CRM123 | 4.63 | 0.265 | (6.30) | 29.82 | (7.4) | 0.018 | 6.6 | 1.69 | 0.41 | 2.75 | 0.041 | 7.84 | 0.095 | 35.56 | 0.202 | 0.010 | 0.021 |

RM PORTLAND CEMENT WITH EXTENSIVE ANALYSIS analysis listed in mass %

| Number | Al ₂ O ₃ | BaO | CaO | T.Fe ₂ O ₃ | K ₂ O | MgO | MnO | Ni | P ₂ O ₅ | SiO ₂ | Sr | TiO ₂ | Zr | | | | |
|--|--------------------------------|--------|-------|----------------------------------|------------------|--------|-------|-----------|-------------------------------|------------------|---------|------------------|---------|-------|--------|------|------|
| IAG OPC-1 | 4.55 | 0.0512 | 62.9 | 3.19 | 0.344 | 2.58 | 0.404 | (0.00870) | (0.044) | 21.85 | 0.01182 | 0.318 | 0.00812 | | | | |
| continued analysis listed in mg/kg ~35 g units | | | | | | | | | | | | | | | | | |
| Number | As | Be | Ce | Co | Cs | Cu | Dy | Er | Eu | Ga | Gd | Hf | Ho | La | Li | Lu | |
| IAG OPC-1 | (4.6) | (0.82) | 48.9 | 21.4 | 1.00 | (23.7) | 2.87 | 1.52 | 1.00 | 7 | 3.75 | 2.12 | 0.55 | 25.9 | (13.1) | 0.20 | |
| Number | Nb | Nd | Pb | Pr | Rb | Sb | Sm | Ta | Tb | Th | Tm | U | V | W | Y | Yb | Zn |
| IAG OPC-1 | 4.9 | 24.7 | (7.2) | 6.2 | 14.7 | 0.26 | 4.5 | (0.35) | 0.52 | 3.93 | 0.21 | 0.83 | (64.0) | (0.7) | 15.5 | 1.34 | 27.8 |

CEMENT chart 1 of 2

= class, where 1 = CRM and 2 = RM

analysis listed in mass %

| # | Number | CaO | SiO ₂ | Al ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SO ₃ | SrO | TiO ₂ | LOI | Units |
|---|--------------|--------|------------------|--------------------------------|--------------------------------|------------------|--------|-------------------|-------------------------------|-----------------|----------|------------------|-----------|---------|
| 1 | BCS 354 | 70.0 | 21.8 | 4.84 | 0.30 | 0.11 | 0.42 | 0.10 | 0.12 | 2.25 | 0.11 | (0.04) | . | 100 g |
| 2 | FLX 138 | 68.6 | 19.0 | 4.39 | 1.78 | 0.77 | 1.09 | 0.15 | 0.114 | 3.44 | 0.189 | 0.220 | . | 30 g |
| 1 | FLX CRM110 | 68.13 | 22.01 | 4.70 | 0.18 | 0.94 | 0.65 | 0.05 | 0.037 | 2.88 | 0.041 | 0.170 | (3.46) | 30 g |
| 1 | SRM 1886a | 67.87 | 22.38 | 3.875 | 0.152 | 0.093 | 1.932 | 0.021 | 0.022 | 2.086 | (0.018) | 0.084 | (1.56) | 4 x 5 g |
| 1 | FLX CRM107 | 67.19 | 21.81 | 4.23 | 1.29 | 0.70 | 0.70 | 0.18 | 0.160 | 3.13 | 0.151 | 0.194 | (6.59) | 30 g |
| 1 | FLX CRM109 | 66.45 | 20.39 | 4.25 | 2.32 | 1.06 | 1.59 | 0.18 | 0.052 | 3.11 | 0.144 | 0.203 | (5.96) | 30 g |
| 2 | JCA RM 611 | 66.25 | 21.84 | 5.41 | 3.20 | 0.34 | 1.08 | 0.40 | 0.59 | 0.25 | 0.28 | 0.30 | (0.51) | 30 g |
| 1 | SRM 1886b | 66.05 | 22.08 | 3.903 | 0.297 | 0.0164 | 1.526 | 0.01682 | 0.0413 | 2.757 | 0.0886 | 0.2054 | (2.174) | 5 x 5 g |
| 1 | FLX CRM106 | 66.05 | 20.29 | 5.70 | 1.98 | 0.86 | 0.96 | 0.12 | 0.111 | 3.01 | 0.206 | 0.271 | (2.06) | 30 g |
| 1 | NCS DC62103h | 65.78 | 22.07 | 4.56 | 3.22 | 0.69 | 2.11 | 0.16 | . | 0.32 | . | 0.37 | 0.60 | 20 g |
| 2 | TL 1Ca | 65.77 | 20.23 | 5.24 | 2.00 | 0.28 | 1.13 | 0.19 | 0.57 | 3.06 | 0.05 | 0.20 | (1.39) | 40 g |
| 1 | NCS DC62117 | 65.71 | 20.49 | 4.61 | 0.26 | 0.05 | 0.14 | 0.05 | . | 1.9 | . | 0.12 | 6.43 | 20 g |
| 1 | FLX CRM105 | 65.24 | 20.84 | 4.27 | 2.50 | 1.24 | 1.57 | 0.21 | 0.053 | 3.37 | 0.146 | 0.179 | (2.46) | 30 g |
| 1 | JCA CRM-1 | 65.21 | 20.99 | 5.26 | 2.67 | 0.56 | 2.13 | 0.26 | 0.28 | 2.05 | 0.05 | 0.35 | (0.63) | 60 g |
| 1 | FLX CRM108 | 65.15 | 20.06 | 4.66 | 2.97 | 0.74 | 2.15 | 0.09 | 0.169 | 3.31 | 0.083 | 0.186 | (2.68) | 30 g |
| 1 | SRM 634a | 65.07 | 20.493 | 5.015 | 3.362 | 0.3572 | 1.0057 | 0.0842 | 0.1767 | 2.780 | (0.0735) | 0.2463 | (1.66) | 100 g |
| 1 | BCS 353 | 64.8 | 20.5 | 3.77 | 4.82 | 0.49 | 2.42 | 0.10 | 0.077 | 2.25 | 0.23 | 0.16 | . | 100 g |
| 2 | FLX 137 | 64.77 | 20.78 | 4.99 | 3.07 | 0.769 | 1.64 | . | 0.171 | 3.17 | 0.076 | 0.221 | . | 30 g |
| 1 | FLX CRM100 | 64.51 | 20.89 | 5.54 | 2.62 | 0.82 | 1.47 | 0.23 | 0.166 | 2.97 | 0.286 | 0.283 | 2.37 | 50 g |
| 1 | CCRL 205 | 64.44 | 20.91 | 4.32 | 3.00 | 0.545 | 1.50 | 0.169 | 0.071 | 2.97 | 0.171 | 0.24 | 1.51 | 30 g |
| 2 | JCA 211S | 64.25 | 20.57 | 5.60 | 2.51 | 0.40 | 1.26 | 0.27 | 0.14 | 2.10 | . | 0.30 | 2.19 | 40 g |
| 1 | SRM 1880b | 64.16 | 20.42 | 5.183 | 3.681 | 0.646 | 1.176 | 0.0914 | 0.2443 | 2.710 | (0.0272) | 0.236 | (1.666) | 4 x 5 g |
| 1 | SRM 633a | 64.129 | 22.38 | 2.911 | 3.738 | 0.391 | 1.1532 | 0.203 | 0.14263 | 2.178 | (0.0507) | 0.2157 | (2.460) | 4 x 5 g |
| 1 | FLX CRM113 | 63.63 | 20.98 | 5.06 | 2.75 | 0.619 | 2.49 | (0.092) | 0.135 | 2.47 | 0.064 | 0.231 | (1.53) | 30 g |
| 1 | SRM 1888b | 63.13 | 20.42 | 4.277 | 3.062 | 0.658 | 3.562 | 0.1364 | 0.07307 | 2.634 | 0.1099 | 0.2316 | (various) | 4 x 5 g |
| 2 | JCA RM 613 | 63.00 | 19.51 | 5.36 | 2.78 | 1.20 | 1.07 | 0.23 | 0.15 | 6.07 | 0.15 | 0.35 | (3.45) | 30 g |
| 2 | JCA RM 612 | 62.95 | 20.12 | 5.19 | 2.81 | 0.90 | 1.52 | 0.52 | 1.02 | 4.51 | 0.045 | 0.28 | (2.52) | 30 g |
| 2 | CCRL 173 | 62.45 | 20.01 | 4.49 | 2.62 | 0.447 | 3.03 | 0.309 | 0.192 | 4.10 | . | 0.27 | 2.02 | 30 g |
| 2 | CCRL 174 | 62.43 | 20.75 | 3.71 | 3.62 | 0.430 | 4.83 | 0.189 | 0.067 | 2.64 | . | 0.21 | 1.14 | 30 g |
| 1 | NCS DC62101c | 62.23 | 20.41 | 4.68 | 3.20 | 0.71 | 2.66 | 0.12 | . | 3.16 | . | 0.27 | 2.18 | 20 g |
| 1 | CCRL 206 | 61.90 | 19.64 | 5.16 | 3.37 | 0.653 | 3.84 | 0.148 | 0.044 | 3.47 | 0.087 | 0.28 | 1.55 | 30 g |
| 1 | SRM 1885b | 61.87 | 20.05 | 4.70 | 3.044 | 0.497 | 3.86 | 0.293 | 0.0737 | 2.832 | 0.0795 | 0.2361 | (2.310) | 5 x 5 g |
| 2 | CCRL 172 | 61.78 | 19.30 | 4.66 | 2.93 | 0.935 | 4.74 | 0.267 | 0.111 | 3.21 | . | 0.26 | 2.09 | 30 g |

| Number | CO ₂ | Free CaO | Cl | Cr ₂ O ₃ | F | Mn | MnO | Mn ₂ O ₃ | S | ZnO | Ins. Res. |
|--------------|-----------------|----------|----------|--------------------------------|----------|--------------|-------|--------------------------------|----------|-----------|-------------------------------|
| BCS 354 | . | . | . | . | . | . | . | 0.058 | . | . | . |
| FLX 138 | . | . | . | . | . | . | . | 0.095 | . | 0.017 | . |
| FLX CRM110 | . | . | (0.008) | 0.004 | . | . | 0.029 | . | . | 0.003 | . |
| SRM 1886a | . | . | (0.0042) | 0.0024 | (0.02) | . | . | 0.0073 | . | (0.001) | (0.23) |
| FLX CRM107 | . | . | 0.043 | 0.006 | . | . | 0.040 | . | . | 0.013 | . |
| FLX CRM109 | . | . | 0.049 | 0.008 | . | . | 0.051 | . | . | 0.042 | . |
| JCA RM 611 | . | . | . | . | . | . | 0.06 | . | . | . | . |
| SRM 1886b | . | (0.24) | 0.00399 | 0.00404 | (0.0118) | . | . | 0.02639 | . | (0.00058) | (0.13) |
| FLX CRM106 | . | . | 0.055 | 0.008 | . | . | 0.161 | . | . | 0.012 | . |
| NCS DC62103h | . | . | . | . | . | . | . | . | . | . | 0.12 |
| TL 1Ca | . | (0.83) | . | . | . | . | . | . | . | . | (0.21) |
| NCS DC62117 | . | . | . | . | . | . | . | . | . | . | . |
| FLX CRM105 | . | . | 0.049 | 0.008 | . | . | 0.040 | . | . | 0.054 | . |
| JCA CRM-1 | . | . | . | . | . | 0.06 | . | . | . | . | . |
| FLX CRM108 | . | . | 0.042 | 0.007 | . | . | 0.219 | . | . | 0.036 | . |
| SRM 634a | . | (1.86) | . | (0.0114) | . | . | . | (0.0229) | . | (0.0222) | (0.21) |
| BCS 353 | . | . | . | . | . | . | . | 0.23 | . | . | . |
| FLX 137 | . | . | . | . | . | . | . | 0.266 | . | 0.029 | . |
| FLX CRM100 | . | . | (0.09) | 0.009 | . | . | . | 0.066 | . | 0.051 | . |
| CCRL 205 | . | 0.74 | 0.008 | 0.011 | . | . | . | 0.121 | . | 0.010 | 0.34 |
| JCA 211S | . | . | 0.022 | . | . | . | 0.05 | . | . | . | 0.23 |
| SRM 1880b | . | (2.227) | 0.01830 | 0.01927 | (0.0539) | . | . | 0.1981 | (0.0131) | (0.01054) | (0.487) |
| SRM 633a | . | (1.60) | 0.0087 | (0.0124) | (0.038) | BaO: (0.256) | . | 0.1176 | (0.049) | 0.123 | (0.23) |
| FLX CRM113 | . | . | . | (0.007) | . | . | . | 0.233 | (0.137) | 0.030 | . |
| SRM 1888b | . | (1.42) | 0.0143 | (0.01021) | (0.048) | . | . | 0.0652 | (0.15) | (0.01253) | (0.32) |
| JCA RM 613 | . | . | . | . | . | . | 0.08 | . | . | . | . |
| JCA RM 612 | . | . | . | . | . | . | 0.06 | . | . | . | . |
| CCRL 173 | 0.6 | 1.65 | 0.023 | 0.009 | . | . | . | 0.060 | . | 0.024 | 0.36 last, expires 12/31/2014 |
| CCRL 174 | . | 1.04 | 0.005 | 0.006 | . | . | . | 0.073 | . | 0.014 | 0.26 last, expires 12/31/2014 |
| NCS DC62101c | . | . | . | . | . | . | . | . | . | . | 1.20 |
| CCRL 206 | . | 0.97 | 0.014 | 0.011 | . | . | . | 0.029 | . | 0.020 | 0.22 |
| SRM 1885b | . | (0.27) | (0.0021) | 0.02709 | (0.0524) | . | . | 0.1282 | (0.042) | 0.0354 | (0.36) BaO: 0.0149 |
| CCRL 172 | 1.38 | 1.04 | 0.008 | 0.009 | . | . | . | 0.088 | . | 0.004 | 0.44 last, expires 12/31/2014 |

CEMENT chart 2 of 2

= class, where 1 = CRM and 2 = RM

analysis listed in mass %

| # | Number | CaO | Ca | SiO ₂ | Al ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SO ₃ | SrO | TiO ₂ | LOI | Units |
|---|--------------|-------|-------|------------------|--------------------------------|--------------------------------|------------------|-------|-------------------|-------------------------------|-----------------|---------|------------------|-----------|-----------|
| 1 | SRM 1884b | 61.31 | . | 19.30 | 4.851 | 2.937 | 0.957 | 4.74 | 0.278 | 0.0965 | 4.034 | 0.0258 | 0.2651 | (various) | 4 x 4.5 g |
| 1 | SRM 1887b | 61.15 | . | 19.59 | 4.911 | 2.471 | 0.961 | 3.624 | 0.288 | 0.1540 | 4.599 | 0.2625 | 0.2034 | 2.121 | 5 x 4 g |
| 1 | NCS DC62118 | 60.99 | . | 21.73 | 4.75 | 4.12 | 0.43 | 4.37 | 0.12 | . | 2.27 | . | 0.23 | 0.81 | 20 g |
| 1 | SRM 1889b | 60.11 | . | 18.39 | 5.79 | 2.891 | 1.115 | 2.82 | 0.36 | 0.297 | 4.3721 | 0.284 | 0.260 | (3.117) | 5 x 5 g |
| 1 | FLX CRM122 | 59.00 | . | 21.94 | 5.60 | 1.67 | 0.900 | 2.02 | 0.204 | 0.066 | (3.27) | 0.131 | 0.353 | (4.86) | 30 g |
| 1 | NCS DC62102d | 57.25 | . | 23.54 | 6.60 | 3.25 | 0.71 | 1.88 | 0.15 | . | 2.74 | . | 0.39 | 3.13 | 20 g |
| 2 | TL 203Ca | 57.09 | . | 18.75 | 4.78 | 2.34 | 0.93 | 4.80 | 0.21 | (0.13) | 3.28 | (0.06) | 0.24 | (7.24) | 40 g |
| 1 | FLX CRM130 | 56.60 | . | 14.35 | 11.62 | 2.88 | 0.682 | 1.84 | 0.277 | 0.067 | 10.91 | 0.052 | 0.563 | (5.12) | 30 g |
| 1 | JCA CRM-2 | 56.33 | . | 25.66 | 8.94 | 2.08 | 0.31 | 3.05 | 0.24 | 0.07 | (2.59) | 0.07 | 0.50 | (0.47) | 60 g |
| 1 | FLX CRM118 | 55.47 | . | 21.57 | 6.94 | 3.88 | 0.983 | 2.08 | 0.181 | 0.163 | (2.97) | 0.087 | 0.301 | (5.02) | 30 g |
| 1 | FLX CRM103 | 54.90 | . | 26.95 | 7.75 | 1.78 | 0.77 | 4.44 | 0.33 | 0.09 | 2.73 | 0.070 | 0.372 | (0.59) | 50 g |
| 1 | TL 201C | 54.48 | . | 25.63 | 6.81 | 2.08 | 0.73 | 3.35 | 0.32 | . | 3.16 | . | . | 1.96 | 40 g |
| 1 | FLX CRM117 | 54.22 | . | 28.05 | 7.68 | 1.37 | 0.916 | 2.95 | 0.223 | (0.035) | (3.49) | 0.098 | 0.564 | (0.402) | 30 g |
| 1 | FLX CRM115 | 53.93 | . | 27.29 | 8.13 | 0.972 | 0.612 | 2.93 | 0.194 | 0.073 | (3.90) | 0.117 | 0.610 | (1.56) | 30 g |
| 1 | FLX CRM119 | 53.78 | . | 24.04 | 7.36 | 3.04 | 1.23 | 1.34 | 0.292 | 0.153 | (2.43) | 0.158 | 0.347 | (6.42) | 30 g |
| 1 | FLX CRM120 | 53.63 | . | 26.51 | 7.64 | 1.11 | 0.692 | 3.55 | 0.184 | 0.086 | (4.35) | 0.124 | 0.598 | (1.71) | 30 g |
| 1 | FLX CRM114 | 51.29 | . | 28.61 | 6.94 | 1.37 | 0.954 | 4.93 | 0.277 | 0.027 | (4.07) | 0.108 | 0.525 | (1.28) | 30 g |
| 1 | FLX CRM116 | 50.05 | . | 30.81 | 9.11 | 0.86 | 0.693 | 4.47 | 0.201 | 0.034 | (3.14) | 0.091 | 0.690 | (-0.394) | 30 g |
| 1 | TL 200Ca | 49.97 | . | 26.55 | 8.72 | 4.07 | 1.10 | 2.06 | 0.21 | 0.45 | 2.84 | 0.13 | 0.46 | (3.30) | 40 g |
| 1 | SRM 1881b | 49.27 | . | 29.045 | 8.812 | 3.365 | 0.721 | 2.741 | 0.790 | 0.0510 | 2.72 | 0.0836 | 0.3011 | (1.699) | 4 x 5 g |
| 1 | FLX CRM101 | 48.43 | . | 30.29 | 8.86 | 3.54 | 2.10 | 1.72 | 0.657 | 0.189 | 3.12 | 0.248 | 0.463 | (5.13) | 50 g |
| 2 | DH X0210 | 46.72 | 33.39 | 30.30 | 9.99 | 1.66 | 0.541 | 4.96 | 0.236 | 0.066 | . | 0.077 | 0.421 | . | 100 g |
| 1 | TL 202C | 45.12 | . | 29.61 | 10.14 | 3.27 | 1.05 | 4.46 | 0.32 | . | 3.17 | . | . | 1.51 | 40 g |
| 1 | NCS DC62116a | 44.73 | . | 27.12 | 10.99 | 3.43 | 1.10 | 2.82 | 0.23 | . | 2.49 | . | 0.57 | 6.00 | 20 g |
| 1 | FLX CRM121 | 44.45 | . | 30.81 | 8.78 | 3.04 | 2.73 | 1.41 | 1.23 | 0.112 | 2.55 | 0.140 | 0.344 | (4.24) | 30 g |
| 1 | FLX CRM131 | 42.89 | . | 8.73 | 23.10 | 3.24 | 0.287 | 1.62 | 0.466 | 0.060 | 18.19 | 0.067 | 1.15 | (2.03) | 30 g |
| 1 | NCS DC62105g | 42.40 | . | 13.76 | 3.12 | 2.13 | 0.50 | 2.11 | 0.14 | . | 0.21 | . | 0.18 | 34.94 | 20 g |
| 1 | NCS DC62104b | 39.65 | . | 12.77 | 3.66 | 2.02 | 0.93 | 1.20 | 0.18 | . | 0.72 | . | 0.18 | 38.32 | 20 g |
| 1 | SRM 1882a | 39.29 | . | 4.01 | 39.14 | 14.67 | 0.051 | 0.51 | 0.021 | (0.070) | . | (0.024) | 1.786 | (0.20) | 4 x 5 g |
| 1 | SRM 1883a | 29.52 | . | 0.24 | 70.04 | 0.078 | 0.014 | 0.19 | 0.30 | (0.003) | . | (0.019) | (0.020) | (0.59) | 4 x 5 g |
| 2 | DH X0209 | . | 48.78 | 21.95 | 4.63 | 0.204 | 1.01 | 0.717 | 0.078 | 0.043 | . | 0.051 | 0.095 | . | 100 g |
| 2 | DH X0212 | . | 46.48 | 21.16 | 4.41 | 3.94 | 0.495 | 0.945 | 0.084 | 0.191 | . | 0.086 | 0.242 | . | 100 g |
| 2 | DH X0211 | . | 40.63 | 25.04 | 6.86 | 2.98 | 0.524 | 2.79 | 0.156 | 0.137 | . | 0.083 | 0.319 | . | 100 g |

| Number | BaO | Free CaO | Cl | Cr ₂ O ₃ | F | Mn | Mn ₂ O ₃ | S | Unignited SO ₃ | V ₂ O ₅ | ZnO | Ins. Res. |
|--------------|---------|----------|----------|--------------------------------|----------|------|--------------------------------|----------------------------------|--|-------------------------------|---------|-----------|
| SRM 1884b | . | (0.418) | (0.0065) | 0.00791 | (0.0394) | . | 0.0750 | . | . | . | . | (0.159) |
| SRM 1887b | (0.022) | 0.21 | 0.01001 | 0.01551 | 0.101 | . | 0.0957 | 0.025 | . | . | 0.01560 | 0.26 |
| NCS DC62118 | . | . | . | . | . | . | . | . | . | . | . | 1.18 |
| SRM 1889b | . | (0.52) | 0.0101 | 0.0083 | (0.10) | . | 0.084 | . | . | . | 0.0770 | (0.30) |
| FLX CRM122 | . | . | . | (0.004) | . | . | 0.111 | (0.288) | SO ₄ as SO ₃ : 2.62 | . | 0.027 | . |
| NCS DC62102d | . | . | . | . | . | . | . | . | . | . | . | . |
| TL 203Ca | . | (1.23) | . | . | . | . | . | . | . | . | . | (0.58) |
| FLX CRM130 | . | . | . | 0.021 | . | . | 0.062 | . | SO ₄ as SO ₃ : (10.18) | . | 0.018 | . |
| JCA CRM-2 | . | . | . | . | . | 0.15 | . | (0.32) | (1.91) | . | . | . |
| FLX CRM118 | . | . | . | 0.009 | . | . | 0.176 | (0.131) | SO ₄ as SO ₃ : 2.89 | . | 0.053 | . |
| FLX CRM103 | . | . | (0.040) | 0.007 | . | . | 0.170 | . | . | . | 0.014 | . |
| TL 201C | . | . | 0.06 Cl- | . | . | . | . | 0.31 S ₂ ⁻ | . | . | . | . |
| FLX CRM117 | . | . | . | (0.004) | . | . | 0.193 | (0.602) | SO ₄ as SO ₃ : 1.80 | . | 0.024 | . |
| FLX CRM115 | . | . | . | (0.004) | . | . | 0.214 | (0.611) | SO ₄ as SO ₃ : 2.41 | . | 0.008 | . |
| FLX CRM119 | . | . | . | 0.008 | . | . | 0.040 | (0.029) | SO ₄ as SO ₃ : 2.38 | . | 0.040 | . |
| FLX CRM120 | . | . | . | (0.005) | . | . | 0.188 | (0.632) | SO ₄ as SO ₃ : 3.18 | . | 0.008 | . |
| FLX CRM114 | . | . | . | (0.005) | . | . | 0.154 | (0.676) | SO ₄ as SO ₃ : 2.50 | . | 0.022 | . |
| FLX CRM116 | . | . | . | (0.005) | . | . | 0.238 | (0.737) | SO ₄ as SO ₃ : 1.22 | . | 0.015 | . |
| TL 200Ca | . | (0.34) | . | . | . | . | . | . | . | . | . | (17.42) |
| SRM 1881b | 0.191 | (1.16) | 0.0081 | 0.00949 | (0.09) | . | 0.1175 | . | . | . | 0.1198 | (16.4) |
| FLX CRM101 | . | . | . | 0.011 | . | . | 0.116 | . | . | . | 0.046 | . |
| DH X0210 | 0.071 | . | . | . | . | . | 0.327 | 1.77 | . | 0.011 | . | . |
| TL 202C | . | . | 0.01 Cl- | . | . | . | . | 0.20 S ₂ ⁻ | . | . | . | . |
| NCS DC62116a | . | . | . | . | . | . | . | . | . | . | . | . |
| FLX CRM121 | . | . | . | 0.008 | . | . | 0.111 | (0.025) | SO ₄ as SO ₃ : 2.37 | . | 0.031 | . |
| FLX CRM131 | . | . | . | 0.038 | . | . | 0.029 | . | SO ₄ as SO ₃ : (17.85) | . | 0.006 | . |
| NCS DC62105g | . | . | . | . | . | . | . | . | . | . | . | . |
| NCS DC62104b | . | . | . | . | . | . | . | . | . | . | . | . |
| SRM 1882a | . | . | . | (0.113) | . | . | (0.060) | . | . | . | (0.004) | . |
| SRM 1883a | . | . | . | (0.006) | . | . | (0.003) | . | . | . | . | . |
| DH X0209 | 0.028 | . | . | . | . | . | 0.025 | 1.19 | . | . | . | . |
| DH X0212 | . | . | . | . | . | . | 0.062 | 1.18 | . | . | . | . |
| DH X0211 | 0.041 | . | . | . | . | . | 0.172 | 1.48 | . | 0.014 | . | . |

RM

CEMENT SET JCA 601B

available in set/15 only

number 1-14 powder 20 g number 15 powder 30 g

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | SrO | TiO ₂ |
|-------------|--------------------------------|-------|--------------------------------|------------------|-------|-------|-------------------|-------------------------------|-----------------|------------------|--------|------------------|
| JCA 601B 1 | 5.10 | 59.07 | 2.37 | 1.191 | 1.31 | 0.404 | 0.809 | 1.479 | 9.83 | 17.61 | (0.38) | 0.282 |
| JCA 601B 2 | 5.98 | 65.23 | 2.90 | 0.370 | 1.37 | 0.097 | 0.272 | 0.119 | 2.34 | 20.75 | (0.06) | 0.299 |
| JCA 601B 3 | 4.45 | 63.18 | 2.37 | 0.696 | 2.36 | 0.050 | 0.373 | 0.758 | 5.97 | 19.11 | (0.19) | 0.222 |
| JCA 601B 4 | 5.23 | 66.63 | 2.77 | 0.281 | 0.87 | 0.090 | 0.166 | 0.198 | 2.70 | 20.36 | (0.04) | 0.333 |
| JCA 601B 5 | 3.28 | 64.96 | 3.37 | 0.355 | 1.02 | 0.058 | 0.199 | 0.111 | 1.79 | 24.41 | (0.03) | 0.189 |
| JCA 601B 6 | 3.87 | 64.15 | 4.13 | 0.413 | 0.83 | 0.105 | 0.223 | 0.126 | 2.20 | 23.18 | (0.03) | 0.201 |
| JCA 601B 7 | 2.82 | 63.36 | 3.14 | 0.398 | 0.80 | 0.037 | 0.153 | 0.096 | 2.57 | 26.27 | (0.03) | 0.164 |
| JCA 601B 8 | 2.73 | 63.49 | 3.05 | 0.304 | 0.59 | 0.104 | 0.186 | 0.169 | 2.48 | 26.61 | (0.02) | 0.143 |
| JCA 601B 9 | 7.94 | 58.91 | 2.07 | 0.320 | 2.49 | 0.081 | 0.328 | 0.110 | . | 24.53 | (0.05) | 0.370 |
| JCA 601B 10 | 9.47 | 54.60 | 1.83 | 0.342 | 3.22 | 0.148 | 0.239 | 0.154 | . | 26.50 | (0.05) | 0.379 |
| JCA 601B 11 | 8.89 | 55.64 | 1.98 | 0.337 | 2.96 | 0.181 | 0.311 | 0.283 | . | 26.11 | (0.06) | 0.399 |
| JCA 601B 12 | 8.75 | 55.78 | 2.14 | 0.319 | 3.13 | 0.620 | 0.189 | 0.085 | . | 24.75 | (0.04) | 1.035 |
| JCA 601B 13 | 11.03 | 50.51 | 1.10 | 0.291 | 4.20 | 0.112 | 0.348 | 0.055 | . | 28.61 | (0.05) | 0.449 |
| JCA 601B 14 | 16.05 | 35.85 | 0.24 | 0.188 | 10.21 | 0.115 | 0.274 | 0.011 | . | 35.03 | (0.05) | 0.435 |
| JCA 601B 15 | 2.07 | 75.62 | 7.01 | 0.000 | 0.01 | 0.002 | 0.030 | 0.003 | 0.02 | 14.87 | (0.02) | 0.004 |

CRM CHLORINE and FLUORINE in CEMENT

| Number | Description | CaF ₂ | F | Cl- | Units |
|--------------|-----------------|------------------|------|-------|-------|
| NCS DC62121a | Cement Raw Meal | . | . | 0.016 | 20 g |
| NCS DC62122a | Cement | . | . | 0.009 | 20 g |
| NCS DC62125a | Cement | (0.37) | 0.18 | . | 20 g |

CRM PORTLAND CEMENT HEAT OF HYDRATION

| Number | Heat of Solution J/g | 7 Days J/g | 28 Days J/g | Units |
|----------|----------------------|------------|-------------|-------|
| JCA 301S | 2,483.5 | 274.0 | 325.2 | 600 g |

CRM COMPRESSIVE STRENGTH N/mm²

| Number | 3 Days | 7 Days | 28 Days | Units |
|----------|--------|--------|---------|--------|
| JCA 401J | 29.6 | 44.4 | 62.1 | 4.8 kg |

CRM CLASSIC CEMENT CHEMISTRIES

20 g units

| Number | P - Pozzolana | S - Slag | D - Limestone | D1 - CO ₂ | R5 - Unsolved Slag (EDTA) | Description |
|--------------|---------------|----------|---------------|----------------------|---------------------------|------------------------------------|
| NCS DC62119a | 9.3 | 4.5 | 2.4 | (1.50) | . | Ordinary Portland Cement |
| NCS DC62120 | 0.5 | 18.5 | 7 | 3.5 | 97.5 | Portland Blast-Furnace Slag Cement |

CRM CEMENT CLINKER PHASE ABUNDANCE

| Number | Alite | Alkali Sulfates | Aluminate | Aphthitalite | Arcanite | Belite | Ferrite | Periclase | Units |
|-----------|-------|-----------------|-----------|--------------|----------|--------|---------|-----------|----------|
| SRM 2686a | 63.53 | 0.86 | 2.46 | (0.74) | (0.27) | 18.80 | 10.80 | 3.40 | 3 x 10 g |
| SRM 2687a | 57.88 | . | 9.56 | . | . | 24.70 | 6.27 | . | 5 x 8 g |
| SRM 2688 | 64.95 | . | 4.99 | . | . | 17.45 | 12.20 | . | 3 x 10 g |

CRM PORTLAND CEMENT FINENESS AND BLAINE STANDARD

| Number | Remaining after passing through 80 micron sieve | Blaine | Density g/cm ³ | Units |
|--------------|---|--------------------------|---------------------------|-------|
| NCS DC62127e | 2.68 % | 356.8 m ² /kg | (3.05) | 150 g |
| TL 201B | . | 4,231 cm ² /g | 3.03 | 40 g |
| TL 202B | . | 4,135 cm ² /g | 2.94 | 40 g |
| JCA 1020 | . | 3,300 cm ² /g | . | 30 g |
| JCA 102N | . | 3,300 cm ² /g | . | 30 g |

TL 201B, TL 202B, JCA 1020, JCA 102N: expires March 2023
JCA 102N: expired, last of stock

CRM CEMENT FINENESS

certified analysis

informational analysis listed in mass %

46H: 10 x 5 g

TL, 114q: powder 20 x 5 g units

| Number ASTM METHOD | Surface Area | | 45 µm Sieve Residue C430-96 | Cement Composition (mass %) | | | | | | | | Other Components (mass %) | | | | | LOI | |
|-----------------------|-------------------------|-------------------------|-----------------------------|-----------------------------|------------------|------------------|-------------------|--------------------------------|------|--------------------------------|------------------|---------------------------|-------------------|-------------------------------|-----------------|------------------|------|------------------|
| | Blaine C204-96a | Wagner C115-96a | | C ₂ S | C ₃ S | C ₃ A | C ₄ AF | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | | TiO ₂ |
| SRM 114q | 3818 cm ² /g | 2183 cm ² /g | 0.79 % | 14 | 60 | 7 | 10 | 4.7 | 64.0 | 3.2 | 0.70 | 2.2 | 0.07 | 0.12 | 2.4 | 20.7 | 0.30 | 1.67 |
| TL 3BGa | 3727 cm ² /g | . | 3.14 % | 15 | 59 | 8 | 8 | 4.9 | 63.9 | 2.8 | 0.68 | 1.9 | 0.19 | 0.21 | 2.9 | 20.6 | 0.30 | 1.5 |
| SRM 46h | . | . | 7.43 % | | | | | | | | | | | | | | | |

CRM CEMENT FINENESS

particle size analysis detailed on certificates

40-100 g units

| Number | Density g/cm ³ | Blaine cm ² /g | C ₂ S | C ₃ S | C ₃ A | C ₄ AF | Al ₂ O ₃ | CaO | F.CaO | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | TiO ₂ | Insol. | LOI |
|-----------|---------------------------|---------------------------|------------------|------------------|------------------|-------------------|--------------------------------|-------|-------|--------------------------------|------------------|------|-------------------|-------------------------------|-----------------|------------------|------------------|--------|------|
| TL 2BGa | 3.15 | 4,206 | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |
| TL 9 | 3.15 | 4,175 | 12 | 62 | 7 | 9 | 4.66 | 64.00 | 1.09 | 3.01 | 0.76 | 2.20 | 0.26 | 0.07 | 2.74 | 20.47 | 0.20 | 0.45 | 1.46 |
| TL 203BGa | 3.05 | 4,329 | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . | . |

RM CEMENT FINENESS - SIEVING METHOD

| Number | 45 µm | 32 µm | 20 µm | 16 µm | 10 µm | units | percent remaining, ordinary portland cement |
|----------|-------|-------|-------|-------|-------|-------|---|
| JCA 701B | 10.4 | 22.8 | 43.0 | 52.4 | 70.1 | 30 g | last of stock |
| JCA 701C | 10.2 | 22.1 | 42.6 | 51.8 | 69.9 | 30 g | |

CRM CEMENT COMPONENT MATERIAL

analysis listed in mass %

NCS DC61106: 50g

others: 20 g units

| Number | Material | CaO | T.CaCO ₃ | Al ₂ O ₃ | SiO ₂ | F | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | S | SO ₃ | TiO ₂ | LOI |
|--------------|---------------------------------|-------|---------------------|--------------------------------|------------------|------|--------------------------------|------------------|-------|-------------------|------|-----------------|------------------|-------|
| NCS DC62110a | Portland Blast Furnace Slag | 55.21 | . | 7.24 | 24.78 | . | 3.00 | 0.71 | 2.64 | 0.18 | . | 2.47 | 0.51 | 2.70 |
| NCS DC62109 | Portland Pozzolanic | 47.57 | . | 6.52 | 32.67 | . | 3.54 | 1.43 | 1.86 | 0.85 | . | 2.59 | 0.16 | 2.44 |
| NCS DC62111 | Portland Fly Ash | 46.52 | . | 8.93 | 24.31 | . | 4.90 | 0.61 | 1.90 | 0.32 | . | 2.47 | 0.33 | 9.09 |
| NCS DC62123 | Sulphoaluminate Cement Clinker | 43.4 | . | 32.6 | 8.56 | . | 2.21 | 0.22 | 1.37 | 0.09 | . | 9.55 | 1.51 | 0.41 |
| NCS DC62126a | Cement Black Raw Meal | 39.28 | 70.3 | . | . | 0.15 | 2.07 | . | . | . | . | . | . | 38.51 |
| NCS DC62113 | Granulated Blast Furnace Slag | 35.62 | . | 12.23 | 34.93 | . | 1.26 | 0.54 | 10.66 | 0.42 | 0.61 | 1.17 | 1.06 | 1.05 |
| NCS DC62112 | Aluminate | 34.56 | . | 51.15 | 7.95 | . | 1.91 | 0.13 | 0.63 | 0.04 | 0.1 | . | 2.03 | 0.68 |
| NCS DC62124 | Sulphoaluminate Cement Raw Meal | 33.05 | . | 22.29 | 5.09 | . | 1.34 | 0.14 | 1.21 | 0.06 | . | 7.07 | 1.07 | 28.21 |
| NCS DC62115a | Fly Ash for Cement | 4.13 | . | 33.07 | 51.38 | . | 4.58 | 0.86 | 1.02 | 0.33 | . | 0.24 | 1.14 | 2.80 |
| NCS DC62114a | Pozzolana for Cement | 2.15 | . | 20.66 | 56.86 | . | 3.52 | 1.95 | 0.86 | 0.83 | . | 0.50 | 0.78 | 11.53 |
| NCS DC61106 | Albite Cement | 0.48 | . | 19.62 | 67.96 | . | 0.10 | 0.098 | 0.015 | 11.26 | . | . | 0.054 | 0.36 |

CRM MAGNETIC CENOSPHERE FROM COAL POWER PLANTS

analysis listed in mass % fine powder concentrate, 100 g units

| Number | Al ₂ O ₃ | Ba | CaO | CO ₂ | F% | FeO | T.Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | S | SiO ₂ | Sr | TiO ₂ | Zr | LOI |
|--------------|--------------------------------|-------|------|-----------------|-------|--------|----------------------------------|------------------|-------|---------|-------------------|-------------------------------|--------|------------------|--------|------------------|--------|--------|
| VS 9234-2008 | 25.9 | 0.084 | 1.19 | (0.43) | . | 1.67 | 2.85 | 2.9 | 0.86 | 0.039 | 0.48 | 0.076 | . | 63.7 | 0.0217 | 0.74 | 0.026 | (0.92) |
| VS 9235-2008 | 19.3 | 0.089 | 2.77 | 0.51 | (0.2) | (2.4) | 4.04 | 3.5 | 1.64 | 0.049 | 1.33 | 0.19 | (0.07) | 64.1 | 0.0480 | 0.78 | 0.0246 | (1.6) |
| VS 9236-2008 | 2.6 | 0.38 | 8.2 | (0.64) | . | (15.4) | (78.4) | 0.12 | (1.3) | (0.082) | 0.27 | (0.036) | 0.86 | 7.7 | 0.17 | 0.21 | 0.0039 | (0.42) |

continued analysis listed in mg/kg except

| Number | Ag | As | Be | Ce | Co | Cr | Cs | Cu | Dy | Er | Eu | Ga | Gd | Ge | Hf | Ho | La | Li | Lu | Mo | Nb |
|--------------|-------|------|-------|-----|-----|----|--------|----------|-------|-------|------|-------|-------|-------|--------|-----|------|------|-------|-------|----|
| VS 9234-2008 | . | . | 4.6 | 115 | 7.6 | 59 | 12 | 33 (6) | (3.5) | 1.5 | (11) | (7.3) | . | 7.1 | (1.2) | 60 | (42) | 0.56 | (16) | 19 | |
| VS 9235-2008 | (0.2) | (11) | 2.8 | 89 | 11 | 66 | 8.6 | 34 (5.5) | (3.0) | 1.4 | 9.3 | (6.7) | (0.9) | (6.9) | (1.1) | 46 | (55) | 0.53 | (3.4) | 15 | |
| VS 9236-2008 | . | . | (0.7) | 11 | 23 | 67 | (0.21) | 28 (1.0) | (0.6) | (0.5) | (3) | (1.2) | . | 0.97 | (0.19) | 5.7 | . | 0.10 | (4.6) | (3.1) | |

| Number | Nd | Ni | Pb | Pr | Rb | Sb | Sc | Sm | Sn | Ta | Tb | Th | Tl | Tm | U | V | Y | Yb | Zn |
|--------------|-------|----|-----------|-----|-------|-----|-----|-------|--------|--------|-----|--------|--------|-------|-----|-----|------|----|----|
| VS 9234-2008 | 48 | 30 | 17 (13) | 116 | (0.5) | 13 | 8.4 | (1.6) | (1.5) | (1.1) | 19 | (0.5) | (0.5) | 4.3 | 60 | 34 | 3.5 | 28 | |
| VS 9235-2008 | 39 | 35 | 20 (10) | 135 | (1.3) | 15 | 7.2 | (2.3) | (1.2) | 0.92 | 14 | (0.7) | (0.5) | 4.6 | 102 | 33 | 3.3 | 50 | |
| VS 9236-2008 | (5.3) | 62 | 3.5 (1.3) | 3.8 | (0.3) | 2.5 | 1.0 | (1) | (0.15) | (0.16) | 1.4 | (0.06) | (0.09) | (1.0) | 26 | 5.6 | 0.50 | 38 | |

COAL chart 1 of 2

#=class, 1=CRM and 2=RM mass % except * for mg/kg ACIRS, AS(C)RM, COCO: 50-250 g SABS: 100-150 g USZ: as shown others: 50 g

| # | Number | S | Ash | Volatile Matter | Heat in J/g or BTU/lb | Density | Moisture | C | Fixed C | Cl | F | H | Hg* | N | O | P |
|---|----------------------|-------|-------|-----------------|-----------------------|-------------|---------------|-------------|-----------|--------|--------|--------|--------|-------------------------|------------|--------------|
| 1 | NCS FC28143 | 6.40 | 31.80 | 14.00 | 20,650 J | 1.82 | . | 53.28 | . | . | . | 2.58 | . | 0.72 | . | . |
| 1 | 502-674-18191 | 4.89 | 8.80 | . | . | . | . | . | . | . | . | . | . | . | . | 17025 |
| 1 | 502-684-18321 | 4.77 | 0.22 | 14.1 | (15,379) BTU | . | . | 88.5 (85.7) | . | . | . | 3.91 | . | 1.40 | . | 17025 |
| 1 | SRM 2685c | 4.72 | . | . | Br:(4.94*) | Mg:(0.0814) | Mn:(0.003684) | . | (0.05540) | . | . | . | 0.1494 | . | . | . |
| 1 | NCS FC28009k | 4.31 | 42.34 | 16.66 | 18,060 J | 1.84 | . | 45.40 | . | . | . | 2.63 | . | 0.80 | . | . |
| 1 | NCS FC28142 | 4.35 | 33.40 | 14.38 | 21,050 J | 1.76 | . | 53.63 | . | . | . | 2.79 | . | 0.81 | . | . |
| 1 | NCS FC28221 | 4.04 | 18.98 | 32.0 | 27,790 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS FC28220 | 4.03 | 16.52 | 11.15 | 28,670 J | . | . | . | . | . | . | . | . | . | . | . |
| 2 | 502-836-10076 | 3.99 | 12.57 | (39.18) | (13,093) BTU | . | . | (70.18) | . | . | . | 4.88 | 0.097 | 1.27 | . | . |
| 1 | NCS FC28012k | 3.39 | 18.15 | 10.90 | 27,740 J | 1.57 | . | 71.11 | . | . | . | 2.97 | . | 0.94 | . | . |
| 2 | 502-837-10182 | 3.26 | 9.00 | (38.43) | (13,862) BTU | . | . | (72.64) | . | . | . | 4.88 | 0.058 | 1.36 | . | . |
| 1 | 502-673-17041 | 3.23 | 8.04 | . | . | . | . | . | . | . | . | . | . | . | . | 17025 |
| 1 | NCS FC28210 | 3.17 | 25.80 | 8.77 | 24,130 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | CZ SF-06-14 | 3.13 | 27.21 | 27.36 | 23,990 J | 10,314 BTU | . | 58.28 | . | . | . | 3.51 | . | 3.80 | . | . |
| 1 | NCS FC28141 | 2.92 | 28.64 | 11.50 | 23,040 J | 1.70 | . | 59.60 | . | . | . | 2.80 | . | 0.80 | . | . |
| 1 | ACIRS S2D | 2.877 | . | . | . | . | . | . | 0.010 | . | . | . | 0.281 | . | . | . |
| 1 | NCS FC28216 | 2.79 | 8.62 | 10.86 | 32,300 J | 1.44 | . | 81.27 | . | . | . | 3.55 | . | 1.16 | . | . |
| 1 | NCS FC28008m | 2.68 | 25.20 | 27.50 | 24,490 J | 1.57 | . | 60.26 | . | . | . | 3.74 | . | 1.06 | . | . |
| 1 | CZ SF-07-14 | 2.52 | 28.73 | 38.80 | 21,337 J | 9,173 BTU | . | 50.97 | . | . | . | 4.26 | . | 1.05 | . | . |
| 2 | ACIRS M1 | 2.39 | 18.68 | (32.8) | (28,360) J | 1.417 | . | (67.4) | . | 0.014 | 0.0160 | (4.73) | 0.122 | (1.34) | . | 0.098 |
| 2 | COCO 036 | 2.28 | 13.38 | 8.28 | 30,550 J | . | . | . | . | . | . | . | . | . | . | 0.088 |
| 1 | NCS FC28215 | 2.17 | 25.2 | 28.79 | 24,830 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | 502-672-17101 | 2.53 | 8.01 | . | . | . | . | . | . | . | . | . | . | . | . | 17025 |
| 1 | NCS FC28011k | 2.24 | 22.20 | 8.57 | 25,460 J | 1.76 | . | 68.80 | . | . | . | 2.07 | . | 0.86 | . | . |
| 1 | NCS FC28112 | 2.07 | 8.08 | 33.70 | 32,620 J | 1.33 | . | 78.64 | . | . | . | 5.01 | . | 1.31 | . | . |
| 2 | 502-687-13233 | 2.06 | . | . | . | . | . | . | . | . | . | . | 0.130 | . | . | . |
| 2 | 502-839-11060 | 2.03 | 7.72 | (37.57) | (13,586) BTU | . | (2.37) | (76.19) | . | 0.0975 | 0.0071 | 5.10 | 0.0132 | 1.46 | Se: 1.0 | ppm |
| 1 | NCS FC28002x | 1.99 | 16.47 | 30.10 | 27,050 J | 1.50 | . | 66.95 | . | . | . | 4.05 | . | 1.19 | . | . |
| 1 | SRM 2683C | 1.95 | . | . | . | . | . | . | (0.1127) | . | . | . | 0.0900 | . | . | . |
| 1 | NCS FC28138b | 1.92 | 43.85 | 11.81 | 18,090 J | 1.82 | . | 46.46 | . | . | . | 2.54 | . | 0.73 | . | . |
| 1 | NCS FC28007p | 1.92 | 16.26 | 30.12 | 27,060 J | 1.49 | . | 67.05 | . | . | . | 4.06 | . | 1.21 | . | . |
| 2 | COCO 002 | 1.89 | 14.25 | 21.16 | 29,200 J | . | . | . | . | . | . | . | . | last of stock | . | 0.037 |
| 1 | NCS FC28217 | 1.79 | 8.68 | 36.06 | 31,330 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS FC28002z | 1.78 | 15.68 | 30.18 | 27,460 J | 1.48 | . | 67.88 | . | . | . | 4.08 | . | 1.19 | . | . |
| 1 | NCS FC28209 | 1.76 | 27.33 | 8.21 | 23,960 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS FC28144b | 1.75 | 85.72 | 9.12 | 2,280 J | 2.54 | . | 5.80 | . | . | . | 1.03 | . | 0.27 | . | . |
| 1 | SABS 062 | 1.74 | 11.46 | 5.34 | 30,170 J | . | . | 82.28 | . | . | . | 2.16 | . | 1.83 | . | 0.040 |
| 1 | NCS FC28005j | 1.73 | 16.90 | 7.65 | 28,110 J | 1.66 | . | 74.90 | . | . | . | 2.36 | . | 0.99 | . | . |
| 1 | NCS FC28106a | 1.72 | 8.25 | 31.68 | 32,420 J | 1.35 | . | 78.58 | . | . | . | 4.84 | . | 1.35 | . | . |
| 2 | COCO 041 | 1.67 | 29.60 | 17.96 | 22,870 J | . | . | . | . | . | . | . | . | . | . | . |
| 2 | COCO 039 | 1.67 | 14.65 | 17.89 | 30,240 J | . | . | . | . | . | . | . | . | . | . | 0.040 |
| 1 | NCS FC28214 | 1.66 | 27.85 | 29.21 | 23,630 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS FC28010L | 1.66 | 19.64 | 31.33 | 27,010 J | 1.48 | . | 65.68 | . | . | . | 4.24 | . | 1.19 | . | . |
| 2 | 502-831-14095 | 1.57 | 6.99 | (34.30) | (12,765) BTU | FSI:0.78 | 5.08 | (73.08) | . | 0.3545 | 0.0048 | 4.72 | 0.117 | 1.69 | Se:0.00013 | . |
| 2 | COCO 046 | 1.49 | 16.65 | 19.36 | 29,310 J | . | . | . | . | . | . | . | . | . | . | 0.039 |
| 1 | NCS FC28213 | 1.46 | 9.87 | 35.24 | 30,620 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS FC28218 | 1.35 | 14.58 | 6.16 | 29,260 J | . | . | . | . | . | . | . | . | . | . | . |
| 2 | COCO 030 | 1.34 | 25.52 | 10.06 | 24,350 J | . | . | . | . | . | . | . | . | . | . | . |
| 1 | CZ SF-01-14 | 1.33 | 44.90 | 31.72 | 14,617 J | 6,284 BTU | . | 36.40 | . | . | . | 3.31 | . | 0.60 | . | . |
| 1 | 502-681-19066 | 1.32 | 9.18 | 33.2 | (13,780) BTU | . | . | 77.5 (57.6) | . | . | . | 4.99 | . | 1.45 | . | 17025 |
| 2 | COCO 044 | 1.30 | 28.20 | 20.08 | 23,000 J | . | . | . | . | . | . | . | . | . | . | 0.047 |
| 1 | NCS FC28139 | 1.30 | 22.70 | 18.37 | 27,040 J | 1.51 | . | 67.18 | . | . | . | 3.68 | . | 1.05 | . | . |
| 1 | NCS FC28140 | 1.28 | 25.88 | 30.43 | 22,500 J | 1.62 | . | 58.12 | . | . | . | 3.40 | . | 1.02 | . | . |
| 1 | NCS FC28111 | 1.26 | 25.50 | 28.50 | 23,850 J | 1.57 | . | 59.75 | . | . | . | 3.73 | . | 1.01 | . | . |
| 2 | USZ TTKN | 1.21 | 15.2 | 29.3 | 35,260 J | . | . | . | . | . | . | 4.88 | . | 80, 150, or 200 g units | . | . |
| 2 | 502-685-14005 | 1.16 | . | . | . | . | . | . | . | . | . | . | 0.040 | . | . | . |
| 1 | 502-685-16140 | 1.15 | . | . | . | . | . | . | . | . | . | . | 0.041 | . | . | 17025 |
| 2 | COCO 040 | 1.15 | 10.52 | 36.37 | 30,050 J | . | . | . | . | . | . | . | . | . | . | 0.008 |
| 1 | 502-671-18363 | 1.14 | 10.21 | . | . | . | . | . | . | . | . | . | . | . | . | 17025 |
| 1 | 502-686-17294 | 1.09 | . | . | . | . | . | . | . | . | . | . | 0.042 | . | . | 17025 |

| # | Number | S | Ash | Volatile Matter | Heat in J/g or BTU/lb | Density | Moisture | C | Fixed C | Cl | F | H | Hg* | N | O | P |
|---|--------|---|-----|-----------------|-----------------------|---------|----------|---|---------|----|---|---|-----|---|---|---|
|---|--------|---|-----|-----------------|-----------------------|---------|----------|---|---------|----|---|---|-----|---|---|---|

COAL chart 2 of 2

| #=class, 1=CRM and 2=RM | | mass % except * for mg/kg | | | | ACIRS, AS(C)RM, COCO: 120-250 g | | | SABS: 100-150 g | | USZ: as shown | | others: 50 g | | | | |
|-------------------------|---------------|---------------------------|---------|-----------------|-----------------------|---------------------------------|----------|--------|------------------------|-----------|---------------|------------|--------------|--------|---------|-------------------|--------------------|
| # | Number | S | Ash | Volatile Matter | Heat in J/g or BTU/lb | Density | Moisture | C | Fixed C | Cl | F | H | Hg* | N | O | P | |
| 1 | SRM 2692C | 1.064 | (7.499) | | | | | | | | | | 0.1790 | | | | |
| 1 | SABS 069 | 1.06 | 8.22 | 6.10 | 32,060 | | | 85.36 | | (0.1338*) | | 2.66 | | 1.89 | | 0.018 | |
| 1 | NCS FC28202 | 1.05 | 8.65 | 33.23 | 30,720 | 1.39 | | 74.78 | | | | 4.78 | | 1.35 | | | |
| 1 | NCS FC28205 | 1.05 | 8.62 | 33.44 | 30,750 | | | | | | | | | | | | |
| 1 | NCS FC28208 | 1.03 | 15.47 | 20.59 | 29,200 | | | 1.46 | | | | 3.98 | | 1.14 | | | |
| 1 | NCS FC28004k | 1.03 | 12.10 | 8.59 | 30,230 | 1.59 | | 79.82 | | | | 2.61 | | 1.20 | | | |
| 1 | NCS FC28105a | 1.03 | 11.87 | 8.99 | 31,190 | 1.46 | | 79.96 | | | | 3.31 | | 1.12 | | | |
| 1 | SABS 035 | 1.02 | 26.81 | 22.18 | | | | 58.01 | | | | 2.98 | | 1.52 | | 0.05 | |
| 1 | SABS 028 | 0.99 | 27.00 | 23.10 | | | | 57.24 | | | | 2.94 | | 1.45 | last | 0.09 | |
| 2 | COCO 049 | 0.98 | 25.23 | 17.23 | 24,400 | | | | | | | | | | | 0.052 | |
| 1 | SABS 065 | 0.97 | 8.0 | 5.7 | 32,070 | | | 85.60 | | | | 2.68 | | 1.84 | | 0.020 | |
| 2 | COCO 054 | 0.96 | 22.87 | 25.36 | 25,120 | | | 63.57 | | | | 3.56 | | 1.43 | | 0.033 | |
| 1 | SABS 068 | 0.95 | 27.09 | 21.22 | 25,600 | | | 63.47 | | | | 3.57 | | 1.46 | | 0.038 | |
| 2 | COCO 028 | 0.95 | 17.86 | volati | 28,660 | | | | | | | | | | | | |
| 1 | NCS FC28204 | 0.95 | 8.11 | 34.01 | 31,340 | 1.36 | | 76.22 | | | | 4.93 | | 1.45 | | | |
| 2 | 502-838-12185 | 0.963 | 6.50 | (33.78) | (14,337) | BTU | | (1.41) | (80.26) | | 0.2161 | 0.0060 | 5.17 | 0.067 | 1.60 | Se: 2.7 ppm | |
| 2 | COCO 024 | 0.96 | 15.51 | 26.76 | 27,920 | | | | | | | | | | | 0.031 | |
| 2 | COCO 011 | 0.93 | 18.65 | 24.31 | 25,020 | | | | | | | | | | | 0.008 | |
| 2 | COCO 052 | 0.92 | 24.00 | 23.09 | 24,540 | | | | | | | | | | | 0.028 | |
| 2 | COCO 047 | 0.92 | 23.59 | 21.56 | 24,420 | | | | | | | | | | | 0.097 | |
| 1 | SABS 055 | 0.88 | 14.1 | 11.5 | 29,120 | | | 75.70 | | | | 2.93 | | 1.93 | | 0.031 | |
| 1 | NCS FC28211 | 0.88 | 13.41 | 9.08 | 30,230 | | | | | | | | | | | | |
| 1 | SABS 029 | 0.86 | 32.97 | 23.96 | | | | 50.86 | | | | 2.86 | | 1.17 | 0.051 | | |
| 2 | COCO 042 | 0.85 | 40.79 | 21.25 | 16,700 | | | 59.32 | | | | 2.98 | | 1.56 | | | |
| 1 | NCS FC28206 | 0.85 | 14.46 | 28.56 | 26,720 | | | | | | | | | | | | |
| 2 | COCO 014 | 0.83 | 21.88 | 8.15 | 27,550 | | | 71.25 | | | | 2.73 | | 1.46 | | 0.016 | |
| 1 | SABS 059 | 0.82 | 16.0 | 11.0 | 28,330 | | | 74.50 | | | | 2.84 | | 1.79 | | 0.031 | |
| 1 | NCS FC28006f | 0.82 | 9.08 | 33.09 | 30,580 | | | 74.96 | | 1.40 | | 4.59 | | 1.37 | | | |
| 1 | NCS FC28110a | 0.81 | 9.62 | 33.64 | 29,830 | | | 74.16 | | 1.42 | | 4.44 | | 1.38 | | | |
| 2 | USZ SOEN | 0.78 | 6.29 | 44.8 | 32,160 | | | 76.62 | | | | 5.41 | | 1.44 | | 80 or 200 g units | |
| 1 | SABS 039 | 0.75 | 24.51 | 23.34 | | | | 59.97 | | | | 2.93 | | 1.56 | | 0.079 | |
| 1 | SABS 043 | 0.74 | 22.31 | 23.82 | | | | 61.69 | | | | 3.03 | | 1.56 | | 0.071 | |
| 2 | COCO 034 | 0.74 | 17.54 | 21.45 | 27,590 | | | | | | | | | | | | |
| 2 | COCO 038 | 0.73 | 15.02 | 23.63 | 28,350 | | | | | | | | | | | 0.090 | |
| 1 | SABS 051 | 0.72 | 39.70 | 20.89 | | | | 44.44 | | | | 2.45 | | 1.10 | | 0.11 | |
| 1 | NCS FC28203 | 0.71 | 10.36 | 20.69 | 31,660 | | | | | | | | | | | | |
| 2 | 502-846-13156 | 0.703 | 12.19 | (5.81) | (13,197) | BTU | | (2.14) | (81.94) | | 0.0202 | 0.0086 | 2.39 | 0.0112 | 1.05 | Se: 2.5 ppm | |
| 1 | NCS FC28107a | 0.70 | 10.96 | 15.62 | 31,120 | 1.45 | | 79.17 | | | | 3.82 | | 1.14 | | | |
| 2 | ACIRS S2B | 0.692 | (9.1) | (34.0) | (31,050) | J | | (75.4) | | | 0.027 | 0.0037 | (4.90) | 0.057 | (1.79) | Se: (0.4 ppm) | |
| 1 | 502-682-19115 | 0.676 | 7.36 | 40.8 | (12,369) | BTU | | 71.4 | (51.8) | | | 4.76 | | 1.44 | | 17025 | |
| 2 | 502-845-13158 | 0.665 | 8.02 | (3.83) | (13,573) | BTU | | (3.15) | (87.33) | | 0.0154 | 0.0051 | 1.76 | 0.0169 | 0.878 | Se: 2.0 ppm | |
| 1 | SABS 046 | 0.66 | 11.86 | 26.87 | | | | 74.21 | | | | 3.77 | | 1.76 | | | |
| 1 | SABS 026 | 0.65 | 37.83 | 22.07 | | | | 46.63 | | | | 2.59 | | 1.11 | 0.066 | | |
| 1 | NM 150 | 0.64 | 35.88 | 22.02 | 4,967 | kcal/kg | | | | | | | | | | | |
| 1 | SABS 047 | 0.60 | 13.58 | 25.45 | | | | 71.85 | | | | 3.81 | | 1.66 | | 0.06 | |
| 2 | COCO 025 | 0.59 | 23.93 | 22.60 | 24,300 | | | | | | | | | | | 0.036 | |
| 2 | ACIRS G8 | 0.590 | 10.14 | 20.23 | 32,357 | J | | 1.382 | 79.44 | | 0.052 | Se:0.00008 | 4.40 | 0.02 | 1.79 | 0.027 | |
| 2 | ACIRS G6 | 0.59 | 10.72 | 19.93 | 32,150 | J | | 1.372 | CO ₂ :0.055 | 79.36 | last | 0.041 | 0.0067 | 4.33 | 0.026 | 1.75 | Se:0.8ppm |
| 2 | COCO 045 | 0.59 | 10.70 | 27.61 | 29,630 | J | | | | | | | | | | 0.035 | |
| 1 | NCS FC28109 | 0.58 | 11.98 | 11.42 | 30,560 | J | | 1.49 | 79.02 | | | 3.28 | | 1.03 | | | |
| 1 | SABS 066 | 0.57 | 15.02 | 24.11 | 27,730 | | | | 71.97 | | | 3.62 | | 1.71 | | 0.110 | |
| 1 | NCS FC28108 | 0.57 | 13.68 | 30.55 | 29,530 | J | | 1.42 | 72.65 | | | 4.46 | | 1.23 | | | |
| 2 | ACIRS G5 | 0.568 | 10.68 | 20.03 | 32,186 | J | | 1.381 | (79.19) | | 0.046 | (4.34) | | (1.77) | | 0.0196 | |
| 2 | COCO 055 | 0.54 | 18.96 | 23.85 | 25,820 | J | | | 66.37 | | | 3.54 | | 1.57 | | 0.100 | |
| 2 | USZ BNN | 0.54 | 9.04 | 46.28 | 28,950 | J | | | 70.49 | | | 5.27 | | | | 150 g units | |
| 2 | COCO 043 | 0.54 | 14.80 | 23.94 | 28,490 | J | | | | | | | | | | 0.107 | |
| 1 | 502-670-17354 | 0.53 | 8.30 | | | | | | | | | | | | | 17025 | |
| 1 | NCS FC28212a | 0.51 | 9.36 | 23.27 | 32,380 | J | | | | | | | | | | | |
| 1 | SABS 058 | 0.50 | 31.1 | 23.1 | 20,180 | J | | | 53.64 | | | 3.06 | | 1.33 | | 0.023 | |
| 1 | NCS FC28116a | 0.50 | 11.06 | 23.77 | 30,000 | J | | 1.44 | 75.34 | | | 4.00 | | 1.18 | | | |
| 1 | SABS CCS 008b | 0.49 | 15.8 | 25.6 | 27,550 | | | | 71.19 | | | 3.85 | | 1.70 | | 0.111 | |
| 1 | SABS 050 | 0.48 | 16.36 | 24.58 | | | | 70.79 | | | | 3.80 | | 1.62 | | 0.095 | |
| 1 | SABS 049 | 0.48 | 16.34 | 24.31 | | | | 70.96 | | | | 3.61 | | 1.63 | last | 0.087 | |
| 1 | SABS 057 | 0.48 | 15.49 | 25.42 | | | | 71.19 | | | | 3.85 | | 1.70 | | 1.093 | |
| 1 | SABS 045 | 0.49 | 15.45 | 24.50 | | | | 71.24 | | | | 3.73 | | 1.71 | | | |
| 1 | SABS 037 | 0.48 | 15.26 | 24.84 | | | | 71.17 | | | | 3.67 | | 1.72 | 0.102 | | |
| 1 | CZ SF-04-14 | 0.48 | 4.43 | 23.67 | 34,618 | J | 14,883 | BTU | 85.53 | | | 4.59 | | 1.35 | | | |
| 2 | COCO 053 | 0.47 | 14.37 | 27.71 | 27,980 | J | | | 71.03 | | | 4.00 | | 1.68 | | 0.129 | |
| 2 | 502-835-10060 | 0.469 | 11.08 | (38.20) | (12,604) | BTU | | | (70.05) | | | 4.85 | 0.015 | 1.75 | | | |
| 1 | SRM 2693 | 0.4571 | | | | | | | | | 0.03696 | | | | | | |
| 2 | COCO 050 | 0.45 | 13.19 | 27.10 | 28,350 | J | | | | | | | 0.0373 | | | 0.141 | |
| 1 | NCS FC28001i | 0.45 | 10.75 | 19.46 | 32,080 | J | | 1.39 | 79.24 | | | 4.20 | | 1.44 | | | |
| 2 | ACIRS S2A | 0.433 | | | | | | | | | | | | | | | |
| 1 | NCS FC28207 | 0.43 | 16.26 | 7.26 | 26,100 | J | | | | | | | | | | | |
| 1 | NCS FC28115 | 0.42 | 6.38 | 31.74 | 30,570 | J | | 1.41 | 77.28 | | | 4.47 | | 1.19 | | | |
| 1 | NCS FC28104 | 0.40 | 10.09 | 11.00 | 31,860 | J | | 1.45 | 81.45 | | | 3.52 | | 1.31 | | | |
| 1 | NCS FC28201a | 0.40 | 9.70 | 18.48 | 31,580 | J | | 1.42 | 78.94 | | | 3.98 | | 1.36 | | | |
| 2 | COCO 001 | 0.39 | 14.66 | 24.58 | 26,820 | J | | | | | | | | | last | 0.079 | |
| 2 | 502-833-14255 | 0.383 | 8.96 | (18.20) | (14,175) | BTU | | | (0.862) | (82.03) | | 0.0241 | 0.0098 | (4.19) | 0.00595 | 1.16 | Se:1.22* FSI:(3.6) |
| 1 | 502-680-18277 | 0.38 | 8.35 | 18.6 | (14,209) | BTU | | | 82.3 | (73.1) | | 4.21 | | 1.21 | | 17025 | |
| 1 | SABS 023 | 0.37 | 16.5 | 25.6 | | | | | 68.25 | | | 3.63 | | 1.62 | | 0.102 | |
| 1 | NCS FC28001h | 0.37 | 8.88 | 31.87 | 31,010 | J | | 1.39 | 76.30 | | | 4.53 | | 1.44 | | | |
| 1 | NCS FC28103 | 0.35 | 10.51 | 9.70 | 31,740 | J | | 1.47 | 81.44 | | | 3.40 | | 1.26 | | | |
| 1 | NCS FC28003m | 0.33 | 10.10 | 10.33 | 32,030 | J | | 1.45 | 8 | | | | | | | | |

RM FUSIBILITY OF COAL

| analysis listed in mass % | | | | | | | | | 250 g units | | Reducing Temperature °C | | | |
|---------------------------|-------|------|------|-------|------|--------|------------|-------|-------------|-----------|-------------------------|------|--|--|
| Number | C | H | N | P | S | MJ/Kg | Volatility | Ash | Deformation | Softening | Hemisphere | Flow | | |
| COCO 005 | 81.70 | 4.57 | 1.44 | 0.015 | 1.05 | 32.90 | 27.19 | 7.49 | 1402 | 1425 | 1443 | 1473 | | |
| COCO 007 | 72.55 | 3.69 | 1.83 | 0.036 | 1.55 | 28.71 | 21.60 | 14.51 | 1329 | 1353 | 1381 | 1420 | | |
| COCO 016 | 56.64 | 2.64 | 1.53 | 0.030 | 1.86 | 21.92 | 18.26 | 31.32 | 1284 | 1317 | 1346 | 1387 | | |
| COCO 035 | 62.17 | 2.98 | 1.71 | 0.039 | 1.92 | 24.16 | 15.46 | 27.55 | 1369 | 1394 | 1419 | 1452 | | |
| COCO 051 | 55.80 | 2.91 | 1.30 | 0.075 | 1.12 | 21,520 | 21.73 | 30.98 | 1359 | 1373 | 1390 | 1448 | | |

SULFUR IN COAL

= class, where 1=CRM and 2=RM

analysis listed in mass %

| # | Number | S | Units | # | Number | S | Units | # | Number | S | Units |
|---|---------------|------|------------|---|-----------|-------|-----------|---|---------------|-------|-----------|
| 1 | IARM HC20800A | 8.6 | 50 g | 1 | BCR 336 | 3.290 | 20 g | 2 | VS1-0.96 | 0.96 | 50g, last |
| 1 | ASCRM 012 D | 5.21 | 125g, last | 2 | VS1-1.91 | 1.91 | 50g, last | 1 | IARM HC20075C | 0.76 | 50g, last |
| 2 | COCO 037 | 4.74 | 50 g | 2 | ACIRS S1C | 1.46 | 250 g | 1 | BCR 331 | 0.499 | 20 g |
| 2 | VS1-4.18 | 4.18 | 50g, last | 1 | BCR 333 | 1.344 | 20 g | 2 | ACIRS S1A | 0.42 | 250 g |
| | | | | 1 | BCR 332 | 0.961 | 20 g | | | | |

RM COAL

typical analysis listed in mass %

50 g units, last of stock

| Number | S | DRY ANALYSIS | | | | IGNITED ANALYSIS | | | | | | | | | | |
|---------|------|--------------|-------------|-------|-----------------|--------------------------------|------|--------------------------------|------------------|------|------------------|-------------------|-------------------------------|-----------------|------------------|------------------|
| | | C | Heat BTU/lb | Ash | Volatile Matter | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO ₂ | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | TiO ₂ |
| VS6-016 | 1.41 | (47.64) | (12,293) | 16.71 | 35.59 | 27.07 | 0.81 | 7.96 | 3.56 | 1.18 | 0.02 | 0.38 | 0.12 | 0.77 | 55.62 | 1.20 |

CRM COAL

analysis listed in mass %

(T) = Total

SARM 20 also contains Ta: 0.00012, Y: 0.0029

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Mn | Na ₂ O | P | P ₂ O ₅ | S | SiO ₂ | Sr | TiO ₂ | Zr | LOI | Units |
|----------|--------------------------------|------|--------------------------------|------------------|--------|----------|-------------------|--------|-------------------------------|--------|------------------|--------|------------------|----------|-------|-------|
| SARM 20 | 11.27 | 1.87 | 1.17 | 0.14 | 0.43 | 0.0080 | 0.27 | . | 0.14 | 0.51 | 17.66 | 0.0330 | 0.63 | (0.0180) | 64.66 | 120 g |
| SARM 19 | 8.01 | 1.39 | 1.75 | 0.24 | 0.20 | 0.0157 | 0.29 | 0.0130 | . | 1.49 | 15.00 | 0.0126 | 0.341 | 0.0351 | 71.28 | 120 g |
| SARM 18 | 2.57 | 0.18 | 0.29 | 0.145 | 0.11 | 0.0022 | . | 0.0030 | . | 0.56 | 6.20 | 0.0044 | 0.114 | 0.0067 | 90.11 | 120 g |
| US CLB-1 | (1.51) | 0.22 | 1.25(T) | 0.0760 | 0.0470 | (0.0008) | 0.0230 | . | (0.0700) | (1.49) | (2.51) | . | (0.0780) | . | . | 40 g |

analysis listed in mg/kg US CLB-1 also contains informational Ash, Li, Mo, Nb, Nd and Sb

| Number | As | Ba | Be | Ce | Co | Cr | Cs | Cu | Ga | Ge | Hf | Hg | La | Ni | Pb | Rb | Sc | Se | Sm | Th | U | V | Zn |
|----------|------|-----|-----|----|-----|------|-----|------|-----|-----|-----|--------|-----|------|-----|-----|-----|-----|-----|-------|--------|----|-----|
| SARM 20 | 4.7 | 372 | 2.5 | 87 | 8.3 | (67) | (2) | 18 | 16 | . | 4.8 | 0.25 | 43 | 25 | 26 | 10 | 10 | 0.8 | 6.3 | 18 | 4 | 47 | 17 |
| SARM 19 | 7 | 304 | 2.8 | 56 | 5.6 | 50 | 1.4 | 13 | 14 | 13 | 5.4 | (0.2) | 27 | 16 | 20 | 9 | 7.6 | . | 4.9 | 12 | 5 | 35 | 12 |
| SARM 18 | . | 78 | 4.1 | 22 | 6.7 | 16 | (1) | 5.9 | (8) | (8) | 1.7 | (0.04) | 10 | 10.8 | (5) | 8.1 | 4.3 | . | 2.0 | 3.4 | 1.5 | 23 | 5.5 |
| US CLB-1 | (13) | 34 | . | 10 | 7.0 | 9.7 | . | (10) | (3) | . | . | (0.2) | (5) | 18 | 5.1 | 5.2 | 2.0 | (2) | . | (1.4) | (0.55) | 12 | 48 |

CRM COAL

analysis in mass %

50g units

analysis in mg/kg

| Number | C | Al | Cl | Fe | H | K | N | Na | S | Hg | Mn | V | Zn |
|-------------|---------|----------|----------|--------|--------|----------|---------|----------|---------|--------|---------|--------|------|
| SRM 2684c | (76.82) | (0.8730) | (0.0975) | . | (5.17) | (0.0981) | (1.395) | (0.0606) | 3.027 | 0.0688 | (20.51) | (16.3) | . |
| SRM 1635a * | (68.97) | 0.5437 | (0.0051) | 0.2472 | 3.92 | 0.01874 | (0.946) | 0.1031 | (0.294) | 0.0836 | 6.69 | 13.34 | 7.3 |
| SRM 1632e * | (76) | 0.960 | (0.0963) | 1.42 | (4.97) | 0.1248 | (1.4) | 0.0374 | 2.738 | 0.1351 | (18.4) | (29.2) | 13.0 |

analysis in mg/kg * SRM 1632d, 1632e and 1635a also detail many other elements, see certificate

| Number | As | Ba | Br | Cd | Ce | Co | Cr | Cu | F | Ni | Pb | Rb | Sb | Se | Sr | Th | Ti | U |
|-------------|--------|--------|--------|-------|-----------|---------|---------|--------|------|----------|------|-------|---------|---------|------|-------|-----|---------|
| SRM 2684c | . | . | (11.1) | . | Ca:(3220) | . | . | . | (64) | Mg:(494) | . | . | . | (1.08) | . | . | . | . |
| SRM 1635a * | 0.860 | 357.8 | (1) | 0.282 | 5.45 | 2.004 | 3.56 | 11.42 | (63) | 5.37 | 2.85 | 1.226 | 0.251 | 0.662 | 160 | 1.299 | 254 | 0.4792 |
| SRM 1632e * | (8.55) | (62.8) | (11.9) | . | (12.24) | (3.622) | (16.57) | (5.70) | . | (11.08) | . | 8.49 | (0.428) | (1.525) | 84.1 | . | 519 | (0.636) |

CRM COAL

analysis listed in mass % except * which is mg/kg

powder 50 g

| Number | Al% | Ca% | Cd* | Co* | Cr* | Cu* | Fe% | K% | Mg% | Mn% | Na% | Ni* | P% | Pb% | Si% | Ti% | V* | Zn% |
|-------------|------|------|------|-----|-----|-----|------|-------|-------|--------|-------|-----|--------|--------|------|-------|----|----------|
| NCS FC28127 | 3.47 | 1.88 | 2 | 9 | 23 | 23 | 1.02 | 0.29 | 0.28 | 0.019 | 0.052 | 16 | 0.010 | . | 5.61 | 0.18 | 60 | 0.0040 |
| NCS FC28125 | 2.27 | 0.28 | (<1) | 11 | 5 | 17 | 0.24 | 0.090 | 0.050 | 0.0009 | 0.048 | 18 | 0.013 | 0.0016 | 2.69 | 0.090 | 33 | . |
| NCS FC28123 | 1.88 | 0.74 | (<1) | 4 | 10 | 12 | 0.35 | 0.026 | 0.081 | 0.0030 | 0.11 | 8 | 0.066 | 0.0016 | 1.86 | 0.096 | 12 | (0.001) |
| NCS FC28124 | 1.75 | 0.79 | (<1) | 4 | 7 | 12 | 0.34 | 0.020 | 0.071 | 0.0016 | 0.13 | 8 | 0.044 | 0.0016 | 1.77 | 0.079 | 11 | . |
| NCS FC28128 | 1.22 | 0.19 | . | 4 | 8 | 12 | 0.86 | 0.043 | 0.059 | 0.0026 | 0.026 | 8 | 0.0044 | . | 1.64 | 0.059 | 28 | (<0.001) |
| NCS FC28126 | 0.83 | 0.65 | (<1) | 3 | 5 | 8 | 0.32 | 0.010 | 0.060 | 0.008 | 0.034 | 5 | 0.019 | . | 1.01 | 0.046 | 11 | . |
| NCS FC28122 | 0.25 | 0.85 | . | 8 | 2 | 2 | 1.79 | 0.016 | 0.24 | 0.022 | 0.081 | 8 | 0.0029 | 0.002 | 0.47 | 0.010 | 1 | . |

CRM COAL

analysis listed in mass %

| Number | Al ₂ O ₃ | C | Fe ₂ O ₃ | H | N | P ₂ O ₅ | S | SO ₃ | SiO ₂ | TiO ₂ | Units |
|--------|--------------------------------|-------|--------------------------------|------|------|-------------------------------|------|-----------------|------------------|------------------|-------|
| NM 151 | 26.09 | 52.71 | 6.86 | 3.30 | 1.13 | 0.91 | 0.64 | 0.32 | 59.45 | 1.70 | 50 g |

CRM COAL

BCR: 40 g units GBW: 50 g units

| Number | As mg/kg | P mass % | Cl mass % | F mg/kg |
|-----------|-------------|-------------|--------------|------------|
| GBW 11115 | 15 | 0.031 | . | . |
| GBW 11116 | 34 | 0.007 | . | . |
| GBW 11117 | 51 | 0.092 | . | . |
| GBW 11118 | . | . | 0.010 | . |
| GBW 11119 | . | . | 0.057 | . |
| GBW 11120 | . | . | 0.110 | . |
| GBW 11121 | . | . | . | 248 |
| GBW 11122 | . | . | . | 864 |
| GBW 11123 | . | . | . | 1496 |
| BCR 460 | . | . | (0.0059) | 225 |

CRM COAL AIR DRIED vs. HEATED DRIED ANALYSIS

20 g powder

| Number | Heat J/g | Volatile Matter% | Ash% | Moisture% | S% | Expiry |
|--------------|--------------|------------------|--------------|-------------------|------------|---------------------|
| NCS FC62002a | 24190, 24980 | 4.29, 4.43 | 23.90, 24.68 | 3.15 (air dried) | 0.36, 0.37 | August 2019 |
| NCS FC62001c | 22840, 25680 | 28.37, 31.90 | 11.30, 12.71 | 11.06 (air dried) | 0.63, 0.71 | May 2018 H: (3.32%) |

CRM FUSIBILITY OF COAL ASH

analysis listed in °C MRed = Mildly Reducing, Oxi = Oxidizing, SRed = Strongly Reducing

| Atmosphere Number | Initial Deformation | | | Softening | | | Hemishpering | | | Fluid | | | Units |
|----------------------|---------------------|-----|------|------------------|-----|------|------------------|-----|------|------------------|-----|------|-------|
| | MRed | Oxi | SRed | MRed | Oxi | SRed | MRed | Oxi | SRed | MRed | Oxi | SRed | |
| NCS FS91001d | 1057, 1176, 1208 | | | 1072, 1202, 1253 | | | 1098, 1236, 1328 | | | 1148, 1320, 1401 | | | 30 g |
| NCS FS28001 | 1161, 1211, . | | | 1190, 1230, . | | | 1198, 1239, . | | | 1204, 1252, . | | | 5 g |
| NCS FS28002 | 1217, 1356, . | | | 1340, 1408, . | | | 1357, 1420, . | | | 1369, 1445, . | | | 5 g |
| NCS FS28003 | 1285, 1314, . | | | 1314, 1345, . | | | 1322, 1360, . | | | 1340, 1381, . | | | 5 g |

COAL ASH

= class, where 1 = CRM and 2 = RM

analysis listed in mass %

| # | Number | SiO ₂ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Mn | Na ₂ O | P ₂ O ₅ | SO ₃ | TiO ₂ | V ₂ O ₅ | Units | Other |
|---|---------------|------------------|--------------------------------|-------|--------------------------------|------------------|------|--------------------------|-------------------|-------------------------------|-----------------|------------------|-------------------------------|-------|-------|
| 2 | COCO ASH 015 | 55.67 | 31.44 | 2.42 | 3.99 | 0.65 | 0.74 | MnO ₂ : 0.036 | 0.13 | 0.37 | 1.89 | 1.78 | . | 20 g | |
| 2 | COCO ASH 013 | 54.89 | 32.44 | 2.16 | 4.43 | 0.84 | 0.54 | MnO ₂ : 0.041 | 0.06 | 0.33 | 2.05 | 1.82 | . | 20 g | |
| 1 | NCS FC28154 | 53.17 | 32.02 | 2.28 | 6.47 | 1.37 | 0.90 | MnO: 0.035 | 0.41 | 0.19 | 0.78 | 1.34 | 0.027 | 5 g | |
| 2 | COCO ASH 014 | 52.29 | 31.67 | 2.53 | 6.46 | 0.85 | 0.51 | MnO ₂ : 0.050 | 0.07 | 0.44 | 2.73 | 1.86 | . | 20 g | |
| 2 | COCOASHSRM017 | 51.24 | 33.96 | 2.90 | 4.34 | 0.55 | 0.41 | MnO ₂ : 0.044 | 0.11 | 2.31 | 0.99 | 1.70 | . | 20 g | |
| 1 | JCFA-1 * | 50.56 | 24.25 | 8.91 | 4.22 | 1.27 | 2.12 | MnO: 0.068 | 2.24 | 0.586 | . | 1.31 | . | 100 g | |
| 2 | COCO ASH 016 | 50.18 | 29.05 | 6.35 | 5.82 | 0.93 | 1.39 | MnO ₂ : 0.061 | 0.21 | 0.57 | 3.23 | 1.56 | . | 20 g | |
| 1 | NCS FC28148 | 48.03 | 35.80 | 3.27 | 2.81 | 0.90 | 0.69 | MnO: 0.0073 | 0.54 | 0.25 | . | 1.25 | 0.049 | 5 g | |
| 1 | NCS FC28150 | 47.64 | 26.03 | 10.44 | 5.79 | 1.41 | 1.87 | MnO: 0.097 | 0.28 | 0.091 | . | 1.21 | 0.042 | 5 g | |
| 1 | NCS FC28151 | 43.42 | 28.53 | 3.33 | 15.18 | 0.64 | 1.21 | MnO: 0.042 | 0.43 | 0.12 | . | 1.22 | 0.062 | 5 g | |
| 1 | NCS FC28146 | 37.86 | 33.71 | 9.90 | 4.74 | 0.30 | 1.27 | MnO: 0.037 | 1.45 | 1.44 | . | 1.52 | 0.020 | 5 g | |
| 1 | NCS FC28147 | 37.52 | 32.78 | 10.97 | 4.81 | 0.24 | 1.17 | MnO: 0.020 | 1.75 | 1.00 | . | 1.31 | 0.019 | 5 g | |
| 1 | NCS FC28149 | 35.54 | 25.92 | 14.92 | 7.56 | 0.20 | 1.63 | MnO: 0.17 | 0.75 | 0.72 | . | 1.26 | 0.032 | 5 g | |
| 1 | NCS FC28145 | 15.66 | 7.34 | 18.37 | 39.61 | 0.30 | 6.05 | MnO: 0.44 | 1.69 | 0.10 | . | 0.26 | 0.0042 | 5 g | |

* JCFA-1 also contains (in mg/kg) Be: 4.06, Co: 37.4, Cr: 75, Cs: 8.6, Cu: 122, Li: 91, Ni: 32.2, Pb: 47.2, Rb: 54.1, S: 1960, Sb: 2.1, V: 243, and Zn: 63.
 * JCFA-1 also contains (in mass %): FeO: 0.88, TFe₂O₃: 5.2, C: 1.35, H₂O: 0.18, H₂O+: 0.37, Sr: 0.110

CRM COAL ASH

analysis listed in mass %

100 g units

| Number | Al ₂ O ₃ | Ba | CaO | CO ₂ | Org.C | FeO | T.Fe ₂ O ₃ | H ₂ O | H ₂ O+ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | S | SO ₃ | SiO ₂ | TiO ₂ | LOI |
|--------------|--------------------------------|-------|-------|-----------------|--------|-------|----------------------------------|------------------|-------------------|------------------|------|-------|-------------------|-------------------------------|------|-----------------|------------------|------------------|---------|
| VS 7177-95 | 27.07 | 0.028 | 4.88 | . | . | 1.59 | 5.48 | (0.13) | . | 0.59 | 1.48 | 0.059 | 0.14 | (0.064) | . | (0.018) | 58.68 | 0.60 | (0.56) |
| VS 9237-2008 | 9.7 | 0.86 | 43.8 | (2.9) | . | (0.3) | 5.1 | . | . | 0.36 | 5.5 | 0.22 | 0.67 | (0.024) | 3.48 | . | 15.6 | 0.59 | 8.5 |
| VS 7125-94 | 6.79 | 0.225 | 20.91 | 13.20 | (1.33) | 0.40 | 6.28 | (2.41) | (5.74) | 0.51 | 6.70 | 0.094 | 0.22 | 0.059 | 0.17 | (0.40) | 35.80 | 0.35 | (21.29) |

continued analysis listed in mg/kg except % which is mass %

| Number | Ag | As | Au | B% | Be | Bi | Ce | Co | Cr | Cs | Cu | Dy | Er | Eu | F% | Ga | Gd | Ge | Hf | Hg | Ho | La |
|--------------|--------|-------|---------|---------|-----|-------|-----|----|----|-------|-----|-------|-------|-------|----------|------|-------|-------|------|--------|--------|----|
| VS 7177-95 | (0.15) | . | . | (0.044) | 11 | (2.1) | 138 | 25 | 99 | (3.3) | 176 | . | . | (2.6) | (0.017) | (21) | . | (2.7) | (13) | . | . | 70 |
| VS 9237-2008 | (0.2) | (8.0) | . | . | 2.9 | . | 37 | 26 | 42 | (1.2) | 51 | (2.7) | (1.4) | (0.8) | . | 15 | (3.3) | (6) | 3.3 | (0.1) | (14) | 20 |
| VS 7125-94 | 0.16 | (1.0) | (0.003) | 0.0097 | 2.9 | . | 38 | 16 | 45 | (1.1) | 45 | (3.8) | (2.4) | 0.9 | (0.0230) | 9 | (4.2) | 2.3 | 2.6 | (0.03) | (0.87) | 20 |

| Number | Li | Lu | Mo | Nb | Nd | Ni | Pb | Pr | Rb | Sb | Sc | Sm | Sr% | Ta | Tb | Th | Tl | Tm | U | V | W | Y | Yb | Zn | Zr% | |
|--------------|-----|------|-------|-----|----|----|----|-------|----|-----|-----|------|-------|--------|--------|------|------|-------|--------|-----|-------|-------|-----|-----|-------|--------|
| VS 7177-95 | 96 | . | 7.4 | 34 | . | 66 | 35 | . | 22 | . | 27 | (15) | 11 | 0.0403 | . | . | (45) | . | (15) | 145 | (3.7) | 87 | 7.8 | 77 | 0.033 | |
| VS 9237-2008 | (9) | 0.26 | (3.0) | 7.8 | 17 | 68 | 22 | (4.2) | 13 | 3.0 | 8.6 | 3.2 | (3.2) | 0.83 | (0.61) | 0.45 | 7.0 | (0.4) | (0.25) | 3.1 | 63 | . | 15 | 1.5 | 76 | 0.013 |
| VS 7125-94 | 32 | 0.40 | 1.4 | 8.4 | 20 | 49 | 13 | (4.3) | 15 | . | 11 | 4.1 | 2.7 | 0.33 | (0.53) | 0.68 | 5.8 | . | (0.38) | 3.3 | 61 | (1.1) | 29 | 2.6 | 65 | 0.0119 |

CRM COAL ASH

| Number | Ash% | C% | S% | Units |
|--------------|-------|------|-------|-------|
| CZ SFA-01-14 | 96.60 | 3.10 | 0.029 | 50 g |

CRM COAL WASTE ROCK analysis listed in mass % 50 g units

| Number | Al | Ca | Fe | K | Mg | Mn | Na | P | Si | Ti | V |
|-------------|-------|------|------|------|------|-------|------|-------|-------|------|-------|
| NCS FC28152 | 10.76 | 0.34 | 2.57 | 1.27 | 0.53 | 0.023 | 0.15 | 0.026 | 20.59 | 0.44 | 0.012 |

CRM ASH OF COAL WASTE ROCK analysis in mass % 5 g units

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SiO ₂ | TiO ₂ | V ₂ O ₅ |
|-------------|--------------------------------|------|--------------------------------|------------------|------|-------|-------------------|-------------------------------|------------------|------------------|-------------------------------|
| NCS FC28153 | 27.71 | 0.65 | 5.01 | 2.09 | 1.20 | 0.041 | 0.27 | 0.082 | 60.03 | 1.01 | 0.028 |

CRM COAL FLY ASH analysis listed in mass %

| Number | As | Al | Ba | Ca | Fe | K | Mg | Mn | Na | Ni | P | S | Si | Ti | Zn | LOI |
|----------------|----------|-------|----------|-------|--------|-------|-------|----------|----------|----------|---------|---------|---------|-------|----------|--------|
| SRM 2689 | (0.0200) | 12.94 | (0.0800) | 2.18 | 9.32 | 2.20 | 0.61 | (0.0300) | 0.25 | (0.0122) | 0.10 | . | 24.06 | 0.75 | (0.0240) | (1.76) |
| SRM 1633c | 0.01862 | 13.28 | 0.1126 | 1.365 | 10.49 | 1.773 | 0.498 | 0.02402 | 0.1707 | 0.0132 | (0.192) | (0.110) | (21.30) | 0.724 | (0.0235) | . |
| IRNT 12-1-02 * | 0.00739 | 16.1 | 0.1090 | 1.49 | 5.17 | 0.651 | 0.581 | 0.0441 | 0.361 | 0.0105 | . | . | 22.9 | 3.61 | 0.0221 | . |
| BCR 176R | 0.0054 | . | (0.4650) | . | 1.3100 | . | . | (0.0730) | (3.4800) | 0.0117 | . | . | . | . | 1.6800 | . |
| IRNT 12-1-03 * | 0.00505 | 14.7 | 0.0731 | 1.86 | 5.59 | 1.36 | 0.788 | 0.0383 | 0.303 | 0.0117 | . | . | 25.0 | 1.36 | 0.0248 | . |
| BCR 038 | 0.00480 | . | . | . | 3.3800 | . | . | 0.0479 | 3.740 | (0.0194) | . | . | . | . | 0.0581 | . |
| SRM 2691 | (0.0030) | 9.81 | (0.5900) | 18.45 | 4.42 | 0.34 | 3.12 | (0.0200) | 1.09 | (0.0053) | 0.51 | 0.83 | 16.83 | 0.90 | (0.0120) | (0.23) |
| SRM 2690 | (0.0026) | 12.35 | (0.5800) | 5.71 | 3.57 | 1.04 | 1.53 | (0.0300) | 0.24 | (0.0046) | 0.52 | 0.15 | 25.85 | 0.52 | (0.0120) | (0.53) |

continued analysis listed in mg/kg except % which is mass %

| Number | Ag | Au | B | Be | Br | Cd | Ce | Co | Cr | Cs | Cu | Eu | F | Hf | Hg | La |
|----------------|--------|---------|------|----|-------|-------|--------|-------|--------|--------|--------|---------|------|----------|--------|--------|
| SRM 2689 | . | . | (21) | . | (3) | . | (48) | (170) | (11) | . | (3) | . | (7) | (<0.003) | . | . |
| SRM 1633c | . | . | (16) | . | 0.758 | (180) | 42.9 | (258) | (9.39) | 173.7 | (4.67) | . | . | 1.005 | (87.0) | . |
| IRNT 12-1-02 * | . | . | . | . | . | 322 | 52.4 | 189 | 20.1 | 231 | 4.99 | . | . | . | 164 | . |
| BCR 176R | (33.1) | (0.604) | . | . | (836) | 226 | (47.7) | 26.7 | 810 | (8.27) | 1050 | (0.868) | . | (4.85) | (1.60) | (30.2) |
| IRNT 12-1-03 * | . | . | . | . | . | 183 | 48.6 | 191 | . | 155 | . | . | . | . | 84.4 | . |
| BCR 038 | . | . | . | . | 4.6 | . | 53.8 | (178) | . | 176 | . | . | . | . | . | . |
| SRM 2691 | . | . | (8) | . | (0.9) | . | (26) | (68) | (1) | . | (2) | . | (10) | (<0.003) | . | . |
| SRM 2690 | . | . | (8) | . | (0.7) | . | (19) | (67) | (8) | . | (2) | . | (8) | (<0.003) | . | . |

| Number | Pb | Rb | Sb | Sc | Se | Sr | Ta | Th | Tl | U | V | W | Yb | Units |
|----------------|------|--------|------|--------|--------|--------|--------|--------|------|--------|-------|--------|-------|----------|
| SRM 2689 | (52) | . | (9) | (32) | (7) | (700) | . | (25) | . | . | . | . | . | 3 x 10 g |
| SRM 1633c | 95.2 | 117.42 | 8.56 | (37.6) | (13.9) | 901 | (1.58) | (23.0) | . | (9.25) | 286.2 | . | (7.7) | 75 g |
| IRNT 12-1-02 * | 40.8 | 69.0 | . | . | . | 583 | . | 23.9 | . | . | 558 | . | . | 50 g |
| BCR 176R | 5000 | (102) | 850 | (2.91) | 18.3 | . | (2.02) | (5.28) | 1.32 | . | (35) | (28.3) | . | 40 g |
| IRNT 12-1-03 * | 66.9 | 141 | . | 29.2 | . | 407 | . | 22.1 | . | . | 381 | . | . | 50 g |
| BCR 038 | 262 | . | . | . | . | . | . | . | . | . | . | . | . | 5 to 6 g |
| SRM 2691 | (29) | . | (3) | (24) | (17) | (2700) | . | (26) | . | . | . | . | . | 3 x 10 g |
| SRM 2690 | (39) | . | (6) | (17) | (0.8) | (2000) | . | (25) | . | . | . | . | . | 3 x 10 g |

* IRNT certificates expired, however use and sales continue without problems worldwide

COAL FLY ASH analysis listed in mass %

ACIRS: RM, 80g SABS: CRM, 20g NCS: CRM. 30g

| Number | Al ₂ O ₃ | BaO | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | SiO ₂ | SO ₃ | SrO | TiO ₂ | Units | REDUCING, OXIDIZING TEMPERATURES °C | | | |
|--------------|--------------------------------|------|------|--------------------------------|------------------|------|--------------------------------|-------------------|-------------------------------|------------------|-----------------|------|------------------|-------|-------------------------------------|------------|---------------|------------|
| | | | | | | | | | | | | | | | Deformation | Spherical | Hemispherical | Flow |
| SABS 109 | 37.96 | . | 2.84 | 2.54 | 1.19 | 1.19 | . | 2.07 | 0.31 | 45.66 | 2.57 | . | 1.97 | 20 g | . | . | . | . |
| NCS FC82012a | 29.45 | . | 8.73 | 7.95 | 1.16 | 1.03 | . | 0.79 | 0.28 | 46.99 | 0.69 | . | 1.28 | 30 g | . | . | . | . |
| ACIRS Al * | 28.9 | 0.18 | 6.05 | 14.6 | 0.46 | 1.25 | 0.22 | 0.43 | 1.26 | 44.1 | 0.32 | 0.16 | 1.56 | 80 g | 1257, 1349 | 1287, 1383 | 1309, 1398 | 1367, 1429 |
| NCS FC82013a | 28.54 | . | 5.30 | 13.43 | 1.61 | 0.71 | . | 0.25 | 0.11 | 41.11 | 5.10 | . | 2.62 | 30 g | . | . | . | . |
| NCS FC82014a | 28.09 | . | 5.15 | 6.04 | 1.44 | 1.05 | . | 0.52 | 0.38 | 54.68 | 0.49 | . | 1.14 | 30 g | . | . | . | . |
| NCS FC82015a | 22.90 | . | 4.01 | 6.55 | 1.53 | 1.31 | . | 0.78 | 0.55 | 59.66 | 0.40 | . | 1.13 | 30 g | . | . | . | . |

* ACIRS Al also contains Co:(0.0043) Cr:(0.0058) Cu:(0.0099) Ni:(0.0047) Pb:(0.0047) V:(0.0176) Zn:(0.0090)

INDUSTRIAL FLY ASH

analysis listed in mass % except * which is mg/kg

| Number | Al | As | Ca | Cd | Cr | Hg* | Na | Ni | Pb | Cu | Fe | K | Sb | V | Zn | ZnO |
|--------------|-------|--------|-------|--------|-------|-------|-------|--------|-------|--------|-------|-------|--------|--------|------|-------|
| ECRM 882-1 | 0.375 | 0.0054 | 10.11 | 0.0183 | 0.490 | 0.75 | 0.697 | 0.0263 | 1.324 | 0.218 | 22.20 | 0.960 | 0.0116 | 0.0090 | . | 28.49 |
| JK 43 | (0.2) | . | (12) | 0.0023 | (8) | 3.9 | (0.5) | (2) | 0.21 | (0.2) | (20) | (0.3) | . | (0.02) | 4.96 | . |
| JK 44 | (0.2) | . | (5) | 0.0469 | (0.2) | 2.8 | (1) | (0.02) | 2.74 | (0.2) | (27) | (1.3) | . | (0.02) | 27.3 | . |
| JK 45 | (0.1) | . | (7) | 0.0047 | (0.3) | 0.25 | (7) | (0.05) | 0.11 | (0.01) | (40) | (0.4) | . | (0.1) | 1.53 | . |
| 502-843-1000 | . | . | . | . | . | 0.827 | . | . | . | . | . | . | . | . | . | . |

| Number | Bi | C | Cl | F | Mg | Mn | S | Si | Sn | Units, Class |
|--------------|--------|-------|--------|--------|--------|-----|-------|--------|--------|--------------|
| ECRM 882-1 | 0.0026 | (1.0) | (2.35) | (0.07) | (0.48) | (2) | (0.5) | (1.05) | (0.02) | 100 g, CRM |
| JK 43 | . | . | . | . | . | . | . | . | . | 15 g, CRM |
| JK 44 | . | . | . | . | . | . | . | . | . | 25 g, CRM |
| JK 45 | . | . | . | . | . | . | . | . | . | 15 g, CRM |
| 502-843-1000 | . | 42.4 | . | . | . | . | 0.29 | . | . | 20 g, RM |

RM**COAL-TAR PITCH**

analysis listed in mg/kg except as noted

60 g units

| Number | %C | S% | %H | Ash | Al | As | Br | Ca | Cd | Cl | Cr | Fe | I | K | Mg | Mn | Na | Ni | P | Pb |
|--------------|------|------|-------|------|-----|------|------|----|-------|-----|------|-----|------|-----|-----|------|-----|-----|-----|----|
| DOMTAR CTP A | 94.0 | 0.49 | 4.0 | 0.27 | 245 | . | 1.7 | 91 | . | 118 | 0.87 | 200 | 0.33 | 43 | 17 | 2.7 | 257 | 2.5 | 10 | 91 |
| DOMTAR CTP B | 93.4 | 0.52 | 4.3 | 0.22 | 228 | 9 | 4.8 | 41 | 2.5 | 122 | 1.1 | 280 | 0.6 | 34 | <30 | 3.3 | 150 | . | 3 | 80 |
| DOMTAR CTP C | 83.4 | 4.46 | 10.31 | 0.19 | 9 | 0.18 | 0.25 | 3 | <0.05 | 18 | 0.4 | 14 | 1.4 | 2.2 | <16 | 0.21 | 10 | 76 | 236 | 1 |

continued informational values listed in mg/kg except as noted

| Number | Sb | Si | Sn | Ti | V | Zn | Soft Point °C |
|--------------|------|-----|------|----|------|----|---------------|
| DOMTAR CTP A | . | 358 | . | 18 | 1.2 | 88 | 115 |
| DOMTAR CTP B | 0.57 | 408 | 3.7 | 16 | 0.89 | 90 | 118 |
| DOMTAR CTP C | 0.03 | 20 | <0.7 | 19 | 170 | 1 | 129 |

CRM COATING THICKNESSNumber nominal μm coating thickness

| | | | | |
|-----------|-----|------|------|------|
| SRM 1361b | 6 | 12 | 25 | 48 |
| SRM 1358b | 20 | 80 | 255 | 1000 |
| SRM 1362b | 40 | 80 | 140 | 205 |
| SRM 1359b | 48 | 140 | 505 | 800 |
| SRM 1363b | 255 | 385 | 505 | 635 |
| SRM 1364b | 800 | 1000 | 1525 | 1935 |

These samples are designed for calibrating thickness gauges using magnetic principles. Each sample is a set of four 45 mm x 45 mm plates of coated 1010 sheet steel substrate coated with copper and a thin protective layer of chromium.

CONTINUOUS CASTING POWDER

analysis listed in mass %

IRSID: RM, 100 g units NCS: CRM, 50 g units

| Number | SiO ₂ | Al ₂ O ₃ | C | C.Free | CO ₂ | Ca | F | Fe | K | MgO | Mn | Na | Na ₂ O | P | S | TiO ₂ | LOI |
|-------------|------------------|--------------------------------|-------|--------|-----------------|-------|---------|---------|---------|--------|-------|------|-------------------|---------|---------|------------------|--------|
| NCS HC26805 | 41.31 | 6.93 | 3.06 | 1.57 | . | 21.46 | (4.79) | . | . | 3.26 | . | . | 4.07 | . | . | . | . |
| NCS HC26804 | 34.95 | 5.30 | 15.86 | 14.49 | . | 19.13 | (5.15) | . | . | 0.78 | . | . | 4.99 | . | . | . | . |
| IRSID 2701 | 32.70 | 6.10 | 3.37 | (1.78) | (5.59) | 22.90 | 7.58 | (0.145) | 0.159 | 2.19 | . | 9.42 | . | (0.014) | (0.055) | (0.048) | (2.08) |
| NCS HC26803 | 30.10 | 2.14 | 5.98 | 4.06 | . | 30.78 | (10.59) | . | . | 1.30 | . | . | 0.52 | . | . | . | . |
| IRSID 2702 | 28.70 | 12.60 | 16.54 | 15.80 | (2.53) | 17.80 | 6.08 | 1.260 | (0.750) | (1.47) | 0.071 | 3.61 | . | (0.180) | (0.490) | 0.564 | (1.26) |
| NCS HC26802 | 23.08 | 14.14 | 12.71 | 9.94 | . | 17.93 | (3.86) | . | . | 5.86 | . | . | 2.94 | . | . | . | . |
| NCS HC26801 | 18.96 | 16.99 | 19.97 | 18.14 | . | 12.89 | (4.47) | . | . | 1.39 | . | . | 9.86 | . | . | . | . |

RM**CONTINUOUS CASTING POWDER**

typical analysis listed in mass %

100 g units

| Number | SiO ₂ | Al ₂ O ₃ | Ca | F | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | S | SrO | TiO ₂ | Other |
|----------|------------------|--------------------------------|-------|-------|--------------------------------|------------------|-------|-------|-------------------|-------------------------------|-------|-------|------------------|------------|
| DH X2802 | 57.50 | 3.09 | 25.15 | 0.074 | 0.488 | 0.830 | 0.981 | 0.030 | 1.097 | 0.060 | 0.132 | 0.020 | 0.055 | ZnO: 0.004 |
| DH X2801 | 55.0 | 3.58 | 23.08 | 0.047 | 0.467 | 1.092 | 4.80 | 0.033 | 1.33 | 0.044 | 0.245 | 0.019 | 0.069 | BaO: 0.019 |
| DH 3005 | 43.10 | 5.14 | 27.35 | 6.23 | 0.233 | 0.376 | 0.39 | 0.007 | 7.93 | 0.091 | 0.019 | . | 0.086 | . |
| DH 3010 | 38.56 | 5.05 | 27.06 | . | 2.63 | 0.155 | 4.13 | 0.059 | 5.479 | 0.454 | 0.131 | . | 0.055 | . |
| DH 3013 | 37.70 | 5.95 | 30.73 | 5.72 | 0.437 | 0.288 | 1.93 | 0.045 | 6.43 | 0.047 | 0.077 | . | 0.064 | BaO: 0.121 |

RM**COVER POWDER**

analysis listed in mass %

100 g units

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | S | SiO ₂ | SrO | TiO ₂ |
|---------|--------------------------------|-------|--------------------------------|------------------|-------|--------------------------------|-------------------|-------------------------------|-------|------------------|-------|------------------|
| DH 5905 | 19.32 | 46.50 | 0.435 | 0.321 | 9.17 | 0.051 | . | 0.039 | 0.074 | 22.93 | . | 0.035 |
| DH 5906 | 14.34 | 33.29 | 0.598 | 0.210 | 19.38 | 0.052 | 0.32 | 0.037 | 0.061 | 30.78 | 0.015 | 0.037 |

CRM COKE

analysis listed in mass % except * which is mg/kg

| Number | Al% | Ca% | Cd* | Co* | Cr* | Cu* | Fe% | K% | Mg% | Mn% | Na% | Ni* | P% | Pb* | Si% | Ti% | V* | Zn* |
|-------------|------|------|-----|-----|-----|-----|------|-------|-------|-------|-------|-----|-------|-----|------|-------|----|-----|
| NCS FC28131 | 2.72 | 0.29 | <1 | 7 | 11 | 16 | 0.51 | 0.094 | 0.046 | 0.008 | 0.050 | 13 | 0.015 | . | 3.22 | 0.12 | 27 | 18 |
| NCS FC28129 | 2.34 | 0.60 | . | 7 | 15 | 21 | 0.75 | 0.093 | 0.11 | 0.021 | 0.13 | 15 | 0.020 | 14 | 2.97 | 0.12 | 41 | 11 |
| NCS FC28130 | 1.96 | 0.52 | <1 | 6 | 12 | 17 | 0.63 | 0.061 | 0.11 | 0.015 | 0.063 | 12 | 0.022 | . | 2.35 | 0.099 | 34 | 11 |

COKE ASH

analysis listed in mass %

| Number | Al ₂ O ₃ | CaO | Co ₃ O ₄ | Fe | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | SrO | TiO ₂ | V ₂ O ₅ | LOI |
|-------------|--------------------------------|-------|--------------------------------|------|--------------------------------|------------------|------|-------|-------------------|-------------------------------|-----------------|------------------|-------|------------------|-------------------------------|------|
| NCS FC28137 | 35.62 | 2.82 | . | . | 5.02 | 0.78 | 0.53 | 0.070 | 0.47 | 0.24 | . | 47.81 | . | 1.38 | 0.033 | . |
| NCS FC28136 | 30.66 | 6.00 | . | . | 7.51 | 0.61 | 1.50 | 0.16 | 0.70 | 0.41 | . | 41.61 | . | 1.37 | 0.050 | . |
| NCS FC28135 | 29.95 | 5.67 | . | . | 7.23 | 0.76 | 1.25 | 0.18 | 1.18 | 0.31 | . | 42.87 | . | 1.41 | 0.049 | . |
| DH 3713 | 26.33 | 4.16 | 0.011 | . | 8.57 | 4.120 | 2.94 | 0.172 | 0.568 | 0.215 | S:0.228 | 49.74 | 0.056 | 1.095 | 0.060 | . |
| DH 3711 | 13.79 | 11.60 | 0.007 | 7.79 | . | 3.29 | 8.69 | 0.189 | 3.08 | 0.607 | 0.091 | 43.5 | 0.103 | 2.78 | 0.058 | 0.52 |

| Number | C.tot | CO ₂ | Cr ₂ O ₃ | CuO | NiO | ZnO | ZrO ₂ | Units |
|-------------|-------|-----------------|--------------------------------|-------|-------|-------|------------------|-----------|
| NCS FC28137 | . | . | . | . | . | . | . | CRM, 5 g |
| NCS FC28136 | . | . | . | . | . | . | . | CRM, 5 g |
| NCS FC28135 | . | . | . | . | . | . | . | CRM, 5 g |
| DH 3713 | . | . | 0.046 | 0.037 | 0.026 | 0.232 | 0.028 | RM, 100 g |
| DH 3711 | 0.039 | 0.045 | 0.036 | 0.009 | 0.030 | 0.010 | 0.041 | RM, 100 g |

BaO: 0.170 PbO: 0.209

CRM DUST

| Number | Type | Al | Al ₂ O ₃ | As | C | CaO | Co | Cr | Cr ₂ O ₃ | Cu | CuO | F | Fe | FeO | K |
|--------|------------------|--------|--------------------------------|---------|-------|------|---------|-------|--------------------------------|--------|-------|-------|------|------|-------|
| VS E5 | Blast Furnace | . | 2.87 | . | 13 | 7.9 | 0.013 | . | 0.085 | . | 0.013 | 0.049 | 44.3 | . | . |
| VS E4 | Blast Furnace | . | 2.33 | 0.0018 | 13.2 | 8.8 | . | . | . | . | 0.034 | 0.023 | 44.6 | . | . |
| VS E2 | Converter | (0.07) | . | (0.002) | 1.383 | 7.97 | (0.003) | (0.1) | . | (0.04) | . | (0.5) | 56.4 | 6.2 | (0.2) |
| VS E1 | Electric Furnace | . | 3.06 | (0.004) | 0.684 | 5.85 | (0.03) | . | 20.3 | (0.1) | . | (0.7) | 29.7 | (21) | (0.1) |

| Number | MgO | MnO | Na | Ni | NiO | P | Pb | S | SiO ₂ | Sn | TiO ₂ | V | V ₂ O ₅ | Zn | Units |
|--------|------|------|-------|--------|-------|--------|--------|-------|------------------|-----------|------------------|--------|-------------------------------|-------|-------|
| VS E5 | 2.26 | 0.5 | . | . | 0.022 | 0.041 | . | 0.26 | 7.17 | . | 1.63 | . | 0.39 | 0.27 | 150 g |
| VS E4 | 0.82 | 0.47 | . | . | 0.033 | 0.015 | 0.44 | 7.46 | . | . | 0.2 | . | 0.041 | 1.52 | 150 g |
| VS E2 | 1.64 | 1.41 | (0.1) | (0.03) | . | 0.065 | 0.276 | 0.116 | 1.76 | (<0.0005) | . | (0.01) | . | 0.59 | 100 g |
| VS E1 | 9.3 | 1.56 | (0.1) | . | 3.68 | (0.02) | (0.05) | 0.072 | 10.3 | (<0.0005) | 2.79 | (0.04) | . | (0.2) | 150 g |

RM DUST

typical analysis listed in mass %

* samples list Cu as CuO and Ni as NiO

DH 6203-6205: 20 g

all others: 100 g

| Number | Type | Al ₂ O ₃ | C | CO ₂ | CaO | Cl | Cr ₂ O ₃ | CuO | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | PbO | SiO ₂ | TiO ₂ | ZnO |
|-----------|------------------|--------------------------------|------|-----------------|-------|------|--------------------------------|-------|------------------|-------|-------------------|-------------------------------|-------|------------------|------------------|-------|
| DH X2901 | Blast Furnace | 0.961 | . | . | 5.28 | . | 0.038 | . | 0.778 | 1.147 | 0.119 | 0.153 | 0.006 | 4.28 | 0.068 | 0.267 |
| DH X2902 | Blast Furnace | 0.823 | . | . | 3.12 | . | 0.037 | . | 0.84 | 0.678 | 0.138 | 0.165 | 0.017 | 3.28 | 0.053 | 0.271 |
| DH X2903 | Blast Furnace | 0.701 | . | . | 2.00 | . | 0.040 | 0.006 | 0.705 | 0.502 | 0.111 | 0.158 | 0.018 | 2.44 | 0.058 | 1.19 |
| DH 6205 | Cupola | 1.30 | 6.80 | 3.84 | 4.91 | 2.88 | 0.041 | 0.163 | 3.68 | 1.85 | 2.26 | 0.147 | 2.43 | 34.52 | 0.060 | 21.01 |
| DH 6204 * | Cupola | 1.06 | 8.08 | 2.02 | 2.54 | 3.62 | 0.072 | 0.079 | 4.16 | 1.53 | 2.63 | 0.051 | 3.48 | 26.94 | 0.184 | 30.65 |
| DH 6206 * | Cupola | 0.220 | 2.57 | . | 0.090 | . | 0.048 | 2.021 | 0.086 | 0.020 | 0.085 | 0.191 | . | 0.430 | 0.014 | 91.1 |
| DH 6203 | Electric Furnace | 2.57 | 4.22 | 1.01 | 1.23 | 2.00 | 0.004 | 0.311 | 2.51 | 3.10 | 5.12 | 0.52 | 1.05 | 15.65 | 0.517 | 12.32 |

continued

| Number | F | Fe | Fe ₂ O ₃ | Mn | Mn ₃ O ₄ | MoO ₃ | NiO | S | SO ₃ | SnO ₂ | V ₂ O ₅ | -H ₂ O |
|-----------|-------|-------|--------------------------------|-------|--------------------------------|------------------|--------|-------|-----------------|------------------|-------------------------------|-------------------|
| DH X2901 | . | 59.37 | . | 0.367 | . | . | 0.015 | 0.488 | . | . | 0.020 | . |
| DH X2902 | . | 61.67 | . | 0.341 | . | . | 0.016 | 0.577 | 1.44 | . | . | . |
| DH X2903 | . | 63.01 | . | 0.425 | . | . | 0.012 | 0.392 | . | . | 0.020 | . |
| DH 6205 | 0.096 | . | 9.49 | . | 2.57 | 0.013 | . | . | 2.70 | 0.018 | 0.019 | 0.107 at 900°C |
| DH 6204 * | 0.247 | . | 6.29 | 0.97 | . | . | 0.0162 | 1.09 | . | . | . | 0.055 at 500°C |
| DH 6206 * | . | . | 0.572 | 0.04 | 0.061 | . | 0.297 | 0.305 | . | 0.047 | . | 1.17 at 900°C |
| DH 6203 | 0.570 | . | 36.85 | . | 4.97 | . | . | . | 5.29 | . | 0.004 | 0.214 at 900°C |

analysis for DH 6204 continues in mg/kg

| Number | Cd | Li | Mo | Sr | V |
|---------|-----|-----|-----|-----|------|
| DH 6204 | 197 | 0.4 | 510 | 5.4 | 2.92 |

CRM FLUE DUST

informational analysis listed in mass %

30 g units

| Number | Type | Al ₂ O ₃ | CO ₂ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | SO ₃ | SiO ₂ | TiO ₂ |
|------------|--------------|--------------------------------|-----------------|-------|--------------------------------|------------------|------|------|-------------------|-----------------|------------------|------------------|
| BL 12-1-11 | Sinter Plant | 4.00 | 8.60 | 6.77 | 3.18 | 1.23 | 2.22 | 0.03 | 4.11 | 1.15 | 65.58 | 0.23 |
| BL 12-1-10 | Foundry | 1.64 | 5.39 | 12.80 | 60.95 | 0.28 | 7.59 | 0.16 | 0.15 | 2.22 | 9.80 | 0.075 |

continued certified analysis listed in mg/kg

| Number | Ag | As | Ba | Cd | Co | Cr | Cu | Mo | Ni | Pb | Sb | Sn | Sr | V | Zn |
|------------|-----|-----|-------|-----|----|------|----|------|----|------|-----|------|------|------|----|
| BL 12-1-11 | . | (8) | 160 | (3) | 8 | 3910 | 27 | (10) | 36 | (25) | . | (43) | (58) | 56 | 50 |
| BL 12-1-10 | (1) | (8) | (150) | 5 | 31 | 189 | 76 | (4) | 47 | 56 | (3) | (40) | (50) | (33) | 86 |

CRM FURNACE DUST

analysis listed in mass %

100 g units

| Number | Ag | Al | As | Bi | C | Ca | Cd | Cl | Co | Cr | Cu | F | Fe | H ₂ O | Hg |
|------------|--------|-------|--------|--------|--------|------|--------|-------|--------|-------|--------|-------|-------|------------------|----------|
| ECRM 876-1 | . | 0.034 | 0.023 | . | . | 3.43 | . | . | . | 0.17 | 0.42 | . | 24.85 | . | . |
| ECRM 880-1 | . | 1.28 | . | . | . | 3.15 | . | 0.085 | . | 0.027 | 0.005 | 0.034 | 31.0 | . | . |
| ECRM 884-1 | 0.0028 | 0.379 | 0.0054 | 0.0280 | (0.82) | 5.22 | 0.0045 | 0.991 | 0.0046 | 1.86 | 0.1569 | 0.411 | 31.67 | (0.30) | (0.0002) |

continued

| Number | K | Mg | Mn | Mo | Na | Ni | P | Pb | S | Si | Sn | Ti | V | Zn | LOI |
|------------|-------|-------|-------|-------|-------|-------|-------|-------|--------|-------|--------|--------|--------|-------|--------|
| ECRM 876-1 | 1.63 | 1.31 | 2.84 | . | 1.98 | 0.034 | 0.128 | . | 0.87 | 1.72 | . | 0.048 | . | 23.29 | . |
| ECRM 880-1 | 0.108 | 0.714 | 0.218 | . | 0.041 | 0.014 | 0.038 | 0.017 | 0.425 | 3.34 | . | 0.081 | . | 0.064 | . |
| ECRM 884-1 | 0.979 | 1.848 | 5.85 | 0.208 | 0.585 | 0.197 | 0.079 | 0.442 | (0.49) | 2.100 | 0.0186 | 0.0230 | 0.0303 | 17.50 | (2.94) |

CRM INDOOR DUST

analysis listed in mg/kg

8 g units

| Number | As | Cd | Cr | Hg | Pb |
|----------|------|------|-------|------|------|
| SRM 2584 | 17.4 | 10.0 | 135.0 | 5.20 | 9761 |
| SRM 2583 | 7.0 | 7.3 | 80 | 1.56 | 85.9 |

CRM ROAD DUST

analysis listed in mg/kg except * which is µg/kg

| Number | Pd* | Pt* | Rh* | Cd | Co | Hf | Mo | Rb | Sb | Th | V | Y | Units |
|---------|-----|------|------|-------|--------|-------|--------|------|--------|-------|--------|--------|-------|
| BCR 723 | 6.1 | 81.3 | 12.8 | (2.5) | (29.8) | (2.2) | (40.0) | (75) | (28.2) | (4.8) | (74.9) | (12.5) | 25 g |

continued analysis listed in mass %

| Number | Al | Ba | Cr | Fe | Mn | Ni | Pb | Sr | Ti | Zn | Moisture |
|---------|--------|---------|----------|--------|---------|----------|----------|----------|---------|---------|----------|
| BCR 723 | (3.75) | (0.046) | (0.0440) | (3.29) | (0.128) | (0.0171) | (0.0866) | (0.0254) | (0.258) | (0.166) | (~3%) |

CRM USED AUTOMOBILE EXHAUST CATALYST

mg/kg

| Number | Pb | +/- | Pd | +/- | Pt | +/- | Rh | +/- | Units |
|------------|-------|-----|-------|-----|-------|-----|-------|-----|-------|
| SRM 2557 | 13931 | 97 | 233.2 | 1.9 | 1131 | 11 | 135.1 | 1.9 | 70 g |
| SRM 2556 | 6228 | 49 | 326.0 | 1.6 | 697.4 | 2.3 | 512 | 0.5 | 70 g |
| FLX CRM132 | . | . | 1673 | 43 | 1770 | 18 | 295 | 6 | 30 g |
| ERM-EB504a | . | . | 1596 | 11 | 1414 | 9 | 210.0 | 2.2 | 200 g |
| FLX CRM133 | . | . | 1075 | 33 | 465 | 32 | 242 | 4 | 30 g |

RM ELECTRODE CARBON

| Number | Size Analysis | Ash | Bulk Density | Relative Density | Moisture | Volatile Matter | Units |
|----------|---------------|------|--------------|------------------|----------|-----------------|-------|
| ACIRS EC | <1% @ +212 µm | 2.8% | 1028 g/L | 2.08 | 0.2% | 0.8% | 1 kg |

ELECTRONIC SCRAP POWDER

analysis listed in mass %

200 g units

| Number | Ag | As | Au | Be | Cd | Cr | Cu | Hg | In | Ni | Pb | Pd | Pt | Sn |
|-----------|--------|--------|---------|---------|---------|-------|-------|-----------|----------|-------|------|---------|---------|-------|
| BAM M505a | 0.0633 | 0.0372 | 0.00524 | 0.00068 | 0.00164 | 0.980 | 16.76 | (<0.0005) | (0.0043) | 0.694 | 1.13 | 0.00480 | 0.00057 | 0.468 |

FERROBORON

= class, where 1 = CRM and 2 = RM typical analysis

| # | Number | B | Al | C | Cr | Cu | Fe | Mn | Ni | P | S | Si | Sn | Ti | V | W | Zn | Units |
|---|-------------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|-------|--------|-------|--------|-------|--------|-------|
| 1 | VS F21/2 | 20.91 | 1.546 | 0.047 | . | 0.012 | . | . | . | 0.0119 | . | 0.73 | 0.0055 | . | . | . | 0.0055 | 100 g |
| 1 | NCS HC11612 | 20.82 | 0.018 | 0.086 | . | . | . | 0.305 | . | 0.020 | 0.0018 | 0.353 | . | . | . | . | . | 70 g |
| 1 | NCS HC25658 | 20.58 | 0.99 | 0.022 | . | . | . | . | . | 0.017 | 0.0016 | 1.68 | . | . | . | . | . | 50 g |
| 2 | DH 1705 | 20.23 | 0.085 | 0.62 | 0.925 | 0.072 | 75.71 | 0.416 | 0.097 | 0.061 | . | 1.011 | . | 0.025 | 0.006 | 0.127 | . | 50 g |
| 1 | NCS HC28632 | 19.33 | 0.185 | 0.19 | 0.26 | 0.050 | . | 0.575 | 0.056 | 0.027 | 0.0020 | 0.65 | . | 0.030 | 0.010 | . | . | 50 g |
| 2 | DH 1704 | 19.13 | 0.064 | 0.322 | 0.038 | 0.013 | 79.74 | 0.299 | 0.012 | 0.015 | . | 0.229 | . | 0.011 | . | 0.015 | . | 50 g |
| 1 | NCS HC28631 | 18.92 | 0.036 | 0.40 | 0.025 | 0.015 | . | 0.310 | 0.013 | 0.025 | 0.0023 | 0.18 | . | 0.017 | 0.0090 | . | . | 50 g |
| 1 | NCS HC93623 | 18.69 | 0.083 | 0.45 | . | . | . | . | . | 0.025 | 0.0044 | 0.44 | . | 0.019 | . | . | . | 50 g |
| 1 | ECRM 587-1 | 18.67 | 0.047 | 0.738 | 0.104 | . | . | 0.272 | . | 0.020 | . | 0.129 | . | 0.039 | 0.004 | . | . | 100 g |
| 1 | NCS HC11613 | 17.65 | 0.035 | 0.181 | . | . | . | 0.35 | . | 0.030 | 0.0023 | 0.549 | . | . | . | . | . | 70 g |
| 1 | VS F22/3 | 8.95 | 7.78 | 0.161 | . | 3.43 | . | . | . | 0.021 | 0.018 | 7.82 | . | . | . | . | . | 100 g |

FERROCHROMIUM

= class, where 1 = CRM and 2 = RM

chips as noted

all others: powder

| # | Number | Cr | Al | C | Co | Cu | Fe | Mn | N | Ni | P | S | Si | Ti | V | Units |
|---|-----------------|-------|-------|---------|-----------|------------|-------|----------|-----------|------------|---------|----------|--------|----------|--------|--------------------|
| 1 | IRSID 509-1 | 72.85 | . | 0.012 | . | . | . | . | 0.029 | . | (0.019) | . | 0.230 | . | . | 100 g |
| 1 | ECRM 580-1 | 72.18 | . | 0.019 | 0.047 | . | . | . | 0.035 | . | 0.011 | . | 0.306 | . | 0.083 | 100 g |
| 1 | VS F11/4 | 71.41 | . | 0.119 | . | . | 27.20 | 0.161 | 0.035 | 0.362 | 0.023 | 0.0021 | 0.68 | . | 0.084 | 100 g |
| 2 | DH 1602 | 71.20 | 0.017 | 8.12 | Nb:0.006 | . | 19.38 | 0.103 | W:0.038 | 0.304 | 0.014 | Mg:0.049 | 0.153 | 0.009 | 0.102 | 50 g |
| 1 | IPT 65 | 71.2 | 9.2 | 0.051 | 0.016 | . | 17.9 | 0.128 | . | 0.077 | 0.006 | 0.016 | 0.71 | . | . | 100 g |
| 1 | BCS 203/6 | 71.01 | . | 0.0270 | 0.0442 | . | . | 0.153 | . | 0.218 | 0.0195 | 0.004 | 0.381 | . | 0.0729 | 100 g |
| 1 | SRM 196 | 70.83 | . | 0.035 | . | . | . | . | . | 0.020 | 0.003 | 0.003 | 0.373 | . | . | 100 g |
| 1 | IRSID 507-1 | 70.30 | . | 5.40 | . | . | 0.270 | 0.049 | . | 0.017 | . | 1.20 | . | . | . | 100 g |
| 1 | VS F47 | 69.8 | . | 8.80 | . | . | . | 0.020 | . | 0.025 | 0.036 | 0.103 | . | . | . | 100 g |
| 2 | DFS 1 | 69.8 | . | 0.0549 | . | H:(0.0024) | . | O:(0.56) | 0.0973 | . | (0.027) | (0.0017) | (0.87) | . | . | 50 g |
| 2 | DFS 2 | 69.6 | . | 9.03 | . | . | . | O:(0.2) | (0.0178) | . | (0.021) | (0.0174) | (0.2) | . | . | 50 g |
| 1 | JSS 733-1 | 69.24 | . | 8.35 | . | . | . | . | . | 0.0219 | 0.0248 | (0.262) | . | . | . | 150 g |
| 1 | NCS HC37618 | 69.12 | . | 8.44 | . | . | . | . | . | 0.026 | 0.026 | 0.30 | . | . | . | 50 g |
| 1 | VS F45 | 69.0 | 0.041 | 0.071 | . | . | . | 0.082 | . | 0.027 | 0.0024 | 1.10 | . | . | . | 100 g |
| 1 | VS F35/2 | 68.9 | . | 0.79 | . | . | 28.16 | . | 0.185 | . | 0.027 | 0.0022 | 0.181 | . | . | 100 g chips |
| 1 | BCS 204/6 | 68.77 | . | 8.57 | 0.0359 | . | . | 0.225 | . | 0.289 | 0.0164 | 0.0190 | 0.715 | 0.0565 | 0.094 | 100 g |
| 1 | NCS HC93604 | 68.13 | . | 0.038 | . | . | . | 0.103 | . | 0.260 | 0.022 | 0.015 | 0.466 | 0.016 | . | 50 g |
| 1 | VS F15/2 | 68.1 | 0.30 | 0.080 | . | . | . | . | 1.79 | . | 0.036 | 0.0019 | 2.08 | . | . | 100 g chips |
| 1 | SRM 64c | 68.00 | . | 4.68 | 0.051 | 0.005 | 24.98 | 0.16 | 0.045 | 0.43 | 0.020 | 0.067 | 1.22 | 0.002 | 0.15 | 100 g |
| 1 | NCS HC26607b | 65.86 | . | 7.23 | . | . | . | 0.38 | . | 0.29 | 0.017 | 0.020 | 2.18 | 0.15 | . | 50 g |
| 1 | NCS HC25603b | 65.27 | . | 7.37 | 0.044 | . | 24.90 | 0.31 | . | 0.39 | 0.020 | 0.015 | 1.27 | 0.104 | 0.138 | 50 g |
| 1 | NCS HC37607 | 64.32 | . | 0.086 | . | . | . | 0.35 | . | . | 0.034 | 0.005 | 1.39 | . | . | 40 g chips |
| 1 | VS F50 | 64.3 | 0.168 | 0.079 | . | . | . | . | 11.2 | . | 0.0048 | 0.0058 | 1.29 | . | . | 100 g |
| 2 | DFS 3 | 64.2 | . | (0.068) | . | . | . | O:(3) | 9.71 | . | (0.023) | (0.016) | (0.54) | . | . | 50 or 100 g |
| 1 | NCS HC11606 | 64.17 | . | 6.37 | . | . | . | 0.32 | . | . | 0.023 | 0.013 | 4.29 | . | . | 100 g |
| 1 | NCS HC37608 | 64.06 | . | 0.243 | . | . | . | 0.34 | . | . | 0.039 | 0.019 | 0.99 | . | . | 40 g chips |
| 1 | NCS HC37615 | 63.50 | . | 8.78 | . | . | . | . | . | 0.023 | 0.025 | 0.78 | 0.097 | . | . | 50 g |
| 1 | NCS HC25635a | 63.44 | . | 0.218 | . | . | . | 0.44 | . | 0.035 | 0.0033 | 1.07 | . | 0.093 | . | 50 g chips |
| 1 | NCS HC25651 | 63.31 | . | 2.55 | . | . | . | 0.47 | . | 0.023 | 0.047 | 2.04 | . | . | . | 50 g |
| 1 | NM 230 | 62.89 | . | 7.86 | . | . | . | 0.31 | . | . | 0.035 | 0.037 | 3.46 | . | . | 100 g |
| 1 | NCS HC25636a | 62.81 | . | 0.078 | 0.16 | 0.032 | 35.56 | 0.39 | . | 0.34 | 0.028 | 0.0052 | 0.32 | . | 0.098 | 50 g chips |
| 1 | NCS HC25644 | 62.57 | . | 0.0064 | . | . | . | 0.207 | 8.69 | . | 0.024 | 0.029 | 0.75 | . | . | 50 g |
| 1 | NCS HC28621 | 62.54 | . | 7.78 | . | . | 27.09 | 0.307 | . | . | 0.026 | 0.033 | 1.45 | 0.166 | 0.138 | 50 g |
| 1 | NCS HC25653 | 62.49 | . | 8.70 | . | . | . | 0.11 | . | 0.025 | 0.024 | 0.15 | 0.016 | . | . | 50 g |
| 1 | NCS HC37609 | 61.54 | . | 8.49 | . | . | . | 0.20 | . | . | 0.024 | 0.015 | 2.15 | 0.24 | . | 40 g |
| 1 | NCS HC93611 | 60.42 | . | 8.13 | . | . | . | 0.21 | . | . | 0.022 | 0.059 | 0.92 | . | . | 50 g |
| 1 | NCS HC28622 | 60.00 | . | 7.72 | . | . | 28.65 | 0.340 | . | 0.025 | 0.033 | 2.43 | 0.261 | 0.153 | . | 50 g |
| 1 | NCS HC93605 | 59.71 | . | 8.00 | . | . | . | 0.308 | . | 0.312 | 0.037 | 0.037 | 2.94 | 0.410 | . | 50 g |
| 1 | NCS HC37617 | 59.35 | . | 8.32 | . | . | . | . | . | 0.021 | 0.040 | 1.91 | . | . | . | 50 g |
| 1 | NCS HC11610 | 58.83 | . | 0.040 | . | . | . | 0.313 | 4.76 | . | 0.020 | 0.042 | 0.525 | . | . | 70 g |
| 1 | NCS HC28630 | 58.65 | . | 0.45 | Ca:0.0012 | 0.020 | . | 0.35 | As:0.0027 | 0.242 | 0.030 | 0.0015 | 1.43 | Zn:0.010 | 0.154 | 50 g |
| 1 | NCS HC93609 | 58.28 | . | 8.36 | . | . | . | 0.207 | . | . | 0.023 | 0.068 | 1.15 | . | . | 50 g |
| 1 | NCS HC14615a | 57.44 | . | 8.00 | . | . | 31.46 | 0.24 | . | 0.23 | 0.017 | 0.015 | 1.96 | . | 0.28 | 50 g |
| 1 | NCS HC28619 | 56.76 | . | 7.28 | . | . | 30.22 | 0.300 | . | 0.023 | 0.024 | 4.25 | 0.412 | 0.203 | . | 50 g |
| 1 | NCS HC14615 | 56.16 | . | 8.07 | . | . | . | 0.23 | . | . | 0.017 | 0.045 | 2.30 | . | . | 50 g |
| 1 | NCS HC14614 | 55.81 | . | 7.67 | . | . | . | 0.28 | . | . | 0.018 | 0.032 | 2.57 | . | . | 50 g |
| 1 | NCS HC28620 | 55.77 | . | 7.60 | . | . | 31.41 | 0.382 | . | 0.022 | 0.031 | 3.95 | 0.423 | 0.175 | . | 50 g |
| 1 | NCS HC14613 | 54.04 | . | 7.56 | . | . | . | 0.18 | . | 0.020 | 0.032 | 2.58 | . | . | . | 50 g |
| 1 | BS 130/2 | 52.6 | . | 7.70 | . | (0.012) | . | 0.44 | . | 0.015 | 0.045 | 2.12 | (0.18) | (0.36) | . | 100 g |
| 1 | BS 130/1 | 51.6 | . | 7.12 | . | (0.015) | . | 1.25 | . | Mo:(0.005) | 0.014 | 0.034 | 4.41 | (0.38) | (0.32) | 100 g 17025 |
| 1 | BS 130/3 | 49.01 | . | 6.47 | . | (0.01) | . | 0.76 | . | Mo:(0.003) | 0.015 | 0.028 | 6.20 | (0.41) | (0.38) | 100 g 17025 |
| 1 | SARM 144 | 49 | . | . | . | . | 36.9 | 0.33 | . | 0.21 | . | . | 5.31 | 0.67 | 0.37 | 100 g |
| 1 | VS F38 | . | . | 4.62 | . | . | . | . | . | . | . | . | . | . | . | 100 g last |
| 1 | VS F12/3 | . | . | 0.289 | . | . | . | . | . | . | . | . | . | . | . | 100 g chips |
| 1 | VS F10/2 | . | . | 0.018 | . | . | . | . | . | . | . | . | . | . | . | 100 g chips |
| 1 | VS F9/2 | . | . | 0.012 | 0.042 | . | . | . | . | . | . | . | . | . | . | 100 g chips |
| # | Number | Cr | Al | C | Co | Cu | Fe | Mn | N | Ni | P | S | Si | Ti | V | Units |

CRM FERRONICKEL

| Number | Ni | N | C | Co | Cr | Cu | Fe | Mn | P | S | Si | Units |
|-------------|-------|-------|--------|-------|------|-------|------|-------|-------|-------|------|-------|
| VS F41 | 91.4 | 0.058 | 0.0124 | 2.04 | . | 0.47 | 5.68 | . | . | 0.132 | . | 100 g |
| NCS HC11617 | 16.45 | . | 1.85 | 0.241 | 1.87 | 0.021 | . | 0.041 | 0.037 | 0.213 | 3.11 | 60 g |
| NCS HC11616 | 13.34 | . | 2.12 | 0.247 | 1.98 | 0.022 | . | 0.051 | 0.039 | 0.283 | 3.25 | 60 g |
| NCS HC25656 | 12.16 | . | 3.06 | . | 3.62 | . | . | . | 0.046 | 0.245 | 1.04 | 50 g |
| NCS HC11618 | 10.70 | . | 1.65 | 0.198 | 1.56 | 0.021 | . | 0.053 | 0.032 | 0.211 | 2.54 | 60 g |

FERRONIUBIUM

= class, where 1 = CRM and 2 = RM

* notes the total of Nb+Ta

| # | Number | Nb | Fe | Si | Al | C | Cr | Cu | Mn | P | Pb | Sn | Ta | Ti | V | W | Zr |
|---|-------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|---------|-------|-------|-------|-------|
| 1 | NCS HC25650 | 66.34 | . | 1.01 | 0.89 | 0.074 | . | 0.023 | . | 0.085 | . | . | (0.081) | 0.49 | . | . | . |
| 1 | NCS HC18606 | 66.24 | . | 1.09 | 1.35 | 0.070 | . | . | 0.29 | 0.159 | . | . | 0.084 | 0.78 | . | . | . |
| 1 | NCS HC11609 | 64.89 | . | 1.34 | 0.711 | 0.114 | . | 0.059 | 0.37 | 0.172 | . | . | 0.087 | 0.870 | . | . | . |
| 1 | NCS HC93607 | 64.60 | . | 1.04 | 1.50 | 0.101 | . | 0.038 | . | 0.194 | . | . | 0.097 | 0.585 | . | . | . |
| 1 | VS F20/3 | 63.5* | 33.3 | 0.67 | 0.35 | 0.136 | . | . | . | 0.039 | . | 0.0014 | 63.5* | 0.292 | . | . | . |
| 1 | ECRM 579-1 | 62.87 | . | 1.03 | 1.86 | 0.037 | . | . | . | 0.064 | . | 0.344 | 3.85 | 0.567 | . | . | . |
| 2 | DH 2815 | 60.15 | 28.77 | 1.580 | 4.82 | 0.043 | 0.028 | 0.209 | 0.842 | 0.065 | 0.140 | 1.38 | 0.856 | 0.185 | 0.013 | 0.064 | 0.105 |
| 1 | ECRM 576-1 | 43.90 | . | 1.79 | 2.53 | 0.201 | . | . | . | . | . | 0.195 | 0.306 | 1.32 | . | . | . |

| Number | Co | Mg | Mo | N | Ni | S | Units |
|-------------|--------|-------|-------|-------|-------|--------|-------|
| NCS HC25650 | . | . | . | . | . | 0.028 | 50 g |
| NCS HC18606 | . | . | . | . | . | 0.008 | 50 g |
| NCS HC11609 | . | . | . | . | . | 0.014 | 70 g |
| NCS HC93607 | . | . | . | . | . | 0.013 | 50 g |
| VS F20/3 | 0.0056 | . | . | 0.067 | . | 0.0091 | 100 g |
| ECRM 579-1 | 0.005 | . | . | . | . | 0.021 | 100 g |
| DH 2815 | . | 0.076 | 0.020 | . | 0.019 | 0.056 | 50 g |
| ECRM 576-1 | . | . | . | . | . | . | 100 g |

FERROPHOSPHORUS

= class, where 1 = CRM and 2 = RM

analysis listed in mass %

| # | Number | P | Fe | Al | C | Ca | Cr | Cu | Mn | Nb | Ni | S | Si | Ti | V | Units |
|---|-------------|-------|------|-------|-------|-------|-------|-------|-------|-------|-------|--------|-------|-------|-------|-------|
| 1 | NCS HC93622 | 27.50 | . | . | 0.228 | . | 0.226 | . | 0.70 | . | . | 0.017 | 0.156 | 0.53 | . | 50 g |
| 1 | SRM 90 | 26.2 | . | . | . | . | . | . | . | . | . | . | . | . | . | 75 g |
| 1 | NCS HC11614 | 25.81 | . | . | 0.032 | . | . | . | 0.638 | . | . | 0.0038 | 0.60 | 2.14 | . | 70 g |
| 2 | DH 2204 | 25.69 | 69.2 | 0.056 | 0.043 | 0.681 | 0.077 | 0.214 | 0.964 | 0.038 | 0.047 | 0.023 | 1.238 | 0.836 | 0.150 | 50 g |
| 1 | NCS HC11615 | 21.49 | . | . | 0.130 | . | . | . | 1.07 | . | . | 0.061 | 0.382 | 0.62 | . | 70 g |
| 1 | VS F28/2 | 16.05 | . | . | . | . | . | . | 1.20 | . | . | 0.021 | 1.11 | . | . | 100 g |

FERROTITANIUM

= Class, where 1 = CRM and 2 = RM

| # | Number | Ti | Al | Sol.Al | C | Co | Cr | Cu | Fe | Mn | P | S | Si | V | Zr |
|---|------------------|-------|-------|--------|-------|--------|-------|---------|-------|-------|---------|---------|-------|-------|-------|
| 2 | DH 2409 | 72.74 | 2.93 | . | . | . | 0.384 | 0.074 | 19.27 | 0.192 | 0.004 | . | 0.180 | 1.167 | 0.383 |
| 1 | VS F30/4 | 70.3 | 3.83 | . | 0.154 | . | 0.154 | 0.065 | 21.51 | 0.189 | 0.0030 | 0.0054 | 0.163 | 2.29 | 0.231 |
| 1 | NCS HC15601 | 70.02 | 0.3 | . | 0.057 | . | 0.039 | 0.037 | 26.57 | 0.106 | 0.0071 | 0.0047 | 1.47 | 0.011 | . |
| 1 | NCS HC19604 | 43.81 | 10.64 | . | 0.041 | . | . | . | . | 1.59 | 0.051 | 0.011 | 3.46 | 0.158 | . |
| 1 | NCS HC19605 | 38.81 | 8.61 | . | 0.032 | . | . | 0.025 | . | 0.81 | 0.032 | 0.009 | 4.20 | 0.303 | . |
| 1 | ECRM 584-1 | 37.17 | 7.19 | (6.0) | 0.044 | . | . | . | . | 1.13 | 0.032 | 0.030 | 1.80 | . | . |
| 1 | NCS HC93608 | 32.22 | 3.00 | . | 0.095 | . | . | 0.281 | . | 0.255 | 0.014 | 0.015 | 0.30 | . | . |
| 1 | VS F43 | 31.9 | 11.11 | . | 0.098 | . | 0.354 | 0.336 | . | 1.22 | 0.038 | 0.0058 | 2.50 | 0.152 | 0.059 |
| 1 | NCS HC26613 | 30.24 | 8.13 | . | 0.019 | . | . | (0.005) | . | 1.11 | 0.020 | 0.012 | 1.84 | 0.19 | . |
| 1 | NCS HC18604 | 27.93 | 5.38 | . | 0.065 | . | . | 0.117 | . | 2.67 | 0.043 | 0.013 | 4.68 | . | . |
| 1 | NCS HC28638 | 27.34 | 7.82 | . | 0.033 | . | 0.055 | . | . | 0.362 | 0.015 | 0.0048 | 4.51 | 0.15 | . |
| 1 | VS F42 | 27.13 | 11.41 | . | 0.55 | . | 2.22 | 1.32 | . | 1.1 | 0.05 | 0.023 | 6.74 | . | . |
| 1 | IRSID 510-1 | 26.95 | (4.9) | . | 0.058 | . | . | . | . | . | (0.035) | . | 4.65 | . | . |
| 1 | NM 341 | 24.91 | 5.54 | . | . | . | . | . | . | . | . | . | 2.55 | . | . |
| 1 | BS FeTi-1 | 20.0 | 12.5 | . | 0.57 | (0.03) | 0.33 | 0.60 | . | 7.7 | (0.05) | (0.009) | 2.8 | 0.69 | 3.7 |
| 1 | BS FeTi-2 | 19.6 | 12.6 | . | 0.455 | 0.037 | 0.30 | 0.43 | . | 7.9 | (0.05) | (0.01) | 3.2 | 0.76 | 3.8 |

| Number | B | Ca | Mg | Mo | N | Nb | Ni | Pb | Sn | W | Zn | Units |
|------------------|------|------|-------|--------|--------|--------|-------|----|-------|-------|--------|-------|
| DH 2409 | . | . | 0.070 | 0.814 | . | 0.072 | 0.047 | . | 0.246 | . | . | 50 g |
| VS F30/4 | . | . | . | 0.60 | 0.38 | . | 0.053 | . | 0.077 | . | . | 100 g |
| NCS HC15601 | . | . | . | 0.028 | . | . | 0.29 | . | . | . | . | 50 g |
| NCS HC19604 | . | . | . | . | . | . | . | . | 0.056 | . | . | 100 g |
| NCS HC19605 | . | . | . | . | . | . | . | . | 0.061 | . | . | 100 g |
| ECRM 584-1 | . | . | . | . | . | . | . | . | . | . | . | 100 g |
| NCS HC93608 | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| VS F43 | . | . | . | 0.0036 | 0.085 | . | . | . | 0.013 | 0.032 | . | 100 g |
| NCS HC26613 | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| NCS HC18604 | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| NCS HC28638 | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| VS F42 | . | . | . | 0.106 | . | . | . | . | 0.33 | 0.129 | . | 100 g |
| IRSID 510-1 | . | . | . | . | . | . | . | . | . | . | . | 100 g |
| NM 341 | . | . | . | . | . | . | . | . | . | . | . | 100 g |
| BS FeTi-1 | 0.60 | 1.14 | (0.4) | 0.058 | 0.144 | (0.05) | 0.17 | . | 0.11 | . | (0.03) | 100g |
| BS FeTi-2 | 1.10 | 0.98 | (0.4) | 0.15 | (0.15) | 0.036 | 0.156 | . | 0.160 | . | (0.03) | 100g |

17025

17025

CRM FERROTUNGSTEN

| Number | W | Si | Al(tot) | As | C | Cu | Fe | Mn | Mo | P | Pb | S | Sb | Sn | Units |
|--------------|-------|------|---------|-------|--------|--------|--------|-------|-------|--------|--------|---------|-------|-------|-------|
| ECRM 555-1 | 79.9 | 1.75 | 0.14 | . | 0.025 | . | (15.2) | . | . | (0.02) | . | (0.018) | . | 0.034 | 100 g |
| ECRM 590-1 | 79.55 | 1.05 | (0.36) | . | 0.0250 | 0.0484 | . | 0.136 | 0.101 | . | . | . | . | 0.045 | 100 g |
| NCS HC25606a | 76.24 | 0.34 | . | 0.041 | 0.036 | 0.079 | . | 0.102 | . | 0.033 | . | 0.052 | . | 0.041 | 50 g |
| VS F48 | 71.0 | 0.47 | 0.64 | 0.037 | 0.074 | 0.096 | . | 0.695 | 0.047 | 0.035 | 0.0048 | 0.211 | 0.014 | 0.031 | 100 g |

FERROVANADIUM

= Class, where 1 = CRM and 2 = RM

| # | Number | V | Fe | Si | Al | C | Cr | Cu | Mg | Mn | Mo | N | Ni | P | S |
|---|------------------|-------|-------|-------|----------|-------|-------|--------|-------|-------|--------|------|-------|--------|---------|
| 1 | NCS HC93629 | 80.90 | . | 0.86 | 1.33 | 0.032 | . | . | . | 0.046 | . | . | . | 0.036 | 0.014 |
| 2 | DH 2510 | 80.85 | 14.25 | 0.894 | 0.783 | 0.120 | 0.201 | 0.038 | 0.010 | 0.154 | 0.029 | . | 0.009 | 0.051 | 0.016 |
| 1 | VS F40 | 80.1 | . | 1.31 | 2.12 | 0.096 | 0.185 | 0.81 | . | 1.49 | . | . | . | 0.022 | 0.014 |
| 1 | NCS HC11608 | 79.27 | . | 0.653 | 1.41 | 0.109 | . | 0.0089 | . | 0.106 | . | . | 0.010 | 0.021 | 0.035 |
| 1 | NCS HC28633 | 54.02 | . | 0.682 | 0.0026 | 0.285 | 0.110 | 0.054 | . | 0.663 | . | . | 0.011 | 0.056 | 0.0044 |
| 1 | NCS HC26608c | 53.78 | . | 0.81 | (0.0025) | 0.17 | 0.71 | . | . | 2.00 | . | . | . | 0.043 | 0.0040 |
| 1 | NM 351 | 52.10 | . | . | . | . | . | . | . | . | . | . | . | . | . |
| 1 | NCS HC19606 | 51.14 | . | 0.68 | 0.084 | 0.565 | 0.32 | . | . | 0.43 | . | . | . | 0.087 | 0.010 |
| 1 | NCS HC26608b | 50.57 | . | 0.84 | (0.002) | 0.22 | 0.70 | . | . | 1.64 | . | . | . | 0.051 | 0.0044 |
| 1 | NCS HC93628 | 50.24 | . | 0.730 | 6.10 | 0.130 | . | . | . | 0.474 | . | . | . | 0.042 | 0.016 |
| 1 | ECRM 577-1 | 50.16 | . | 1.79 | 0.414 | 0.089 | . | 0.054 | . | 0.158 | . | . | 0.053 | 0.035 | 0.034 |
| 1 | NCS HC93628a | 50.09 | . | 0.730 | 6.03 | 0.152 | . | . | . | 0.475 | . | . | . | 0.043 | 0.017 |
| 1 | NCS HC37616 | 49.72 | . | 0.50 | 5.18 | 0.081 | . | . | . | 0.58 | . | . | . | 0.016 | 0.012 |
| 1 | NCS HC11607 | 49.40 | . | 1.67 | . | 0.235 | . | 0.022 | . | 0.321 | . | . | . | 0.121 | 0.010 |
| 1 | NCS HC18608 | 48.93 | . | 0.76 | 0.158 | 0.40 | . | . | . | 0.26 | . | . | . | 0.049 | 0.043 |
| 1 | NCS HC28634 | 47.32 | . | 1.89 | 0.0061 | 0.475 | 0.289 | 0.064 | . | 0.365 | . | . | 0.067 | 0.093 | 0.014 |
| 1 | BS FeV 45 | 45.1 | 33.7 | 4.90 | 0.017 | 0.241 | 5.82 | 0.40 | 0.014 | 4.12 | 0.0079 | 0.26 | 4.32 | (0.13) | 0.334 |
| 1 | VS F19/3 | 42.6 | . | 1.47 | (0.005) | 0.418 | 1.21 | 0.204 | . | 3.30 | . | . | . | 0.059 | 0.0102 |
| 1 | BS FeV 42 | 42.2 | 39.2 | 3.77 | (0.05) | 0.297 | 5.18 | 0.31 | 0.059 | 3.37 | 0.023 | 0.19 | 3.87 | 0.127 | 0.31 |
| 1 | VS F32/3 | 40.2 | (40) | (1.2) | (<0.05) | (0.4) | . | (0.2) | . | 3.14 | . | 7.51 | . | (0.05) | (0.008) |

| Number | As | Ca | O | Pb | Ti | Units |
|------------------|----------|---------|---|--------|-------------|-------|
| NCS HC93629 | . | . | . | . | . | 25 g |
| DH 2510 | . | . | . | . | 0.071 | 50 g |
| VS F40 | . | . | . | . | . | 100 g |
| NCS HC11608 | 0.0024 | . | . | . | . | 70 g |
| NCS HC28633 | 0.0017 | 0.022 | . | 0.0006 | . | 50 g |
| NCS HC26608c | . | . | . | . | Zn:(0.004) | 50 g |
| NM 351 | . | . | . | . | . | 100 g |
| NCS HC19606 | . | . | . | . | . | 50 g |
| NCS HC26608b | . | . | . | . | Zn:(0.0024) | 50 g |
| NCS HC93628 | . | . | . | . | . | 30 g |
| ECRM 577-1 | . | . | . | . | . | 100 g |
| NCS HC93628a | . | . | . | . | . | 30 g |
| NCS HC37616 | . | . | . | . | . | 50 g |
| NCS HC11607 | 0.021 | . | . | . | . | 70 g |
| NCS HC18608 | . | . | . | . | . | 50 g |
| NCS HC28634 | 0.024 | 0.115 | . | 0.0004 | . | 50 g |
| BS FeV 45 | (0.013) | 0.010 | . | . | 0.021 | 100 g |
| VS F19/3 | 0.0009 | . | . | . | . | 100 g |
| BS FeV 42 | (0.01) | (0.052) | . | . | 0.033 | 100 g |
| VS F32/3 | (<0.001) | . | . | . | . | 100 g |

Co: 0.008 Nb: 0.013 W: 0.025

17025, 34

17025, 34

CRM RARE EARTH FERROSILICON

* VS F31/2 lists Rare Earth Oxides

| Number | RE | Si | Fe | Ca | Mg | Mn | Ti | Al | C | Ce | Cu | La | Units |
|-------------|-------|-------|---------|------|-------|-------|-------|------|-------|--------|------|--------|-------|
| VS F31/3 * | 39.0 | 39.6 | 16.26 | 1.76 | 0.320 | . | . | 7.60 | 0.032 | 15.65 | 0.51 | . | 100 g |
| NCS HC39602 | 21.20 | 37.18 | 22.18 | 1.98 | 10.56 | 3.43 | 1.92 | . | . | . | . | . | 100 g |
| NCS HC39601 | 20.09 | 40.31 | 20.81 | 3.21 | 9.50 | 2.72 | 1.50 | . | . | . | . | . | 100 g |
| NCS HC28615 | 20.00 | 41.02 | . | 5.60 | . | 0.390 | 0.235 | . | . | . | . | . | 100 g |
| NCS HC28609 | 8.66 | 43.90 | (31.67) | 1.01 | 10.20 | 0.70 | 0.54 | . | . | . | . | . | 80 g |
| NCS HC28612 | 6.42 | 43.44 | (36.43) | 0.90 | 8.25 | 0.63 | 0.435 | . | . | . | . | . | 80 g |
| NCS HC28611 | 5.10 | 43.22 | (40.7) | 0.84 | 5.70 | 0.55 | 0.362 | . | . | . | . | . | 80 g |
| NCS HC28610 | 3.71 | 42.05 | (43.4) | 0.76 | 5.52 | 0.46 | 0.275 | . | . | (1.86) | . | (0.88) | 80 g |

last

FERROSILICOALUMINUM, FERROSILICOCALCIUM, FERROSILICOCHROMIUM, and FERROSILICOTITANIUM

= class, where 1 = CRM and 2 = RM

DH, NCS: 50 g units VS: 100 g units

| # | Number | Si | Fe | Ca | Cr | Ti | Al | C | Cu | Mg | Mn | Mo | Ni | P | S | V | Zr |
|---|-------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|------|-------|-------|-------|--------|-------|-------|
| 2 | DH 2902 | 59.25 | 24.80 | 0.220 | 0.059 | 11.21 | 0.613 | 0.284 | 0.022 | 0.234 | 1.64 | 0.126 | 0.043 | 0.010 | 0.005 | 0.154 | 0.046 |
| 2 | DH 2901 | 56.73 | 26.58 | 0.200 | 0.062 | 12.03 | 0.597 | 0.183 | 0.021 | 0.210 | 1.72 | 0.149 | 0.044 | 0.013 | 0.005 | 0.161 | 0.046 |
| 1 | VS F25/3 | 51.5 | 23.06 | 21.3 | . | . | 0.67 | . | . | . | . | . | . | 0.011 | 0.0056 | . | . |
| 2 | DH 5403 | 40.46 | 20.93 | 0.028 | 36.93 | 0.124 | 0.579 | 0.034 | 0.020 | . | 0.41 | . | 0.190 | 0.022 | . | 0.074 | 0.005 |
| 1 | VS F24/2 | 49.9 | . | . | 29.2 | . | 0.9 | 0.02 | . | . | . | . | . | 0.03 | 0.002 | . | . |
| 1 | NCS HC93635 | 27.36 | 26.23 | . | . | . | 38.51 | 1.90 | . | . | 0.18 | . | . | 0.072 | 0.015 | . | . |
| 1 | NCS HC93636 | 26.11 | . | . | . | . | 36.22 | 0.11 | . | . | 1.70 | . | . | 0.021 | 0.0071 | . | . |

CRM IRON IN FLAT SODA LIME GLASS

| Number | Fe | FeII | FeII as Fe2O3 | FeIII | Units |
|----------|--------|--------|---------------|----------|-------------------------|
| BAM S052 | 0.597 | 0.164 | (0.234) | (0.433) | 100 mm x 50 mm x 3.2 mm |
| BAM S051 | 0.0481 | 0.0158 | (0.0226) | (0.0323) | 100 mm x 50 mm x 5.9 mm |
| BAM S050 | 0.0084 | 0.0027 | (0.0038) | (0.0057) | 100 mm x 50 mm x 3.2 mm |

CRM MULTI-ELEMENT GLASS DISCS

listed in mg/kg each unit contains uncertified 72% SiO₂, 12% CaO, 14% Na₂O, and 2% Al₂O₃ each sample is 4 wafers ~13 mm Ø

| 3 mm | 1 mm | Ag | Au | B | Ba | Cd | Ce | Co | Cu | Dy | Er | Eu | Fe | Ga | Gd | K | La |
|---------|---------|-------|--------|--------|------|--------|------|--------|--------|------|------|--------|--------|--------|------|-------|---------|
| SRM 610 | SRM 611 | (254) | (25) | (351) | . | . | . | (390) | (444) | . | . | . | 458 | . | . | (461) | . |
| SRM 612 | SRM 613 | 22.0 | (5) | (32) | (41) | . | (39) | (35.5) | (37.7) | (35) | (39) | (36) | 51 | . | (39) | (64) | (36) |
| SRM 614 | SRM 615 | 0.42 | (0.5) | (1.30) | . | (0.55) | . | (0.73) | 1.37 | . | . | (0.99) | (13.3) | (1.3) | . | 30 | (0.83) |
| SRM 616 | SRM 617 | . | (0.18) | (0.20) | . | . | . | . | (0.80) | . | . | . | (11) | (0.23) | . | 29 | (0.034) |

| 3 mm | 1 mm | Mn | Nd | Ni | Pb | Rb | Sb | Sc | Sm | Sr | Th | Ti | Tl | U | Yb | Zn |
|---------|---------|--------|------|--------|-------|---------|---------|---------|------|-------|--------|--------|----------|----------|---------|----------|
| SRM 610 | SRM 611 | 485 | . | 458.7 | 426 | 425.7 | . | . | . | 515.5 | 457.2 | (437) | (61.8) | 461.5 | . | (433) |
| SRM 612 | SRM 613 | (39.6) | (36) | 38.8 | 38.57 | 31.4 | . | . | (39) | 78.4 | 37.79 | (50.1) | (15.7) | 37.38 | (42) | . |
| SRM 614 | SRM 615 | . | . | (0.95) | 2.32 | 0.855 | (1.06) | (0.59) | . | 45.8 | 0.748 | (3.1) | (0.269) | 0.823 | SRM 615 | SOLD OUT |
| SRM 616 | SRM 617 | . | . | . | 1.85 | (0.100) | (0.078) | (0.026) | . | 41.72 | 0.0252 | (2.5) | (0.0082) | (0.0721) | . | . |

CRM URANIUM IN GLASS

analysis listed in mg/kg 12 mm Ø x 5 mm

| Number | U |
|-----------|-----------|
| IRMM 540R | 15.0 last |

CRM GLASS SAND

T = Total

SRM 89: 45 g

other SRM: 75 g

all others: 100 g units

| Number | SiO ₂ | Al ₂ O ₃ | BaO | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | PbO | SO ₃ | TiO ₂ | ZrO ₂ | LOI |
|----------|------------------|--------------------------------|---------|--------|--------------------------------|--------------------------------|------------------|---------|---------|--------------------------------|-------------------|-------------------------------|-------------|-----------------|------------------|------------------|-------|
| BCS 531 | 99.74 | 0.0327 | 0.00112 | 0.0040 | . | 0.00636 | 0.0039 | 0.00132 | 0.00014 | . | . | 0.00082 | SrO:0.00017 | 0.0160 | . | . | . |
| UNS SPS | 99.32 | 0.248 | . | 0.029 | . | 0.037 | 0.058 | 0.0071 | . | . | 0.045 | . | . | 0.035 | . | . | 0.167 |
| BCS 516 | 98.73 | 0.513 | 0.0040 | 0.0243 | 0.0081 | 0.0596 | 0.127 | 0.0387 | . | 0.0012 | 0.0195 | (0.013) | 0.0127 | . | 0.175 | (0.075) | 0.24 |
| BCS 528 | 95.62 | 2.447 | 0.0298 | 0.237 | 0.0008 | 0.1111 | 0.875 | 0.0887 | . | . | 0.101 | (0.20) | 0.0006 | . | 0.0486 | (0.014) | 0.271 |
| SRM 1413 | 82.77 | 9.90 | 0.12 | 0.74 | . | 0.24 | 3.94 | 0.06 | . | . | 1.75 | . | . | 0.11 | . | . | . |
| SRM 89 * | 65.35 | 0.18 | 1.40 | 0.21 | . | 0.049 | 8.40 | 0.03 | 0.088 | . | 5.70 | 0.23 | 17.50 | 0.03 | 0.01 | 0.005 | 0.32 |
| SRM 81a | . | 0.66 | . | . | 0.0046 | 0.082 | . | . | . | . | . | . | . | . | 0.12 | 0.034 | . |
| SRM 165a | . | 0.059 | . | . | . | 0.012 | . | . | . | . | . | . | . | . | 0.011 | 0.006 | . |

* SRM 89 also contains As₂O₃: 0.03, As₂O₅: 0.36, Cl: 0.05

RM GRAVEL

typical analysis listed in mass %

100 g units

| Number | SiO ₂ | Al ₂ O ₃ | CO ₂ | CaO | Co ₃ O ₄ | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | S | TiO ₂ | -H ₂ O 900°C |
|---------|------------------|--------------------------------|-----------------|-------|--------------------------------|--------------------------------|--------------------------------|------------------|-------|--------------------------------|-------------------|-------------------------------|-------|------------------|-------------------------|
| DH 3610 | 98.80 | 0.234 | . | 0.008 | . | 0.030 | 0.419 | 0.014 | . | 0.009 | <0.003 | . | 0.009 | . | 0.153 |
| DH 3609 | 96.44 | 1.46 | 0.010 | 0.047 | 0.005 | 0.029 | 0.703 | 0.334 | 0.104 | 0.020 | 0.045 | 0.019 | . | 0.086 | 0.48 |

HARDGROVE GRINDABILITY INDEX

| Class | Set Number | HGI | HGI | HGI | HGI | Units |
|-------|-----------------|-----------------|-----------------|-----------------|-----------------|-------------------------------------|
| CRM | NCS AG82001d-4d | sample 1d: 35 | sample 2d: 56 | sample 3d: 74 | sample 4d: 107 | 250 g of each 1d - 4d last of stock |
| CRM | NCS AG82001g-4g | sample 1g: 35 | sample 2g: 54 | sample 3g: 72 | sample 4g: 119 | 250 g of each 1d - 4d |
| RM | ACIRS H7 | sample A : 30.0 | sample B : 45.9 | sample C : 63.5 | sample D : 85.5 | 1 kg of each A - D |
| RM | ACIRS H5 | sample A : 28 | sample B : 46 | sample C : 63 | sample D : 85 | 1 kg of each A - D last |

RM HARDGROVE GRINDABILITY INDEX

individually available in 1 kg units

| Number | HGI | Number | HGI | Number | HGI | Number | HGI | Number | HGI |
|--------------|-----|--------------|-----|--------------|-----|--------------|-----|--------------|-----|
| COCO HGI 007 | 66 | COCO HGI 021 | 57 | COCO HGI 006 | 54 | COCO HGI 018 | 51 | COCO HGI 022 | 47 |
| COCO HGI 005 | 64 | COCO HGI 019 | 57 | COCO HGI 026 | 54 | COCO HGI 020 | 51 | COCO HGI 023 | 46 |
| COCO HGI 013 | 64 | COCO HGI 004 | 56 | COCO HGI 014 | 53 | COCO HGI 016 | 50 | COCO HGI 010 | 31 |
| COCO HGI 008 | 60 | COCO HGI 003 | 55 | COCO HGI 024 | 52 | COCO HGI 015 | 48 | | |
| COCO HGI 009 | 59 | COCO HGI 011 | 55 | COCO HGI 012 | 51 | COCO HGI 027 | 48 | | |

CRM HARDNESS TEST BLOCKS

for NCS items, please indicate desired hardness when ordering

| Number | Scale | Available Range | Units (mm) |
|-----------------|--|-----------------|-----------------------------|
| NCS HBW | Brinell Hardness W | (8-650) | 100 x 80 x 16 |
| NCS HL | Leeb Hardness | (200-900) | 90 Ø x 55 |
| NCS HLG | Leeb Type G Hardness | (300-750) | 120 Ø x 70 |
| NCS HRA | Rockwell Hardness A | (20-88) | 60 x 40 x 10 |
| NCS HRB | Rockwell Hardness B | (20-100) | 60 x 40 x 10 |
| BS TRM-3 | Rockwell Hardness B | 86.3 | 300 x 300 x -2 17025 |
| NCS HRC | Rockwell Hardness C | (20-70) | 60 x 40 x 10 |
| NCS HR15N | Rockwell Superficial Hardness 15N | (70-94) | 60 x 40 x 10 |
| NCS HR30N | Rockwell Superficial Hardness 30N | (42-86) | 60 x 40 x 10 |
| NCS HR45N | Rockwell Superficial Hardness 45N | (20-77) | 60 x 40 x 10 |
| BS TRM-4 | Rockwell Superficial Hardness 15T | 71.9 | 300 x 300 x -2 17025 |
| NCS HR15T | Rockwell Superficial Hardness 15T | (67-93) | 60 x 40 x 10 |
| NCS HR30T | Rockwell Superficial Hardness 30T | (29-82) | 60 x 40 x 10 |
| NCS HR45T | Rockwell Superficial Hardness 45T | (1-72) | 60 x 40 x 10 |
| NCS HSD | Shore Hardness | (5-105) | 65 x 52 x 15 |
| NCS HV | Vickers Hardness | (5-1000) | 60 x 40 x 10 |
| NCS HVM | Vickers Microhardness | (5-1000) | 25 x 25 x 6 |

CRM INCINERATED WASTE

analysis listed in mg/kg

30 g powder

| Number | As | Ba | Be | Cd | Co | Cr | Cu | Hg | Mo | Ni | Pb | Sb | Se | Sn | Sr | V | Zn |
|------------|----|------|-----|------|----|-----|-----|-----|------|-----|--------|------|----|-------|-------|------|-------|
| BL 12-1-12 | 45 | 3600 | (8) | (60) | 23 | 731 | 375 | 7.8 | (10) | 198 | (1389) | (67) | 4 | (815) | (233) | (69) | 10450 |

informational analysis listed in mass %

| Number | Al ₂ O ₃ | CO ₂ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | TiO ₂ |
|------------|--------------------------------|-----------------|---------|--------------------------------|------------------|--------|--------|-------------------|-------------------------------|-----------------|------------------|------------------|
| BL 12-1-12 | (11.92) | (11.05) | (13.68) | (4.44) | (3.23) | (3.41) | (0.46) | (2.56) | (1.77) | (2.22) | (41.78) | (1.14) |

CRM IMPACT

approximate analysis

| Number | Energy | Uncertainty | Temperature | Units | Type |
|-----------|-------------|--------------|---------------|--------------------------------|-------------------------|
| SRM 2098 | 176 - 244 J | 8.8 - 12.2 J | 21 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| ERM-FA415 | 155.1 J | 4.6 J | 20 °C +/- 2' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| ERM-FA016 | 122.0 J | 3.6 J | 20 °C +/- 2' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| SRM 2097 | 101.9 J | 0.572 | -40 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| SRM 2096 | 88 - 136 J | 4.4 - 6.8 J | -40 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| ERM-FA015 | 79.8 J | 2.4 J | 20 °C +/- 2' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| ERM-FA013 | 28.1 J | 0.8 J | 20 °C +/- 2' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| SRM 2093 | 15.4 J | 0.125 | -40 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| SRM 2092 | 13 - 20 J | 1.4 J | -40 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 55 mm | CHARPY v-notch |
| SRM 2115 | 13 - 25 J | 1.4 J | 21 °C +/- 1' | 5 pcs of 10 mm x 10 mm x 75 mm | IZOD beam last of stock |

CRM LAYER THICKNESS

BCR: 2 sets of 4 Tantalum foils, 5 mm x 10 mm NMIJ: 13-15 mm squares

| Number | Material | Thickness | (+/-) | Layer 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
|------------|---|-----------|---------|-----------------|------|------|------|------|-----------------------|---|---|
| NMIJ 5202a | Si, SiO ₂ multi layer | n/a nm | 0.7 nm | (20.5) | 20.0 | 20.5 | 19.9 | 20.4 | surface oxide: (1.32) | | |
| NMIJ 5203a | GaAs, AlAs multi layer | n/a nm | 0.10 nm | (9.24) | 9.65 | 9.51 | 9.64 | 9.51 | 9.62 | . | . |
| NMIJ 5204b | GaAs, AlAs single layer | 3.26 nm | 0.41 nm | . | . | . | . | . | . | . | . |
| BCR 261T | Ta ₂ O ₅ single layer | 1.72 nm | 0.07 nm | 30 nm material | . | . | . | . | . | . | . |
| BCR 261T | Ta ₂ O ₅ single layer | 5.40 nm | 0.12 nm | 100 nm material | . | . | . | . | . | . | . |

CRM NANOSCALE LAYER THICKNESS

| Number | Certified Values | Informational Data | Units |
|----------|---|----------------------------------|----------------------|
| BAM L200 | 35 certified lengths from 3.5 - 4642 nm | 5 informational lengths 1 - 5 nm | block ~10 x 4 x 5 mm |

CRM LEAD PAINT FILMS

sold in SET/6 only, thin paint film on polyester sheets last of stock ~7cm wide and ~10 cm long

| Number | film, Pb in mg/kg | film, Pb in mg/kg | film, Pb in mg/kg | film, Pb in mg/kg | film, Pb in mg/kg | film, Pb in mg/kg | film, Pb in mg/kg |
|-----------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
| SRM 2579a | 2571 3.58 | 2572 1.527 | 2572 1.527 | 2573 1.040 | 2574 0.714 | 2575 0.307 | 2570 <0.001 |

RM ELECTROLYTIC MANGANESE

typical analysis

50 g units

| Number | Al | C | Co | Cr | Cu | Fe | Mn | Ni | P | S | Si | Zn | -H ₂ O@900°C |
|---------|----------|-------|--------|-------|--------|------|-------|--------|-------|--------|------|--------|-------------------------|
| DH 7701 | (0.0015) | 0.120 | 0.0012 | 0.411 | 0.0070 | 2.07 | 95.85 | 0.0068 | 0.056 | 0.0160 | 1.09 | 0.0011 | 0.019 |

CRM MANGANESE METAL POWDER

analysis listed in mass %

| Number | Mn | C | Fe | N | P | S | Se | Si | Units |
|-------------|-------|-------|-------|------|-------|-------|-------|-------|-------|
| NCS HC25655 | 97.43 | 0.080 | 1.81 | . | 0.018 | 0.016 | . | 0.28 | 50 g |
| NCS HC15604 | 95.52 | 0.154 | 3.37 | . | 0.032 | 0.019 | . | 0.92 | 100 g |
| NCS HC26615 | 91.56 | 0.007 | 0.039 | 7.84 | . | 0.031 | 0.049 | 0.009 | 50 g |

last

MELTING POINT

| Class | Number | Form | Melting point °C |
|-------|---------|--------------------|------------------|
| CRM | 501-951 | 6 inch nickel wire | 1455 |
| RM | 502-496 | 6 inch gold wire | 1063 |

CRM OXIDE

analysis listed in mg/kg except % which is mass %

100 g units

| Number | Notes | Ag | Al | As | B | Ba | Be | C | Ca | Cd | Ce | Cl | Co | Cr |
|-----------|---|----|-------|--------|------|-------|--------|--------|------|--------|--------|-------|------|-------|
| BAM RS 1 | SiO ₂ > 99.99% | . | 8.7 | <0.1 | . | . | . | . | 0.42 | <0.05 | . | . | . | 0.062 |
| BAM RS 2 | Al ₂ O ₃ = 99.76% | . | . | (<0.5) | (<5) | . | (<0.2) | . | 3.1 | (<0.5) | (<0.1) | (<10) | <1 | <1.5 |
| BAM RS 5 | NiO | <1 | (<15) | <0.2 | . | <1 | . | 14 | 2.2 | <0.2 | . | . | <2 | 16.1 |
| BAM RS 6A | MgO 100 - 350 μm | . | 46 | . | . | (<10) | . | (<50) | 994 | . | . | . | (<5) | 9.2 |
| BAM RS 6B | MgO 50 - 100 μm | . | 49 | . | . | (<20) | . | (<210) | 956 | . | . | . | (<5) | 8.1 |

continued

| Number | Cu | Fe | Ga | Ge | Hg | In | K | La | Li | Mg | Mn | Mo | Na | Ni | Pb |
|-----------|------|------|------|----|-------|--------|------|--------|------|--------|------|-------|-----|--------|-------|
| BAM RS 1 | <0.1 | 0.62 | . | <1 | <0.05 | . | 0.48 | . | 0.25 | <0.5 | <0.2 | . | <2 | <0.2 | <0.15 |
| BAM RS 2 | <2.5 | 3.3 | (<2) | . | . | (<0.5) | (<5) | (<0.3) | <1 | <3 | <1.5 | (<1) | <15 | <10 | . |
| BAM RS 5 | 1.53 | 41 | <0.5 | . | . | <1 | <2 | . | (<2) | <1 | <1 | <5 | <2 | 78.57% | <2 |
| BAM RS 6A | (<6) | 72 | . | . | . | . | . | . | . | 60.19% | 5.4 | (<10) | . | 3.9 | (<5) |
| BAM RS 6B | (<6) | 71 | . | . | . | . | . | . | . | 60.17% | 5.2 | (<10) | . | 3.3 | (<5) |

continued

| Number | S | Sb | Se | Si | Sn | Sr | Te | Ti | Tl | V | W | Zn | Zr |
|-----------|-----|--------|----|------|------|------|--------|------|--------|------|------|------|--------|
| BAM RS 1 | . | . | . | . | . | . | . | 1.3 | . | . | . | <1.3 | <0.1 |
| BAM RS 2 | . | . | . | <20 | (<1) | . | . | <2 | . | (<1) | . | <2 | 3.2 |
| BAM RS 5 | (4) | (<0.1) | <1 | (<5) | (<1) | (<1) | (<0.2) | (<2) | (<0.5) | <1 | (<1) | 3.4 | (<1) |
| BAM RS 6A | . | . | . | . | . | 2.0 | . | 1.3 | . | 8.4 | . | (<6) | (<20) |
| BAM RS 6B | . | . | . | . | . | 2.1 | . | 1.2 | . | 7.8 | . | (<6) | (<105) |

CRM IRON OXIDE

analysis listed in mass %

75 g units

| Number | Fe ₂ O ₃ | FeO | Al | C | Ca | Cr | Cu | K | Mg | Mn | Ni | S | Si | Other Impurities |
|----------|--------------------------------|--------|-------|---------|---------|--------|--------|---------|----------|-------|-------|--------|--------|------------------|
| VS P26/2 | 99.49 | (<0.1) | 0.026 | (0.005) | (0.005) | 0.0194 | 0.0090 | (0.001) | (<0.005) | 0.292 | 0.024 | (0.04) | 0.0110 | (0.1) |

CRM IRON OXIDE

analysis listed in mass %

analysis listed in mg/kg

100 g units

| Number | T.Fe | Cl | Mn | Al | Ca | Co | Cr | Cu | K | Mg | Mo | Na | Ni | P | Si | Sn | Ti | Zn |
|------------|-------|-------|-------|-----|----|----|-----|----|----|----|----|----|-----|----|----|----|----|----|
| ECRM 686-1 | 69.44 | 0.095 | 0.231 | 407 | 97 | 19 | 182 | 38 | 24 | 27 | 7 | 58 | 127 | 78 | 83 | 25 | 14 | 4 |

CRM NICKEL OXIDE

certified analysis listed in mass % except * which is mg/kg

25 g units

| Number | Al | Co | Cr | Cu | Fe | Mg | Mn | Si | Ti | Bi* | Pb* | Se* |
|---------|-------|-------|--------|-------|-------|-------|--------|-------|-------|------|-----|-----|
| SRM 673 | 0.001 | 0.016 | 0.0003 | 0.002 | 0.029 | 0.003 | 0.0037 | 0.006 | 0.003 | 0.06 | 3.5 | 0.2 |

continued informational analysis in mg/kg

Certified values show concentrations in nickel oxide. To convert values to the percent concentration in total metal present, multiply the values by 1.29.

| Number | Ag | As | Cd | Ga | Sb | Sn | Te | Tl | Zn | |
|---------|------|-----|------|------|------|------|-----|------|-----|---------------|
| SRM 673 | <0.1 | 0.4 | 0.05 | <0.1 | <0.5 | <0.5 | 0.4 | <0.1 | 1.7 | last of stock |

CRM SILICON OXIDE

analysis listed in mass % except

| Number | SiO ₂ | Al ₂ O ₃ | CaO | Fe | MgO | MnO | TiO ₂ | Units |
|-------------|------------------|--------------------------------|------|------|------|-------|------------------|-------|
| IRSID 608-1 | 60.39 | 9.94 | 8.70 | 4.00 | 1.34 | 0.057 | 0.714 | 100 g |

CRM TITANIUM DIBORIDE

analysis listed in mass %

powder 50 g

| Number | Ti | B | B ₂ O ₃ | Al | C | Ca | Cr | Fe | Mg | Mn | Mo | Ni | O | V | Zr | InsRes |
|----------|------|------|-------------------------------|--------|---------|--------|--------|-------|---------|---------|---------|--------|--------|---------|--------|--------|
| BAM S012 | 68.3 | 30.7 | 0.35 | 0.0012 | (0.169) | 0.0044 | 0.0097 | 0.064 | 0.00016 | 0.00038 | 0.00117 | 0.0023 | (0.48) | 0.00102 | 0.0121 | (0.22) |

CRM TITANIUM DIOXIDE SET

analysis listed in mass %

ONLY available in SET/8 x 20 g units

| Number | Cr | Cu | Fe | Mn | Mo | Ni | Si | Sn | V |
|-------------|---------|--------|---------|--------|--------|---------|---------|---------|--------|
| GSO 2158-81 | 0.0010 | . | . | . | . | . | 0.00054 | 0.00020 | . |
| GSO 2159-81 | 0.00035 | . | 0.00055 | . | . | 0.00046 | 0.0010 | . | . |
| GSO 2160-81 | 0.0013 | 0.0110 | 0.0010 | 0.090 | . | 0.0120 | 0.0015 | 0.0018 | 0.0014 |
| GSO 2161-81 | 0.0013 | 0.0024 | 0.0023 | 0.0010 | 0.0130 | 0.0088 | . | 0.0028 | . |
| GSO 2162-81 | 0.0023 | 0.0043 | 0.0180 | 0.0025 | 0.0048 | 0.0029 | 0.0130 | 0.0047 | 0.1800 |
| GSO 2163-81 | 0.038 | 0.032 | . | 0.0180 | . | 0.0280 | 0.0030 | . | 0.0016 |
| GSO 2164-81 | . | . | 0.0095 | . | 0.0110 | . | 0.0180 | . | . |
| GSO 2165-81 | . | 0.0023 | 0.0082 | 0.0040 | 0.0017 | 0.0014 | . | 0.035 | 0.0040 |

CRM TITANIUM DIOXIDE

| Number | TiO ₂ | Uncertainty | Units |
|----------|------------------|-------------|-----------|
| SRM 154c | 99.591 | +/- 0.062 | 90 g last |

CRM VANADIUM PENTOXIDE

analysis listed in mass %

NCS: 25-50 g units

SARM, VS: 100 g units

| Number | V ₂ O ₅ | V ₂ O ₄ | V | Al ₂ O ₃ | C | CaO | Fe | Fe ₂ O ₃ | K | K ₂ O | Na | Na ₂ O | P | S | Si | SiO ₂ | TiO ₂ | Others |
|--------------|-------------------------------|-------------------------------|-------|--------------------------------|------|------|-------|--------------------------------|---|------------------|----|-------------------|--------|----------|-------|------------------|------------------|--------------|
| NCS HC26612a | 98.99 | . | . | . | . | . | 0.078 | . | . | 0.12 | . | 0.58 | 0.022 | (0.001) | 0.080 | . | . | As: (0.0008) |
| NCS HC19611 | 98.80 | . | . | Cr:0.018 | . | . | 0.061 | . | . | 0.14 | . | 1.03 | 0.010 | 0.011 | 0.102 | . | . | As: (<0.001) |
| NCS HC19610 | 96.68 | . | . | Cr:0.099 | . | . | 0.43 | . | . | 0.18 | . | 0.96 | 0.007 | 0.014 | 0.40 | . | . | As: (<0.001) |
| SARM 38 | 95.52 | 3.07 | 55.84 | 0.14 | . | . | . | 0.119 | . | 0.600 | . | 0.22 | . | (0.0045) | . | 0.11 | . | MgO: 0.0037 |
| VS R30 | 94.3 | . | . | 0.007 | 0.88 | 0.51 | . | 0.053 | . | 0.032 | . | 0.0064 | 0.0072 | . | 0.43 | 0.21 | MnO: 2.58 | |

CRM ZINC OXIDE

| Number | Zn | CaO | Cd | Cl | F | Fe | MgO | Pb | S | SiO ₂ | Units |
|---------|------|------|--------|-------|-------|------|------|------|------|------------------|-------|
| IMN TC9 | 53.4 | 6.96 | 0.0049 | 0.033 | 0.055 | 5.64 | 3.50 | 3.77 | 0.52 | 5.47 | 220 g |

CRM PAPER

AVAILABLE IN SET/20 ONLY

includes software for data processing

5 pages per sample, 8.5 x 11" each

| Number | dry TAPPI analysis listed in mass % | | | | | | | Total | 400°C | 900°C | Base Weight |
|--------|-------------------------------------|--------|------------------|-------|-----------|--------------------------------|-------------------------------|-------|-------|-------|-------------|
| | CaCO ₃ | Kaolin | TiO ₂ | Talc | Muscovite | Al ₂ O ₃ | P ₂ O ₅ | | | | |
| A | 9.88 | 0.28 | 0.00 | 1.41 | 0.00 | . | . | 11.57 | 11.88 | 7.32 | 75 |
| B | 18.20 | 0.28 | 0.00 | 0.00 | 0.00 | . | . | 18.48 | 18.53 | 10.65 | 75 |
| C | 12.53 | 0.56 | 0.00 | 0.60 | 0.00 | . | . | 13.69 | 13.58 | 8.11 | 75 |
| D | 18.29 | 0.00 | 0.00 | 0.00 | 0.00 | . | . | 18.29 | 18.76 | 10.51 | 75 |
| E | 9.45 | 0.00 | 0.00 | 0.00 | 0.00 | . | . | 9.45 | 10.14 | 5.78 | 75 |
| F | 11.22 | 0.00 | 0.39 | 0.60 | 0.00 | . | . | 12.21 | 12.34 | 7.49 | 75 |
| G | 12.26 | 0.18 | 0.00 | 0.41 | 0.00 | . | . | 12.85 | 13.08 | 7.56 | 75 |
| H | 11.19 | 1.34 | 0.00 | 0.38 | 0.00 | . | . | 12.91 | 11.98 | 8.01 | 75 |
| I | 18.94 | 0.00 | 0.00 | 0.28 | 0.00 | . | . | 19.22 | 19.71 | 11.11 | 80 |
| J | 14.79 | 0.51 | 0.09 | 1.48 | 0.00 | . | . | 16.87 | 17.11 | 10.65 | 75 |
| K | 14.12 | 2.10 | 0.28 | 1.88 | 0.00 | . | . | 18.38 | 18.30 | 12.17 | 75 |
| L | 0.00 | 7.54 | 1.75 | 0.00 | 0.00 | . | . | 9.29 | 8.81 | 8.38 | 75 |
| M | 0.16 | 10.91 | 0.18 | 0.00 | 0.00 | . | . | 11.25 | 11.16 | 10.12 | 75 |
| N | 1.74 | 0.00 | 1.51 | 10.74 | 0.00 | . | . | 13.99 | 14.70 | 13.28 | 75 |
| O | 1.86 | 12.69 | 0.00 | 0.47 | 7.57 | . | . | 22.59 | 22.99 | 20.34 | 80 |
| P | 25.61 | 0.35 | 0.00 | 0.00 | 0.00 | . | . | 25.96 | 26.93 | 15.61 | 105 |
| Q | 0.00 | 0.30 | 38.60 | 0.00 | 0.00 | 2.70 | 1.87 | 43.47 | 43.39 | 43.13 | 85 |
| R | 0.13 | 19.02 | 0.25 | 0.00 | 0.65 | . | . | 20.05 | 20.21 | 17.56 | 45 |
| S | 0.14 | 32.04 | 0.42 | 0.00 | 1.08 | . | . | 33.68 | 33.57 | 29.43 | 60 |
| BLANK | 0.00 | 0.00 | 0.00 | 0.00 | 0.00 | . | . | 0.02 | 0.02 | 0.01 | 75 |

CRM PARTICLE SIZE and MASS VOLUME in ALUMINA

| Number | Permeametry | BET Absorption | Obligatory Porosity | Size Range | Median Size |
|--------|---------------------------|---------------------------|---------------------|--------------|-------------|
| TL AA | 2,300 cm ² /g | 5,000 cm ² /g | 0.57 | 1-64 Ø µm | 12.7 Ø µm |
| TL AB | 10,300 cm ² /g | 31,000 cm ² /g | 0.67 | 1-31.50 Ø µm | 2.1 Ø µm |

CRM PARTICLE SIZE analysis listed in μm

| Number | Weight Percentile | Certified Value | Uncertainty | Units |
|-----------|-------------------|-----------------|-------------|---------------------------|
| SRM 1018b | 1.2 - 98.1 | 220 - 750 | 0.4 - 2.0 | 87 g of 0.2 g glass beads |

CRM PARTICLE SIZE

| Number | Average Diameter, μm | Uncertainty, μm | Material | Units |
|----------|---------------------------------|----------------------------|---------------------|-------|
| SRM 1691 | 0.269 | ± 0.007 | Polystyrene Spheres | 5 mL |

CRM PARTICLE SIZE

| Number | Quartz Form | Certified Property | Size Range in Microns | Unit Size |
|---------|-------------|--------------------|-----------------------|-----------|
| BCR 066 | Powder | Stokes' diameter | 0.35 - 3.50 | 10 g |
| BCR 070 | Powder | Stokes' diameter | 1.2 - 20 | 10 g |
| BCR 067 | Powder | Stokes' diameter | 2.4 - 32 | 10 g |
| BCR 069 | Powder | Stokes' diameter | 14 - 90 | 10 g |
| BCR 130 | Powder | Volume diameter | 50 - 220 | 50 g |
| BCR 068 | Sand | Volume diameter | 160 - 630 | 100 g |
| BCR 131 | Powder | Volume diameter | 480 - 1800 | 200 g |
| BCR 132 | Gravel | Volume diameter | 1400 - 5000 | 700 g |

CRM PARTICLE SIZE

| Number | Percentage of Particles Under 20 Microns | Standard Deviation | Uncertainty @ 95% CL | Units |
|-----------|--|--------------------|----------------------|-------|
| ASCRM 026 | 1.0 | ± 0.1 | ± 0.2 | 210 g |

CRM PARTICLE DENSITY, SURFACE AREA, AND SIZE DISTRIBUTION

| Number | Particle Density Picnometer Method | Blaine Area With EN 196-6 | Particle Size by Laser Diffraction ISO 13320-1 | Air Jet Seiving Alpine Test NF X11-640 | Units |
|---------|---------------------------------------|------------------------------|---|---|----------|
| TL 1BGa | 3.11 g/cm ³ | 3396 cm ³ /g | 11.4% @ 2.0 μm - 99.9% @ 160 μm | 71.4% @ 31.5 μm - 100% @ 160 μm | 20 x 5 g |

CRM pH STANDARDS

| Number | PH | Accreditation | Units |
|--------|-------|---------------|---|
| TL pH9 | 9.180 | ISO GUIDE 34 | 10 x 25 mL or as a SET/15 bottles = 5 x pH4, 5 x pH7, 5 x pH9 |
| TL pH7 | 6.865 | ISO GUIDE 34 | 10 x 25 mL or as a SET/15 bottles = 5 x pH4, 5 x pH7, 5 x pH9 |
| TL pH4 | 4.005 | ISO GUIDE 34 | 10 x 25 mL or as a SET/15 bottles = 5 x pH4, 5 x pH7, 5 x pH9 |

RM PLASTER analysis listed in mass % 100 g units

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ | SrO | TiO ₂ | LOI |
|----------|--------------------------------|------|--------------------------------|------------------|------|-------------------|-------------------------------|-----------------|------------------|------|------------------|-----|
| BCS 202A | 0.33 | 37.4 | 0.10 | 0.10 | 0.39 | <0.03 | <0.01 | 53 | 1.38 | 0.33 | 0.03 | 7.0 |

RM PLASTIC analysis listed in mass % 50 g pellets

| Number | Al | Ca | Cl | Cr | F | Fe | Mg | Na | P | S | Si | Ti | Zn |
|------------|--------|--------|----------|--------|---------|--------|--------|--------|--------|----------|----------|--------|--------|
| JSM P702-1 | 0.0012 | 0.0013 | (0.0017) | 0.0012 | . | 0.0015 | 0.0015 | 0.0012 | 0.0011 | (0.0014) | (0.0008) | 0.0009 | 0.0010 |
| JSM P703-1 | 0.022 | 0.023 | (0.018) | 0.020 | (0.018) | 0.021 | 0.029 | 0.024 | 0.017 | (0.021) | 0.018 | 0.017 | 0.020 |

CRM POROUS MATERIAL

| Number | Description | Units | Specific Pore Volume | Median Pore Diameter | Density |
|-----------|-----------------|-----------------------|------------------------|----------------------|----------------------------|
| BAM P 128 | Alumina Ceramic | 6 Cylinders 7 g total | 220 mm ³ /g | 27.6 μm | (3.6405) g/cm ³ |

CRM POROUS MATERIALS and SURFACE AREA

| Number | Description | Units | A_{BET} (m ² /g) BET Specific Surface Area | V_p (cm ³ /g) Specific Pore Volume | D_1 (nm) Hydraulic Pore Diameter | D_2 (nm) Most Frequent Pore Diameter | D_3 (nm) Most Frequent Pore Diameter | (nm) Median Pore Width |
|-----------|-------------------|-------|--|--|---------------------------------------|---|---|---------------------------|
| BAM P 109 | Activated Carbon | 10g | 1396 | . | . | . | . | . |
| BAM P 105 | Glass Material | 10g | 198.5 | 0.2327 | 4.69 | 4.38 | 5.80 | . |
| BAM P 106 | Titanium Dioxide | 15g | 96.6 | 0.2341 | 9.69 | 8.2 | 11.5 | . |
| BAM FD107 | Faujasite Zeolite | 10g | . | 0.217 cm ³ /g ⁻¹ | . | . | . | 0.86 |

| Number | Description | Units | (nm) Mean Pore Radius | (nm) Most Frequent Pore Radius | (cm ² /g) Specific Surface Area | (mm ³ /g) Pore Volume 100 Mpa | (mm ³ /g) Pore Volume 195 Mpa | (mm ³ /g) Pore Volume 200 Mpa | (mm ³ /g) Pore Volume 395 Mpa |
|------------|--------------------------------------|-------|--------------------------|-----------------------------------|---|---|---|---|---|
| BAM PM 101 | SiO ₂ | 10g | . | . | 0.177 | . | . | . | . |
| BAM PM 102 | Alpha-Al ₂ O ₃ | 10g | . | . | 5.41 | . | . | . | . |
| BAM FD 120 | Alpha-Al ₂ O ₃ | 10g | 228.0 | 232.2 | . | 545.0 | 546.7 | 546.8 | 548.1 |
| BAM FD 121 | Porous glass | 12g | 15.1 | 15.3 | . | 621.8 | 621.9 | 621.9 | 624.6 |
| BAM FD 122 | Porous glass | 15g | 139.0 | 140.2 | . | 919.7 | 922.5 | 922.6 | 924.4 |

ROHS/WEEE SAMPLES

= class, where 1 = CRM and 2 = RM analysis listed in mass %

| # | Type | Units | Number | As | Br | Cd | Cr | Hg | Pb | Se | Sold As |
|---|--------------------------|------------------|---------------|----------|--------|------------|----------|-----------|----------|------------|----------------------------------|
| 1 | ABS resin | pellets 25 g | NMIJ 8112a | . | . | 0.000938 | 0.009447 | 0.009410 | 0.009498 | . | individually |
| 1 | ABS resin | granules 100g | BAM H010 gran | . | 0.0240 | 0.0093 | 0.0470 | (0.0415) | 0.0479 | . | individually |
| 1 | ABS resin | 40 mm Ø x 1 mm | BAM H010 1mm | . | 0.0240 | 0.0093 | 0.0470 | (0.0415) | 0.0479 | . | individually |
| 1 | ABS resin | 40 mm Ø x 2 mm | BAM H010 2mm | . | 0.0240 | 0.0093 | 0.0470 | (0.0415) | 0.0479 | . | individually |
| 1 | ABS resin | 40 mm Ø x 6 mm | BAM H010 6mm | . | 0.0240 | 0.0093 | 0.0470 | (0.0415) | 0.0479 | . | individually |
| 1 | ABS resin | 40 mm Ø x 1-6 mm | BAM H010 set | . | 0.0240 | 0.0093 | 0.0470 | (0.0415) | 0.0479 | . | set of above 3 discs |
| 1 | ABS resin | 30 mm Ø x 2 mm | NMIJ 8105a | . | . | 0.001070 | 0.002751 | . | 0.010828 | . | individually |
| 1 | plastic | 40 mm Ø x 4 mm | JSAC 0631 | . | . | 0.00225 | 0.00258 | 0.00197 | 0.00245 | . | set only |
| 1 | plastic | 40 mm Ø x 4 mm | JSAC 0632 | . | . | 0.00461 | 0.00933 | 0.00594 | 0.00929 | . | |
| 1 | polyester | chips 50 g | JSAC 0602-3 | . | . | 0.00506 | 0.01125 | 0.00121 | 0.01121 | . | individually |
| 1 | polyester | chips 50 g | JSAC 0601-3 | . | . | 0.00050 | 0.0113 | 0.000130 | 0.00121 | . | individually |
| 1 | low density polyethylene | pellets 100 g | ERM-EC681m | 0.00170 | 0.143 | 0.0146 | 0.00451 | 0.00099 | 0.00697 | | also Cl S Sb Sn Zn, individually |
| 1 | low density polyethylene | pellets 100 g | ERM-EC680m | 0.00047 | 0.0181 | 0.00208 | 0.00096 | 0.000256 | 0.00113 | | also Cl S Sb Sn Zn, individually |
| 1 | low density polyethylene | pellets 100 g | ERM EC680k | 0.00041 | 0.0096 | 0.00196 | 0.00202 | 0.000464 | 0.00136 | | also Cl and S, individually |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 1 | . | 0 | 0 | 0 | 0 | 0 | . | set only |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 2 | . | 0.1300 | 0.0029 | 0.0053 | 0.0053 | 0.0058 | . | |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 3 | . | 0.0670 | 0.0059 | 0.0110 | 0.0110 | 0.0120 | . | |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 4 | . | 0.0310 | 0.0089 | 0.0320 | 0.1200 | 0.0360 | . | |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 5 | . | 0.0110 | 0.0120 | 0.0660 | 0.0620 | 0.0730 | . | |
| 2 | low density polyethylene | 48 mm Ø x 5 mm | SCAS PE-5E6 6 | . | 0.0053 | 0.0350 | 0.1300 | 0.0290 | 0.1400 | . | |
| 1 | polypropylene resin | 30 mm Ø x 2 mm | NMIJ 8136a | . | . | 0.00937 | 0.08906 | 0.0952 | 0.09435 | . | individually |
| 1 | PVC resin | pellets 25 g | NMIJ 8123a | . | . | 0.009562 | 0.09490 | 0.09370 | 0.09655 | . | individually |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0611 | . | . | 0.00000 | 0.00000 | . | 0.00000 | . | set only |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0612 | . | . | 0.00086 | 0.00243 | . | 0.00242 | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0613 | . | . | 0.00219 | 0.00488 | . | 0.00485 | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0614 | . | . | 0.00430 | 0.00966 | . | 0.00959 | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0615 | . | . | 0.00866 | 0.01941 | . | 0.01929 | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0621 | . | . | . | . | (<0.0001) | . | . | set only |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0622 | . | . | . | . | 0.00100 | . | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0623 | . | . | . | . | 0.00490 | . | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0624 | . | . | . | . | 0.01211 | . | . | |
| 1 | polyvinyl chloride | 40 mm Ø x 4 mm | JSAC 0625 | . | . | . | . | 0.0244 | . | . | |
| 1 | soil | powder 50 g | JSAC 0403 | 0.0199 | . | 0.0183 | 0.0257 | 0.00111 | 0.0224 | . | individually |
| 1 | soil | powder 25 g | JSAC 0466 | 0.01093 | . | 0.01199 | 0.1483 | 0.01135 | 0.1214 | 0.1175 | set only |
| 1 | soil | powder 25 g | JSAC 0465 | 0.0550 | . | 0.06074 | 0.0738 | 0.00578 | 0.6124 | 0.0587 | |
| 1 | soil | powder 25 g | JSAC 0464 | 0.02711 | . | 0.03010 | 0.0499 | 0.00286 | 0.03027 | 0.02919 | |
| 1 | soil | powder 25 g | JSAC 0463 | 0.01376 | . | 0.01468 | 0.0244 | 0.001476 | 0.01516 | 0.01415 | |
| 1 | soil | powder 25 g | JSAC 0462 | 0.00715 | . | 0.00742 | 0.01496 | 0.000727 | 0.00737 | 0.00716 | |
| 1 | soil | powder 25 g | JSAC 0461 | 0.002153 | . | (0.000030) | 0.00972 | 0.0000075 | 0.00244 | (0.000044) | |
| 1 | zinc | 50 mm Ø x 20 mm | 41X ZSC6A | . | . | 0.215 | <0.0002 | 0.029 | 0.0077 | . | individually |
| 1 | zinc | 50 mm Ø x 20 mm | 41X ZSC3A | . | . | 0.119 | 0.0148 | 0.0021 | 0.0273 | . | |
| 1 | zinc | 50 mm Ø x 20 mm | 41X ZSC1A | . | . | 0.0288 | 0.0039 | 0.026 | 0.06 | . | |
| 1 | zinc | 50 mm Ø x 20 mm | 41X ZSC4B | . | . | 0.0131 | 0.0299 | 0.050 | 0.156 | . | |
| 1 | zinc | 50 mm Ø x 20 mm | 41X ZSC2A | . | . | 0.0016 | 0.0036 | 0.0053 | 0.111 | . | |
| # | Type | Units | Number | As | Br | Cd | Cr | Hg | Pb | Se | Sold As |

REFRACTORIES

= class, where 1 = CRM and 2 = RM

T = Total

IPT: 80 g

VS K6: 75 g

VS K10: 125 g

all others: 100 g

| # | Number | SiO ₂ | Al ₂ O ₃ | C | CO ₂ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | S | SO ₃ | TiO ₂ | ZrO ₂ | LOI |
|---|------------|------------------|--------------------------------|---------|-----------------|--------|--------------------------------|------------------|-------|-------|--------------------------------|-------------------|-------------------------------|-------|-----------------|------------------|------------------|--------|
| 1 | IPT 63 | 96.28 | 0.48 | . | . | 2.21 | 0.52 | 0.043 | 0.18 | 0.008 | . | 0.013 | 0.013 | . | . | 0.030 | (0.002) | 0.17 |
| 2 | DH 2612 | 40.80 | 36.45 | 0.437 | 0.54 | 1.80 | 3.10 | 0.759 | 13.13 | 0.125 | . | 0.242 | 0.279 | 0.034 | . | 1.25 | 0.163 | 0.75 |
| 2 | DH 2613 | 25.83 | 42.78 | 1.779 | 0.53 | 2.31 | 2.57 | 0.404 | 21.03 | . | 0.122 | 0.118 | 0.122 | 0.066 | . | 1.199 | . | . |
| 2 | DH 2609 | 23.41 | 63.82 | 0.739 T | 0.170 | 2.25 | 1.75 | 0.526 | 4.17 | 0.282 | . | 0.220 | 0.339 | . | 0.121 | 1.27 | 0.097 | . |
| 1 | FLX CRM112 | 12.16 | 79.81 | . | . | 0.147 | 0.326 | 0.090 | 0.755 | . | 0.024 | 0.267 | 0.074 | . | (0.04) | 0.273 | 5.95 | (5.42) |
| 1 | VS K6/4 | 2.12 | 0.54 | . | . | 2.95 | 2.26 | . | 92.4 | . | . | . | . | . | . | . | . | . |
| 1 | VS K10/3 | (0.2) | 97 | (0.05) | . | (0.03) | 1.82 | (0.03) | . | . | . | (0.5) | . | . | . | 0.35 | . | . |

continued

| Number | Co ₃ O ₄ | CuO | Cr ₂ O ₃ | HfO ₂ | Li ₂ O | NiO | V ₂ O ₅ | WO ₃ | Units |
|------------|--------------------------------|-------|--------------------------------|------------------|-------------------|-------|-------------------------------|-----------------|-------|
| IPT 63 | . | . | . | . | (0.0005) | . | . | . | 80 g |
| DH 2612 | . | . | 0.385 | . | . | 0.032 | 0.027 | . | 100 g |
| DH 2613 | . | 0.004 | 0.140 | . | . | . | 0.020 | . | 100 g |
| DH 2609 | . | . | . | . | . | . | . | . | 100 g |
| FLX CRM112 | <0.01 | . | 0.017 | 0.099 | . | <0.01 | . | 0.041 | 80 g |
| VS K6/4 | . | . | . | . | . | . | . | . | 100 g |
| VS K10/3 | . | . | . | . | . | . | . | . | 125 g |

CRM ALUMINA-MAGNESIA REFRACTORY SET

| Number | SOLD IN SET/10 ONLY | | | | | | | | | certified values | | | | 20 g units | | | | informational values | | | |
|----------|--------------------------------|-------|------|--------------------------------|------------------|-------------------|-------------------------------|------------------|------------------|--------------------------------|------|------------------|------|------------|--|--|--|----------------------|--|--|--|
| | Al ₂ O ₃ | MgO | CaO | Fe ₂ O ₃ | K ₂ O | Na ₂ O | P ₂ O ₅ | SiO ₂ | TiO ₂ | Cr ₂ O ₃ | MnO | ZrO ₂ | LOI | | | | | | | | |
| JRRM 801 | 93.49 | 3.26 | 0.14 | 2.00 | 0.01 | 0.19 | 0.00 | 0.35 | 0.21 | 0.00 | 0.00 | 0.00 | 0.14 | | | | | | | | |
| JRRM 802 | 84.25 | 6.13 | 2.00 | 1.03 | 0.46 | 0.15 | 0.95 | 3.32 | 1.48 | 0.00 | 0.00 | 0.00 | 0.06 | | | | | | | | |
| JRRM 803 | 74.23 | 16.20 | 0.57 | 4.90 | 0.00 | 0.86 | 0.01 | 0.58 | 2.51 | 0.00 | 0.00 | 0.00 | 0.36 | | | | | | | | |
| JRRM 804 | 64.66 | 20.84 | 4.76 | 4.02 | 0.04 | 0.08 | 0.11 | 5.17 | 0.13 | 0.01 | 0.02 | 0.00 | 0.01 | | | | | | | | |
| JRRM 805 | 58.03 | 36.04 | 0.28 | 0.73 | 0.01 | 0.54 | 0.68 | 2.49 | 1.05 | 0.00 | 0.00 | 0.00 | 0.17 | | | | | | | | |
| JRRM 806 | 48.85 | 49.43 | 0.97 | 0.16 | 0.00 | 0.04 | 0.04 | 0.51 | 0.00 | 0.00 | 0.02 | 0.00 | 0.21 | | | | | | | | |
| JRRM 807 | 39.96 | 55.07 | 2.75 | 0.32 | 0.15 | 0.32 | 0.53 | 0.58 | 0.19 | 0.00 | 0.00 | 0.00 | 0.57 | | | | | | | | |
| JRRM 808 | 28.68 | 67.01 | 0.99 | 0.56 | 0.69 | 0.40 | 0.22 | 0.79 | 0.71 | 0.00 | 0.01 | 0.00 | 0.84 | | | | | | | | |
| JRRM 809 | 19.86 | 70.11 | 4.47 | 0.11 | 0.98 | 0.04 | 1.06 | 0.36 | 2.88 | 0.00 | 0.00 | 0.00 | 0.48 | | | | | | | | |
| JRRM 810 | 10.08 | 78.96 | 0.18 | 3.11 | 0.16 | 0.75 | 0.51 | 4.21 | 1.91 | 0.00 | 0.01 | 0.00 | 0.22 | | | | | | | | |

CRM ALUMINA-ZIRCONIA-SILICA REFRACTORY SET

| Number | SOLD IN SET/10 ONLY | | | | | | | | | | certified values | | | 20 g units | | | informational values | | |
|----------|--------------------------------|------------------|------------------|------|--------------------------------|--------------------------------|------------------|------------------|------|-------------------|------------------|------|-------------------------------|------------|--|--|----------------------|--|--|
| | Al ₂ O ₃ | ZrO ₂ | SiO ₂ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | HfO ₂ | K ₂ O | MgO | Na ₂ O | TiO ₂ | MnO | P ₂ O ₅ | LOI | | | | | |
| JRRM 710 | 82.29 | 2.96 | 5.62 | 0.22 | 1.02 | 1.15 | 1.51 | 0.63 | 0.04 | 1.41 | 3.00 | 0.00 | 0.04 | 0.09 | | | | | |
| JRRM 708 | 79.52 | 12.84 | 0.54 | 1.17 | 0.29 | 0.80 | 1.03 | 0.74 | 1.64 | 0.08 | 1.02 | 0.00 | 0.00 | 0.13 | | | | | |
| JRRM 705 | 64.14 | 27.96 | 1.99 | 0.19 | 2.01 | 0.14 | 0.48 | 0.01 | 0.46 | 0.30 | 2.02 | 0.00 | 0.01 | 0.16 | | | | | |
| JRRM 707 | 55.78 | 18.16 | 21.17 | 1.08 | 0.18 | 1.81 | 0.36 | 0.15 | 0.84 | 0.19 | 0.28 | 0.00 | 0.05 | 0.01 | | | | | |
| JRRM 709 | 50.35 | 8.32 | 34.38 | 0.52 | 2.91 | 0.47 | 0.18 | 0.21 | 1.20 | 1.03 | 0.09 | 0.00 | 0.00 | 0.20 | | | | | |
| JRRM 703 | 46.34 | 37.35 | 14.64 | 0.03 | 0.00 | 0.05 | 0.72 | 0.00 | 0.01 | 0.53 | 0.07 | 0.00 | 0.03 | 0.09 | | | | | |
| JRRM 702 | 38.14 | 42.54 | 9.99 | 1.55 | 0.11 | 0.37 | 2.08 | 0.57 | 1.97 | 2.02 | 0.21 | 0.00 | 0.02 | 0.18 | | | | | |
| JRRM 706 | 25.95 | 22.72 | 39.33 | 1.58 | 0.01 | 0.13 | 1.19 | 0.95 | 0.15 | 3.49 | 3.77 | 0.00 | 0.01 | 0.72 | | | | | |
| JRRM 704 | 19.58 | 33.46 | 42.61 | 0.15 | 0.51 | 0.55 | 0.68 | 1.40 | 0.51 | 0.22 | 1.02 | 0.08 | 0.13 | 0.07 | | | | | |
| JRRM 701 | 10.09 | 48.06 | 28.44 | 2.07 | 1.01 | 2.00 | 0.85 | 0.02 | 0.47 | 1.84 | 4.96 | 0.00 | 0.02 | 0.09 | | | | | |

CRM BURNT REFRACTORIES

IPT: 80 g units SRM: 75 g units

| Number | Al ₂ O ₃ | SiO ₂ | CaO | Fe ₂ O ₃ | K ₂ O | Li ₂ O | MgO | Na ₂ O | P ₂ O ₅ | SrO | TiO ₂ | ZrO ₂ | LOI |
|---------|--------------------------------|------------------|------|--------------------------------|------------------|-------------------|------|-------------------|-------------------------------|-------|------------------|------------------|--------|
| SRM 78a | 71.7 | 19.4 | 0.11 | 1.2 | 1.22 | 0.12 | 0.70 | 0.078 | 1.3 | 0.25 | 3.22 | . | (0.42) |
| IPT 57 | 71.5 | 24.3 | 0.05 | 1.25 | 0.83 | 0.008 | 0.13 | 0.35 | 0.054 | 0.009 | 1.19 | 0.20 | 0.20 |
| SRM 77a | 60.2 | 35.0 | 0.05 | 1.00 | 0.090 | 0.025 | 0.38 | 0.037 | 0.092 | 0.009 | 2.66 | . | (0.22) |
| IPT 51 | 40.3 | 55.0 | 0.06 | 1.19 | 0.69 | 0.018 | 0.20 | 0.09 | 0.09 | . | 2.19 | 0.070 | 0.16 |
| SRM 76a | 38.7 | 54.9 | 0.22 | 1.60 | 1.33 | 0.042 | 0.52 | 0.07 | 0.120 | 0.037 | 2.03 | . | (0.34) |

CRM CHROME-MAGNESIA REFRACTORARY SET

SOLD IN SET/12 ONLY

certified values

informational values

20 g units

| Number | MgO | Cr ₂ O ₃ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | MnO | SiO ₂ | TiO ₂ | NiO | P ₂ O ₅ | V ₂ O ₅ | ZnO | LOI |
|----------|-------|--------------------------------|--------------------------------|------|--------------------------------|------|------------------|------------------|------|-------------------------------|-------------------------------|------|-------|
| JRRM 501 | 87.60 | 2.82 | 2.92 | 0.92 | 4.80 | 0.02 | 0.92 | 0.00 | 0.01 | 0.03 | 0.01 | 0.00 | 0.13 |
| JRRM 502 | 76.28 | 7.49 | 11.98 | 0.20 | 1.02 | 0.01 | 3.11 | 0.01 | 0.02 | 0.02 | 0.02 | 0.00 | 0.06 |
| JRRM 503 | 63.11 | 13.60 | 7.14 | 3.81 | 3.00 | 0.03 | 9.09 | 0.04 | 0.03 | 0.03 | 0.03 | 0.01 | 0.11 |
| JRRM 504 | 54.85 | 18.35 | 17.56 | 2.60 | 4.11 | 0.01 | 2.18 | 0.01 | 0.01 | 0.03 | 0.01 | 0.01 | 0.12 |
| JRRM 505 | 50.14 | 21.74 | 7.76 | 0.49 | 17.76 | 0.10 | 1.82 | 0.11 | 0.07 | 0.02 | 0.07 | 0.02 | 0.08 |
| JRRM 506 | 46.65 | 28.19 | 14.69 | 0.46 | 7.49 | 0.07 | 2.16 | 0.13 | 0.09 | 0.01 | 0.08 | 0.01 | 0.07 |
| JRRM 508 | 30.86 | 38.18 | 3.98 | 1.03 | 22.70 | 0.00 | 3.08 | 0.01 | 0.01 | 0.01 | 0.00 | 0.00 | 0.05 |
| JRRM 512 | 24.81 | 4.98 | 29.25 | 4.06 | 26.01 | 0.02 | 10.57 | 0.04 | 0.01 | 0.01 | 0.01 | 0.01 | 0.02 |
| JRRM 507 | 22.36 | 32.03 | 25.02 | 1.61 | 12.98 | 0.11 | 5.69 | 0.16 | 0.20 | 0.01 | 0.13 | 0.03 | -0.11 |
| JRRM 509 | 20.45 | 42.57 | 20.28 | 2.86 | 10.15 | 0.08 | 1.96 | 1.20 | 0.04 | 0.01 | 0.11 | 0.03 | 0.13 |
| JRRM 510 | 16.86 | 50.38 | 12.21 | 0.29 | 14.99 | 0.17 | 4.91 | 0.13 | 0.19 | 0.01 | 0.11 | 0.04 | -0.25 |
| JRRM 511 | 10.62 | 52.51 | 6.68 | 0.07 | 27.22 | 0.12 | 2.90 | 0.10 | 0.10 | 0.00 | 0.05 | 0.05 | -0.48 |

CRM FIRECLAY REFRACTORARY SET

SOLD IN SET/10 ONLY

20 g units

| Number | SiO ₂ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | TiO ₂ |
|-----------|------------------|--------------------------------|------|--------------------------------|------------------|------|------|-------------------|------------------|
| JRRM 101 | 88.57 | 8.10 | 1.06 | 0.31 | 0.16 | 0.21 | 0.11 | 1.01 | 0.30 |
| JRRM 102 | 80.47 | 13.79 | 0.04 | 3.97 | 0.14 | 0.67 | 0.01 | 0.30 | 0.45 |
| JRRM 103 | 80.32 | 18.07 | 0.07 | 0.40 | 0.35 | 0.01 | 0.00 | 0.12 | 0.37 |
| JRRM 104 | 67.35 | 22.52 | 0.25 | 3.24 | 3.04 | 0.07 | 0.01 | 0.30 | 2.94 |
| JRRM 105a | 69.17 | 25.35 | 0.40 | 0.76 | 0.81 | 0.22 | 0.11 | 0.65 | 2.24 |
| JRRM 106 | 63.61 | 29.91 | 0.14 | 1.92 | 1.81 | 0.98 | 0.02 | 0.59 | 0.67 |
| JRRM 107 | 55.32 | 37.08 | 0.71 | 2.20 | 2.57 | 0.49 | 0.01 | 0.21 | 1.15 |
| JRRM 108 | 55.31 | 40.08 | 0.27 | 1.54 | 0.80 | 0.27 | 0.02 | 0.20 | 1.05 |
| JRRM 109 | 54.23 | 41.24 | 0.14 | 0.89 | 0.79 | 0.12 | 0.01 | 0.30 | 1.96 |
| JRRM 110 | 49.54 | 46.68 | 0.10 | 0.84 | 0.34 | 0.16 | 0.01 | 0.08 | 1.66 |

CRM FIRECLAY REFRACTORARY SET

SOLD IN SET/15 ONLY

20 g units

| Number | SiO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | TiO ₂ | ZrO ₂ | LOI |
|----------|------------------|--------------------------------|------|--------------------------------|--------------------------------|------------------|------|------|-------------------|-------------------------------|------------------|------------------|--------|
| JRRM 121 | 86.3 | 6.07 | 1.96 | 0.01 | 0.40 | 0.23 | 0.12 | 0.02 | 3.20 | 0.32 | 0.05 | 1.11 | (0.05) |
| JRRM 125 | 79.2 | 18.7 | 0.13 | 0.01 | 0.50 | 0.69 | 0.08 | 0.00 | 0.07 | 0.04 | 0.30 | 0.02 | (0.07) |
| JRRM 123 | 79.1 | 13.3 | 0.13 | 0.01 | 4.13 | 0.10 | 1.32 | 0.01 | 0.29 | 0.80 | 0.45 | 0.00 | (0.03) |
| JRRM 122 | 78.2 | 10.2 | 0.43 | 0.81 | 0.24 | 2.05 | 0.65 | 0.20 | 1.04 | 4.89 | 1.03 | 0.20 | (0.12) |
| JRRM 124 | 73.9 | 16.5 | 1.09 | 0.11 | 2.60 | 1.79 | 0.10 | 0.24 | 0.31 | 0.19 | 2.74 | 0.11 | (0.10) |
| JRRM 127 | 68.5 | 23.0 | 0.18 | 0.27 | 0.92 | 0.54 | 0.15 | 0.17 | 1.75 | 1.78 | 2.19 | 0.04 | (0.07) |
| JRRM 126 | 66.9 | 21.3 | 0.45 | 0.65 | 3.34 | 3.13 | 0.12 | 0.03 | 0.28 | 0.49 | 2.84 | 0.04 | (0.17) |
| JRRM 129 | 62.2 | 30.1 | 0.15 | 0.10 | 1.46 | 1.92 | 2.23 | 0.01 | 0.23 | 0.20 | 0.96 | 0.11 | (0.11) |
| JRRM 128 | 54.3 | 26.0 | 2.80 | 0.85 | 4.45 | 1.84 | 3.10 | 0.24 | 0.37 | 3.36 | 1.37 | 1.01 | (0.02) |
| JRRM 130 | 53.4 | 32.7 | 1.95 | 1.05 | 0.53 | 1.42 | 0.61 | 0.37 | 2.32 | 0.91 | 3.35 | 0.83 | (0.11) |
| JRRM 131 | 52.7 | 36.6 | 0.78 | 0.07 | 2.20 | 2.61 | 1.02 | 0.03 | 0.76 | 1.61 | 1.16 | 0.26 | (0.17) |
| JRRM 132 | 50.6 | 39.1 | 1.29 | 0.11 | 1.64 | 0.79 | 0.34 | 0.11 | 2.16 | 2.38 | 0.29 | 0.75 | (0.15) |
| JRRM 133 | 50.1 | 39.0 | 0.10 | 1.27 | 3.69 | 0.91 | 2.03 | 0.01 | 0.33 | 0.34 | 1.93 | 0.57 | (0.08) |
| JRRM 134 | 47.2 | 44.3 | 0.20 | 0.24 | 1.07 | 0.37 | 0.20 | 0.24 | 0.13 | 3.83 | 1.74 | 0.35 | (0.14) |
| JRRM 135 | 37.2 | 48.9 | 2.36 | 0.42 | 3.05 | 2.77 | 1.24 | 0.04 | 2.87 | 0.48 | 0.07 | 0.20 | (0.18) |

CRM MAGNESIA REFRACTORY SET

| Number | SOLD IN SET/10 ONLY | | | | | certified values | | | | | | | informational values | | | | 20 g units | | | |
|----------|---------------------|--------------------------------|------|--------------------------------|------------------|-------------------------------|--------------------------------|------------------|------|-------------------|-------------------------------|------------------|-------------------------------|--------------------------------|------------------|------|-------------------|-------------------------------|------------------|--|
| | MgO | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | SiO ₂ | B ₂ O ₃ | Cr ₂ O ₃ | K ₂ O | MnO | Na ₂ O | P ₂ O ₅ | TiO ₂ | B ₂ O ₃ | Cr ₂ O ₃ | K ₂ O | MnO | Na ₂ O | P ₂ O ₅ | TiO ₂ | |
| JRRM 410 | 99.08 | 0.05 | 0.59 | 0.05 | 0.18 | 0.02 | 0.00 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | 0.02 | 0.00 | 0.00 | 0.01 | 0.00 | 0.02 | 0.00 | |
| JRRM 409 | 98.03 | 0.20 | 0.74 | 0.49 | 0.53 | 0.03 | 0.01 | 0.00 | 0.01 | 0.00 | 0.02 | 0.00 | 0.03 | 0.01 | 0.00 | 0.01 | 0.00 | 0.02 | 0.00 | |
| JRRM 408 | 96.19 | 2.55 | 0.67 | 0.13 | 0.46 | 0.09 | 0.00 | 0.00 | 0.01 | 0.00 | 0.01 | 0.00 | 0.09 | 0.00 | 0.00 | 0.01 | 0.00 | 0.01 | 0.00 | |
| JRRM 407 | 94.55 | 0.10 | 0.67 | 2.14 | 2.43 | 0.02 | 0.08 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | 0.02 | 0.08 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | |
| JRRM 405 | 91.95 | 1.37 | 1.69 | 1.34 | 3.47 | 0.01 | 0.01 | 0.01 | 0.07 | 0.00 | 0.12 | 0.05 | 0.01 | 0.01 | 0.01 | 0.07 | 0.00 | 0.12 | 0.05 | |
| JRRM 406 | 91.85 | 1.13 | 4.80 | 0.87 | 1.19 | 0.01 | 0.00 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | 0.01 | 0.00 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | |
| JRRM 404 | 88.02 | 6.01 | 1.78 | 2.90 | 1.22 | 0.01 | 0.00 | 0.00 | 0.03 | 0.00 | 0.05 | 0.00 | 0.01 | 0.00 | 0.00 | 0.03 | 0.00 | 0.05 | 0.00 | |
| JRRM 403 | 85.48 | 4.06 | 0.61 | 1.55 | 8.14 | 0.03 | 0.01 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | 0.03 | 0.01 | 0.00 | 0.01 | 0.00 | 0.04 | 0.00 | |
| JRRM 402 | 83.77 | 1.99 | 3.57 | 5.05 | 5.46 | 0.12 | 0.00 | 0.00 | 0.01 | 0.01 | 0.07 | 0.02 | 0.12 | 0.00 | 0.00 | 0.01 | 0.01 | 0.07 | 0.02 | |
| JRRM 401 | 81.24 | 8.10 | 0.20 | 3.89 | 6.42 | 0.01 | 0.00 | 0.00 | 0.01 | 0.00 | 0.03 | 0.01 | 0.01 | 0.00 | 0.00 | 0.01 | 0.00 | 0.03 | 0.01 | |

CRM SILICA REFRACTORY SETS

| Number | SOLD IN SETS ONLY, AS GROUPED | | | | | | | | | | | | 20 g units | | | | | | | | | | | | | |
|----------|--------------------------------|-------|--------------------------------|--------------------------------|------------------|-------|-------|-------------------|-------------------------------|------------------|------------------|------------------|--------------------------------|-------|--------------------------------|--------------------------------|------------------|-------|-------|-------------------|-------------------------------|------------------|------------------|------------------|-------|----------|
| | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SiO ₂ | TiO ₂ | ZrO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SiO ₂ | TiO ₂ | ZrO ₂ | | |
| JRRM 221 | 10.03 | 2.78 | 0.02 | 1.57 | 0.27 | 0.68 | 0.15 | 0.46 | 0.01 | 83.8 | 0.04 | 0.01 | 10.03 | 2.78 | 0.02 | 1.57 | 0.27 | 0.68 | 0.15 | 0.46 | 0.01 | 83.8 | 0.04 | 0.01 | 0.01 | new 2017 |
| JRRM 222 | 7.66 | 0.16 | 0.006 | 3.86 | 0.78 | 0.94 | 0.05 | 0.20 | 0.006 | 84.8 | 0.78 | 0.48 | 7.66 | 0.16 | 0.006 | 3.86 | 0.78 | 0.94 | 0.05 | 0.20 | 0.006 | 84.8 | 0.78 | 0.48 | 0.48 | |
| JRRM 223 | 5.22 | 4.14 | 0.03 | 2.04 | 0.37 | 0.27 | 0.20 | 0.69 | 0.01 | 86.0 | 0.04 | 0.67 | 5.22 | 4.14 | 0.03 | 2.04 | 0.37 | 0.27 | 0.20 | 0.69 | 0.01 | 86.0 | 0.04 | 0.67 | 0.67 | |
| JRRM 224 | 4.66 | 1.95 | 0.30 | 2.47 | 0.90 | 0.29 | 0.16 | 0.28 | 0.68 | 87.9 | 0.15 | 0.003 | 4.66 | 1.95 | 0.30 | 2.47 | 0.90 | 0.29 | 0.16 | 0.28 | 0.68 | 87.9 | 0.15 | 0.003 | 0.003 | |
| JRRM 225 | 3.22 | 3.19 | 0.01 | 1.27 | 0.63 | 0.13 | 0.07 | 0.90 | 0.01 | 89.9 | 0.42 | 0.01 | 3.22 | 3.19 | 0.01 | 1.27 | 0.63 | 0.13 | 0.07 | 0.90 | 0.01 | 89.9 | 0.42 | 0.01 | 0.01 | |
| JRRM 226 | 2.63 | 0.97 | 0.24 | 2.99 | 0.47 | 0.09 | 0.02 | 0.19 | 0.23 | 91.2 | 0.29 | 0.32 | 2.63 | 0.97 | 0.24 | 2.99 | 0.47 | 0.09 | 0.02 | 0.19 | 0.23 | 91.2 | 0.29 | 0.32 | 0.32 | |
| JRRM 227 | 1.66 | 2.41 | 0.45 | 0.81 | 0.11 | 0.05 | 0.23 | 0.05 | 0.003 | 92.9 | 0.09 | 0.88 | 1.66 | 2.41 | 0.45 | 0.81 | 0.11 | 0.05 | 0.23 | 0.05 | 0.003 | 92.9 | 0.09 | 0.88 | 0.88 | |
| JRRM 228 | 0.39 | 1.78 | 0.08 | 0.08 | 0.10 | 0.11 | 0.03 | 1.18 | 0.99 | 93.8 | 1.21 | 0.01 | 0.39 | 1.78 | 0.08 | 0.08 | 0.10 | 0.11 | 0.03 | 1.18 | 0.99 | 93.8 | 1.21 | 0.01 | 0.01 | |
| JRRM 229 | 1.17 | 1.41 | 0.37 | 0.19 | 0.07 | 0.46 | 0.07 | 0.07 | 0.01 | 95.7 | 0.12 | 0.20 | 1.17 | 1.41 | 0.37 | 0.19 | 0.07 | 0.46 | 0.07 | 0.07 | 0.01 | 95.7 | 0.12 | 0.20 | 0.20 | |
| JRRM 230 | 0.18 | 0.60 | 0.05 | 0.70 | 0.02 | 0.01 | 0.12 | 0.07 | 0.38 | 97.7 | 0.03 | 0.001 | 0.18 | 0.60 | 0.05 | 0.70 | 0.02 | 0.01 | 0.12 | 0.07 | 0.38 | 97.7 | 0.03 | 0.001 | 0.001 | |
| JRRM 231 | 0.63 | 0.005 | 0.18 | 0.04 | 0.004 | 0.004 | 0.004 | 0.006 | 0.001 | 98.6 | 0.003 | 0.38 | 0.63 | 0.005 | 0.18 | 0.04 | 0.004 | 0.004 | 0.004 | 0.006 | 0.001 | 98.6 | 0.003 | 0.38 | 0.38 | |
| JRRM 232 | 0.05 | 0.004 | 0.002 | 0.05 | 0.004 | 0.001 | 0.005 | 0.005 | 0.001 | 99.7 | 0.002 | (0.001) | 0.05 | 0.004 | 0.002 | 0.05 | 0.004 | 0.001 | 0.005 | 0.001 | 99.7 | 0.002 | (0.001) | (0.001) | | |

CRM SILICON CARBIDE REFRACTORY SET available in SET/9 ONLY 50 g

| Number | SiC | Tot C | Free C | LOI | Al | Ca | Fe | Mg | N | O | Ti | Free Si |
|-----------|-------|-------|--------|-------|-------|--------|-------|--------|-------|-------|--------|---------|
| JRRM 1001 | 99.58 | 29.81 | 0.04 | . | 0.008 | <0.001 | 0.044 | <0.001 | 0.030 | 0.048 | 0.0035 | 0.06 |
| JRRM 1002 | 0.06 | 5.03 | 4.98 | 5.11 | . | . | . | . | . | . | . | . |
| JRRM 1003 | . | 10.06 | 10.01 | 10.11 | . | . | . | . | . | . | . | . |
| JRRM 1004 | . | 20.04 | 19.92 | 20.01 | . | . | . | . | . | . | . | . |
| JRRM 1005 | . | 29.93 | 29.81 | 29.95 | . | . | . | . | . | . | . | . |
| JRRM 1006 | . | 49.99 | 49.97 | 49.95 | . | . | . | . | . | . | . | . |
| JRRM 1007 | 89.29 | 36.75 | 10.01 | . | . | . | . | . | . | . | . | . |
| JRRM 1008 | 29.74 | 14.12 | 5.21 | . | . | . | . | . | . | . | . | . |
| JRRM 1009 | 6.18 | 39.43 | 37.67 | . | . | . | . | . | . | . | . | . |

CRM ZIRCON-ZIRCONIA REFRACTORY SET

| Number | SOLD IN SET/10 ONLY | | | | | | | | | | | | 20 g units | | | | | | | | | | | | |
|----------|---------------------|------------------|------------------|--------------------------------|------|--------------------------------|--------------------------------|------------------|------|-------------------|-------------------------------|------------------|------------------|------------------|------------------|--------------------------------|------|--------------------------------|--------------------------------|------------------|------|-------------------|-------------------------------|------------------|------|
| | ZrO ₂ | HfO ₂ | SiO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | TiO ₂ | ZrO ₂ | HfO ₂ | SiO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | P ₂ O ₅ | TiO ₂ | |
| JRRM 601 | 92.0 | 1.59 | 0.26 | 0.11 | 5.58 | 0.00 | 0.10 | 0.00 | 0.06 | 0.00 | 0.00 | 0.16 | 92.0 | 1.59 | 0.26 | 0.11 | 5.58 | 0.00 | 0.10 | 0.00 | 0.06 | 0.00 | 0.00 | 0.16 | 0.16 |
| JRRM 602 | 88.4 | 1.52 | 0.33 | 0.07 | 0.22 | 0.01 | 1.62 | 0.00 | 5.30 | 0.76 | 1.34 | 0.16 | 88.4 | 1.52 | 0.33 | 0.07 | 0.22 | 0.01 | 1.62 | 0.00 | 5.30 | 0.76 | 1.34 | 0.16 | 0.16 |
| JRRM 603 | 84.8 | 1.45 | 0.96 | 5.29 | 0.95 | 0.02 | 2.86 | 0.65 | 0.96 | 0.18 | 0.83 | 0.93 | 84.8 | 1.45 | 0.96 | 5.29 | 0.95 | 0.02 | 2.86 | 0.65 | 0.96 | 0.18 | 0.83 | 0.93 | 0.93 |
| JRRM 604 | 79.4 | 1.35 | 3.05 | 6.93 | 0.09 | 3.06 | 0.43 | 1.94 | 0.01 | 1.09 | 1.99 | 0.13 | 79.4 | 1.35 | 3.05 | 6.93 | 0.09 | 3.06 | 0.43 | 1.94 | 0.01 | 1.09 | 1.99 | 0.13 | 0.13 |
| JRRM 605 | 75.5 | 1.31 | 10.8 | 4.84 | 1.94 | 1.55 | 0.17 | 0.54 | 1.99 | 0.45 | 0.35 | 0.12 | 75.5 | 1.31 | 10.8 | 4.84 | 1.94 | 1.55 | 0.17 | 0.54 | 1.99 | 0.45 | 0.35 | 0.12 | 0.12 |
| JRRM 606 | 72.5 | 1.26 | 22.1 | 0.53 | 0.02 | 0.00 | 0.93 | 0.01 | 0.32 | 2.03 | 0.01 | 0.11 | 72.5 | 1.26 | 22.1 | 0.53 | 0.02 | 0.00 | 0.93 | 0.01 | 0.32 | 2.03 | 0.01 | 0.11 | 0.11 |
| JRRM 607 | 61.6 | 1.21 | 32.9 | 3.53 | 0.04 | 0.00 | 0.12 | 0.04 | 0.03 | 0.02 | 0.08 | 0.13 | 61.6 | 1.21 | 32.9 | 3.53 | 0.04 | 0.00 | 0.12 | 0.04 | 0.03 | 0.02 | 0.08 | 0.13 | 0.13 |
| JRRM 608 | 58.8 | 1.21 | 34.6 | 0.70 | 0.52 | 0.49 | 0.09 | 0.01 | 3.12 | 0.03 | 0.11 | 0.10 | 58.8 | 1.21 | 34.6 | 0.70 | 0.52 | 0.49 | 0.09 | 0.01 | 3.12 | 0.03 | 0.11 | 0.10 | 0.10 |
| JRRM 609 | 55.6 | 1.12 | 40.5 | 0.88 | 0.30 | 0.01 | 0.15 | 0.02 | 0.15 | 0.94 | 0.08 | 0.15 | 55.6 | 1.12 | 40.5 | 0.88 | 0.30 | 0.01 | 0.15 | 0.02 | 0.15 | 0.94 | 0.08 | 0.15 | 0.15 |
| JRRM 610 | 48.7 | 0.98 | 45.7 | 0.45 | 3.07 | 0.00 | 0.30 | 0.01 | 0.54 | 0.04 | 0.11 | 0.09 | 48.7 | 0.98 | 45.7 | 0.45 | 3.07 | 0.00 | 0.30 | 0.01 | 0.54 | 0.04 | 0.11 | 0.09 | 0.09 |

RM RICE STRAW ASH - THERMOSTIL

| Number | typical analysis | | | | | | | | | | | | | | 100 g units | |
|---------|------------------|--------------------------------|------|-----------------|------|--------------------------------|------------------|-------|-------|-------------------|-------------------------------|-----------------|------------------|-------------------------|-------------|--|
| | SiO ₂ | Al ₂ O ₃ | C | CO ₂ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | SO ₃ | TiO ₂ | -H ₂ O 900°C | | |
| DH 5704 | 92.49 | 0.198 | 3.60 | 0.008 | 0.30 | 0.090 | 0.97 | 0.362 | 0.062 | 0.070 | 0.273 | 0.177 | 0.004 | 1.38 | | |
| DH 5708 | 86.67 | 1.15 | 3.83 | 0.094 | 0.97 | 0.931 | 0.872 | 3.10 | 0.117 | 0.085 | 0.226 | 0.255 | 0.126 | 1.70 | | |
| DH 5705 | 76.31 | 0.363 | 4.33 | 0.265 | 2.51 | 2.89 | 0.653 | 9.60 | 0.245 | 0.116 | 0.123 | 0.409 | 0.217 | 2.32 | | |

| RM | | SAND FOR SLIDING GATES | | | | | | | | | | | | | | | | typical analysis listed in mass % | | 100 g units | |
|---------|------------------|--------------------------------|-------|--------|--------------------------------|-------|------------------|------|--------------------------------|-------------------|-------|-------------------------------|-------|------------------|-------------------------------|-----------------|------------------|-----------------------------------|-------|-------------|--|
| Number | SiO ₂ | Al ₂ O ₃ | C | CaO | Cr ₂ O ₃ | Fe | K ₂ O | MgO | Mn ₃ O ₄ | Na ₂ O | NiO | P ₂ O ₅ | S | TiO ₂ | V ₂ O ₅ | WO ₃ | ZrO ₂ | -H ₂ O | 900°C | | |
| DH 4501 | 72.21 | 4.92 | 0.607 | 0.025 | 11.53 | 5.14 | 0.633 | 2.40 | 0.065 | 0.059 | 0.053 | 0.008 | . | 0.195 | 0.102 | . | . | . | 0.204 | | |
| DH 4502 | 65.97 | 5.69 | 0.47 | 0.038 | 14.75 | 6.31 | 0.693 | 3.24 | 0.074 | 0.062 | 0.033 | 0.007 | 0.010 | 0.203 | 0.110 | . | . | . | 0.177 | | |
| DH 4507 | 27.95 | 11.00 | 0.326 | 0.096 | 33.41 | 14.51 | . | 7.29 | 0.179 | . | 0.090 | CO ₂ | 0.013 | 0.486 | 0.270 | 0.019 | . | . | 0.129 | | |
| DH 4506 | 10.22 | 12.93 | 0.700 | <0.017 | 42.01 | 25.03 | . | 8.18 | 0.703 | . | . | CO ₂ | 0.007 | 0.510 | 0.382 | . | . | . | 0.091 | | |

| RM | | FOUNDRY SAND | | | | | | | | | | | | typical analysis listed in mass % except * which is mg/kg | | 100 g units | |
|---------|-------|------------------|--------------------------------|--------------------------------|-------|------|-----|-----|-------|------|------------------|-------------------------------|-------|---|--|-------------|--|
| Number | BaO | CeO ₂ | Co ₃ O ₄ | Cr ₂ O ₃ | CuO | La* | Li* | Nd* | NiO | Sr* | TiO ₂ | V ₂ O ₅ | ZnO | ZrO ₂ | | | |
| DH 3301 | 0.015 | 0.003 | 0.020 | 0.538 | 0.012 | 14.0 | 6.7 | 8.3 | 0.003 | 35.0 | 0.213 | 0.007 | 0.015 | 0.127 | | | |

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | SO ₃ | SiO ₂ |
|---------|--------------------------------|-------|--------------------------------|------------------|-------|--------------------------------|-------------------|-------------------------------|-----------------|------------------|
| DH 3301 | 2.76 | 0.720 | 3.84 | 0.169 | 0.570 | 0.070 | 0.297 | 0.027 | 0.116 | 90.36 |

| CRM | | ZIRCON SAND | | | | | | |
|-----------|-------------------------------------|--------------------------------|--------------------------------|------------------|------------------|------|-------|--|
| Number | ZrO ₂ + HfO ₂ | Al ₂ O ₃ | Fe ₂ O ₃ | SiO ₂ | TiO ₂ | LOI | Units | |
| JCRM R501 | 66.5 | 0.39 | 0.06 | 32.6 | 0.16 | 0.11 | 100 g | |
| JCRM R502 | 60.3 | 5.87 | 0.10 | 32.8 | 0.24 | 0.26 | 100 g | |

| CRM | | SILICA POWDER SET | | | | | | | SOLD IN SET/3 ONLY | | 100 g units | |
|-----------|--------------------------------|-------------------|--------------------------------|------------------|----------|-------------------|------------------|------------------|--------------------|--|-------------|--|
| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | Na ₂ O | SiO ₂ | TiO ₂ | LOI | | | |
| JCRM R405 | 1.07 | 0.029 | 0.053 | 0.71 | 0.023 | 0.060 | 97.78 | 0.022 | 0.13 | | | |
| JCRM R406 | 1.31 | 0.016 | 0.102 | 0.13 | 0.005 | 0.030 | 96.71 | 0.564 | 0.97 | | | |
| JCRM R404 | 0.0011 | 0.00002 | 0.00006 | 0.00004 | <0.00001 | 0.0001 | >99.99 | 0.0006 | 0.00 | | | |

| SILICA BRICK | | # = class, where 1 = CRM and 2 = RM analysis listed in mass % | | | | | | | | | | | | | NH, VS: 75g | | SRM: 45g | | others: 100g | |
|--------------|------------|---|--------------------------------|-------|-------|--------------------------------|--------------------------------|------------------|-------------------|-------|-------|-------------------|--------|-------------------------------|------------------|------|----------|--|--------------|--|
| # | Number | SiO ₂ | Al ₂ O ₃ | BaO | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | Li ₂ O | MgO | MnO | Na ₂ O | P | P ₂ O ₅ | TiO ₂ | LOI | | | | |
| 1 | VS K1/3 | 96.1 | 0.55 | . | 1.35 | . | 1.36 | . | . | 0.045 | 0.031 | . | 0.0122 | . | . | . | | | | |
| 1 | ECRM 777-1 | 95.06 | 0.795 | . | 2.826 | . | 0.330 | 0.154 | . | 0.071 | . | 0.02 | . | . | 0.444 | . | | | | |
| 1 | ECRM 776-1 | 62.76 | 29.28 | 0.122 | 0.31 | 0.022 | 1.43 | 2.92 | 0.019 | 0.476 | . | 0.488 | . | 0.062 | 1.62 | . | | | | |
| 1 | VS K2/4 | 58.6 | 35.1 | . | 0.4 | . | 2.94 | 0.69 | . | 0.48 | 0.06 | 0.19 | . | . | 1.91 | . | | | | |
| 1 | VS K3/2 | 32.3 | 63.6 | . | 0.44 | . | 1.15 | 0.15 | . | 0.27 | . | 0.17 | . | . | 1.34 | . | | | | |
| 1 | SRM 198 | . | 0.16 | . | 2.71 | . | 0.66 | 0.017 | 0.001 | 0.07 | . | 0.012 | . | 0.022 | 0.02 | 0.21 | | | | |
| 1 | SRM 199 | . | 0.48 | . | 2.41 | . | 0.74 | 0.094 | 0.002 | 0.13 | . | 0.015 | . | 0.015 | 0.06 | 0.17 | | | | |

| SILICEOUS MATERIAL | | # = class, where 1 = CRM, 2 = RM analysis listed in mass % | | | | | | | | | | | | | T = Total | |
|--------------------|--------------|--|--------------------------------|--------|--------------------------------|--------------------------------|------------------|--------|----------|-------------------|-------------------------------|------------------|--------|------------|--------------|--|
| # | Number | SiO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | TiO ₂ | LOI | Units | Other | |
| 1 | BCS 313/2 | 99.73 | 0.068 | 0.0160 | BaO:0.00067 | 0.0229 | 0.0108 | 0.0038 | 0.00032 | 0.0057 | . | 0.0243 | . | 100 g | SrO: 0.00024 | |
| 1 | NCS DC60116a | 98.32 | 1.10 | 0.038 | 0.00030 | 0.076 | 0.15 | 0.026 | 0.0013 | 0.076 | (0.0069) | 0.023 | 0.14 | 50 g | | |
| 1 | GBW 03113 | 95.74 | 2.36 | 0.17 | 0.00054 | 0.21 | 0.67 | 0.098 | (0.0033) | 0.25 | (0.0076) | 0.036 | 0.35 | 60 g | | |
| 1 | SRM 2696 | 95.61 | 0.2080 | 0.426 | . | (0.055) | 0.652 | 0.235 | 0.032 | (0.129) | (0.0863) | . | (2.11) | 70 g | ZnO:0.051 | |
| 1 | NCS DC60117a | 94.41 | 3.20 | 0.094 | 0.00034 | 0.088 | 1.26 | 0.025 | 0.0011 | 0.47 | (0.0070) | 0.019 | 0.27 | 50 g | | |
| 1 | GBW 03114 | 89.59 | 5.48 | 0.34 | 0.0012 | 0.48 | 2.07 | 0.16 | (0.010) | 1.09 | (0.014) | 0.102 | 0.53 | 60 g | | |
| 2 | CERAM CEB1 | 74.0 | 16.2 | 0.52 | BaO:0.05 | 0.48 | 1.75 | 0.16 | . | 0.71 | 0.14 | 0.34 | 5.60 | 25 or 100g | SrO: 0.02 | |
| 1 | GBW 03117 | 71.25 | 2.56 | 6.37 | . | 0.18 | 1.10 | 3.98 | . | 13.77 | . | 0.057 | 0.44 | 50 g | | |

| RM | | CERAMIC POWDER | | | | | | | | | | | | | | |
|---------|------------------|--------------------------------|------|--------------------------------|------------------|------|-------|------------------|--------|--------|--------|--------|--------|--------|--------|-------|
| Number | SiO ₂ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | MgO | MnO | TiO ₂ | Ba | Co | Cr | Cu | Ni | Sc | Zn | Units |
| SARM 69 | 66.6 | 14.4 | 2.37 | 7.18 | 1.96 | 1.85 | 0.129 | 0.777 | 0.0518 | 0.0028 | 0.0223 | 0.0046 | 0.0053 | 0.0020 | 0.0068 | 100 g |

CRM SYNTHETIC SILICATE WITH TRACE ELEMENTSMaterial base: SiO₂ 72%, Al₂O₃ 15%, Fe₂O₃ 4%, CaMg(CO₃)₂ pure dolomite 4%, Na₂SO₄ 2.5%, K₂SO₄ 2.5% analysis listed in mg/kg 70 g units

| Number | Ag | As | B | Ba | Be | Bi | Cd | Ce | Co | Cr | Cu | La | Li | Mn |
|-----------|---------|-----|------|-------|------|------|-------|------|------|------|------|-----|------|-------|
| GBW 07701 | (0.034) | 2.0 | 2.1 | 24 | 0.26 | 0.31 | 0.022 | 2.0 | 2.6 | 2.3 | 2.0 | 2.1 | 15 | 27 |
| GBW 07702 | 0.064 | 5.0 | 5.1 | 54 | 0.56 | 0.61 | 0.052 | 5.0 | 5.6 | 5.3 | 5.0 | 5.1 | 18 | 57 |
| GBW 07703 | 0.11 | 10 | 10.0 | 104 | 1.1 | 1.1 | 0.1 | 10.0 | 10.6 | 10.3 | 10.0 | 10 | 23 | 107 |
| GBW 07704 | 0.21 | 20 | 20 | 204 | 2.1 | 2.1 | 0.2 | 20 | 20.6 | 20.3 | 20.0 | 20 | 33 | 207 |
| GBW 07705 | 0.51 | 50 | 50 | 504 | 5.1 | 5.1 | 0.5 | 50 | 50.6 | 50 | 50 | 50 | 63 | 507 |
| GBW 07706 | 1.0 | 100 | 100 | 1000 | 10 | 10 | 1.0 | 100 | 101 | 100 | 100 | 100 | 113 | 1000 |
| GBW 07708 | 5.0 | 500 | 500 | 5000 | 50 | 50 | 5.0 | 500 | 500 | 500 | 500 | 500 | 513 | 5000 |
| GBW 07709 | 10.0 | . | 1000 | 10000 | 100 | 100 | 10 | 1000 | . | 1000 | 1000 | . | 1010 | 10000 |
| GBW 07710 | 20 | . | . | . | 200 | 200 | 20 | . | . | . | 2000 | . | . | . |
| GBW 07711 | 50 | . | . | . | 500 | . | 50 | . | . | . | 5000 | . | . | . |

continued

| Number | Mo | Nb | Ni | Pb | Sb | Sn | Sr | Ti | V | W | Y | Yb | Zn | Zr |
|-----------|------|------|------|------|------|------|------|-------|------|------|-----|-----|------|------|
| GBW 07701 | 0.21 | 2.3 | 2.6 | 2.5 | 0.28 | 0.28 | 5.0 | 24 | 2.8 | 0.20 | 2.0 | 0.2 | 3.0 | 2.2 |
| GBW 07702 | 0.51 | 5.3 | 5.6 | 5.5 | 0.58 | 0.58 | 8.0 | 54 | 5.8 | 0.50 | 5.0 | 0.5 | 6.0 | 5.2 |
| GBW 07703 | 1.0 | 10.3 | 10.6 | 10.5 | 1.1 | 1.1 | 13 | 104 | 10.8 | 1.0 | 10 | 1.0 | 11.0 | 10.2 |
| GBW 07704 | 2.0 | 20.3 | 20.6 | 20.5 | 2.1 | 2.1 | 23 | 204 | 20.8 | 2.0 | 20 | 2.0 | 21 | 20 |
| GBW 07705 | 5.0 | 50 | 50.6 | 50 | 5.1 | 5.1 | 53 | 504 | 51 | 5.0 | 50 | 5.0 | 51 | 50 |
| GBW 07706 | 10 | 100 | 101 | 100 | 10 | 10 | 103 | 1000 | 101 | 10 | 100 | 10 | 101 | 100 |
| GBW 07708 | 50 | 500 | 500 | 500 | 50 | 50 | 500 | 5000 | 500 | 50 | 500 | 50 | 500 | 500 |
| GBW 07709 | 100 | . | . | 1000 | 100 | 100 | 1000 | 10000 | 1000 | 100 | . | 100 | 1000 | 1000 |
| GBW 07710 | 200 | . | . | 2000 | 200 | 200 | 2000 | 20000 | . | 200 | . | . | 2000 | . |
| GBW 07711 | 500 | . | . | 5000 | 500 | 500 | 5000 | . | . | 500 | . | . | 5000 | . |

CRM SILICON METAL POWDER

analysis listed in mass %

| Number | Al | C | Ca | Cr | Cu | Fe | Mg | Mn | Ni | P | S | Ti | V | Zr |
|-------------|--------|----------|-----------|-----------|-----------|--------|--------|---------|---------|---------|----------|--------|----------|---------|
| NCS DC25007 | 0.24 | . | 0.31 | . | . | 0.39 | . | . | . | . | . | . | . | . |
| SRM 57B | 0.1690 | (0.0200) | (0.00222) | (0.00173) | (0.00172) | 0.3400 | . | 0.00782 | 0.00153 | 0.00163 | (0.0030) | 0.0346 | (0.0025) | 0.00178 |
| IPT 134 | 0.085 | 0.025 | 0.102 | 0.0011 | 0.0014 | 0.29 | 0.0048 | 0.0113 | 0.0006 | 0.0033 | 0.002 | 0.0097 | 0.0004 | . |
| IPT 135 | 0.045 | 0.018 | 0.011 | 0.0006 | 0.0008 | 0.125 | 0.0012 | 0.0070 | 0.0005 | 0.0027 | 0.002 | 0.0113 | 0.0003 | . |
| NCS HC25649 | 0.032 | . | 0.060 | . | . | 0.53 | . | . | . | 0.0067 | . | 0.026 | . | . |
| NCS HC25648 | 0.026 | . | 0.055 | . | . | 0.44 | . | . | . | 0.0065 | . | 0.023 | . | . |

CRM SILICON CARBIDE

analysis listed in mass %

| Number | SiC | Al | Fe | Units |
|---------|------|---------|--------|-------|
| VS K9/2 | 99.6 | (0.002) | (0.06) | 150 g |

CRM SILICON CARBIDE

in the chart below, (F) = Free and (T) = Total analysis listed in mass % except * which is mg/kg

| Number | C (T) | C (F) | Si (T) | Si (F) | SiO ₂ (F) | Al | B | Ca | Cr | Cu | Fe | K | Mg |
|------------|--------|---------|--------|---------|----------------------|----------|----------|----------|----------|----------|----------|----------|----------|
| ECRM 781-1 | 48.251 | (37.22) | 35.56 | (4.66) | . | 4.39 (T) | (0.0149) | (0.0433) | (0.0240) | . | (0.8061) | (0.3765) | (0.0421) |
| NMIJ 8002a | 29.93 | . | 68.01 | . | . | 0.0189 | . | 0.00619 | 0.0115 | 0.0130 | . | . | . |
| BAM S008 | 29.9 | 0.045 | . | (<0.03) | (<0.01) | 0.0047 | 0.00030 | 0.00025 | 0.000016 | 0.000010 | 0.00048 | . | 0.000007 |
| NMIJ 8001a | 29.80 | . | 68.31 | . | . | 0.00832 | . | . | . | . | 0.00467 | . | . |
| ECRM 780-1 | 26.381 | . | 63.5 | . | . | 1.86 (T) | . | 0.84 | . | . | 1.30 (T) | (0.0112) | 0.051 |
| BCS 360 | 23.53 | (0.085) | 60.8 | (0.54) | . | 6.52 | . | 0.115 | . | . | (0.19) | . | . |
| BCS 359 | 23.46 | (0.061) | 67.6 | (0.32) | . | 0.118 | . | 0.108 | . | . | 0.175 | . | . |

| Number | Mn | Mo | N | Na | Ni | O | Ti | V | Y* | Zr | Notes | Units |
|------------|----------|--------|----------|----------|----------|--------|----------|----------|------|---------|--|-------|
| ECRM 781-1 | (0.0274) | . | (0.0282) | (0.0308) | (0.0210) | . | (0.0320) | (0.0216) | . | . | P: (0.0117) Mo: (0.0264) | 100 g |
| NMIJ 8002a | 0.000160 | 0.0109 | . | . | . | . | 0.00477 | . | 0.58 | . | Beta Phase | 50 g |
| BAM S008 | 0.000005 | 0.0018 | . | 0.000017 | 0.00009 | 0.0146 | 0.0067 | 0.0275 | . | 0.00044 | (SiC-6H:99.7, SiC-15R:0.23, SiC-4H:0.06) | 50 g |
| NMIJ 8001a | . | . | . | . | . | . | 0.000637 | . | 0.31 | . | Alpha Phase | 50 g |
| ECRM 780-1 | 0.029 | . | 0.325 | (0.0502) | . | . | . | . | . | . | n/a | 100 g |
| BCS 360 | . | . | (4.77) | . | . | (4.03) | 0.025 | . | . | . | Sialon Bonded | 100 g |
| BCS 359 | . | . | (7.84) | . | . | (0.53) | 0.022 | . | . | . | Nitrogen Bearing | 100 g |

CRM SILICON CARBIDE SET SOLD IN SET/3 ONLY F = Free T = Total 50 g each

| Number | T.Si | F.Si | F.SiO ₂ | T.C | F.C | Al | Ca | Cl | Cr | Cu | F | Fe | Mg |
|-----------|-------|---------|--------------------|-------|---------|--------|--------|----------|-----------|-----------|----------|--------|-----------|
| JCRM R024 | 68.97 | (0.042) | (0.593) | 29.85 | (0.423) | 0.0193 | 0.0019 | (<0.002) | 0.0056 | (<0.0006) | (<0.001) | 0.0219 | 0.0002 |
| JCRM R025 | 68.43 | (0.014) | (0.356) | 30.49 | (1.24) | 0.0184 | 0.0008 | (<0.002) | 0.0097 | 0.0021 | (0.0574) | 0.0233 | (<0.0001) |
| JCRM R026 | 69.03 | (0.012) | (0.311) | 29.85 | (0.598) | 0.0059 | 0.0004 | (<0.002) | (<0.0005) | (<0.0006) | 0.0686 | 0.0011 | (<0.0001) |

| Number | Mn | Mo | N | Ni | O | P | S | Ti | V | Zn | Zr |
|-----------|-----------|----------|---------|----------|------|---------|----------|--------|--------|-----------|-----------|
| JCRM R024 | 0.0004 | (<0.001) | (0.048) | 0.0060 | 0.97 | (<0.01) | (<0.005) | 0.0340 | 0.0013 | (<0.0005) | 0.0047 |
| JCRM R025 | (<0.0003) | 0.0126 | 0.113 | 0.0011 | 0.94 | (<0.01) | (0.0431) | 0.0040 | 0.0053 | (<0.0005) | 0.0012 |
| JCRM R026 | (<0.0003) | (<0.001) | 0.034 | (<0.001) | 0.71 | (<0.01) | (<0.005) | 0.0016 | 0.0018 | (<0.0005) | (<0.0005) |

CRM SILICON NITRIDE

analysis in mass %

analysis in mg/kg

| Number | Si | N | Al | C | Ca | Fe | O | Co | Mg | Cr | Mn | Na | Ni | Ti | W | Zr | ̑-phase of Si ₃ N ₄ | Units |
|------------|--------|--------|---------|-------|---------|---------|--------|-------|-----|-------|----|-------|-------|------|-------|------|---|-------|
| INDIVIDUAL | | | | | | | | | | | | | | | | | | |
| SRM 8983 | . | 39.23 | . | 0.107 | . | . | 1.20 | . | . | . | . | . | . | . | . | . | . | 4.5 g |
| NMIJ 8004a | 59.226 | 38.485 | 0.07397 | . | 0.00727 | 0.01969 | . | 10.29 | . | 2.987 | . | 2.485 | 8.519 | . | 2.146 | . | last | 25 g |
| BAM ED101 | . | 38.1 | 0.0469 | 0.162 | 0.00141 | 0.00795 | (1.91) | 43.5 | 4.3 | . | . | 7.59 | . | 41.3 | . | 7.43 | last | 50 g |
| SET ONLY | | | | | | | | | | | | | | | | | | |
| JCRM R006 | 59.57 | 38.98 | <0.002 | 0.101 | <0.0003 | 0.0012 | 1.18 | . | <2 | <6 | <1 | . | <8 | <4 | . | <7 | . | 20 g |
| JCRM R007 | 59.45 | 39.13 | 0.0707 | 0.136 | 0.0931 | 0.0169 | 0.79 | . | 68 | 49 | 28 | . | <8 | 58 | . | <7 | . | 20 g |
| JCRM R008 | 59.03 | 38.46 | 0.116 | 0.097 | 0.225 | 0.171 | 1.56 | . | 12 | 92 | 86 | . | <8 | 72 | . | 9 | . | 20 g |

CRM SILICON NITRIDE

analysis listed in mass %

SRM 656 is two 10 g powder units, one ̑ and one ̒ phase powder

| Number | Powder | Mass ̑ | Uncertainty ± | Mass ̒ | Uncertainty ± | Amorphous | Uncertainty ± |
|---------|--------|--------|---------------|--------|---------------|-----------|---------------|
| SRM 656 | ̑ | 87.5 | 0.59 | 3.0 | 0.05 | 9.5 | 0.61 |
| SRM 656 | ̒ | 16.3 | 0.81 | 75.1 | 2.54 | 8.6 | 0.60 |

CRM BORON NITRIDE

analysis listed in mass %

T = Total

AO = adherent oxide

50 g units

| Number | B.T | B.AO | N | Al | C | Ca | Co | Cr | Fe | H ₂ O | Mg | Na | O | Si | Ti |
|-----------|------|-------|------|---------|---------|--------|------------|---------|---------|------------------|--------|---------|------|--------|---------|
| BAM ED103 | 43.5 | 0.070 | 55.6 | 0.00070 | (0.018) | 0.0273 | (<0.00001) | 0.00047 | 0.00150 | (<0.1) | 0.0056 | 0.00123 | 0.68 | 0.0017 | 0.00049 |

CRM SILICOALUMINUM

analysis listed in mass %

| Number | Al | Si | Fe | Ba | C | Ca | Co | Cr | Cu | Mg | Mn | Ni | P | S | Sr | Ti | Units |
|-------------|-------|-------|-------|-------|------|-------|--------|--------|--------|-------|-------|--------|-------|--------|-------|-------|-------|
| NCS HC14605 | 36.67 | 25.94 | 24.97 | 9.12 | 0.13 | 1.33 | . | 0.152 | 0.045 | . | 0.12 | 0.167 | 0.018 | 0.012 | . | . | 70 g |
| NCS HC93615 | 34.80 | 29.87 | 30.47 | . | . | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| NCS HC14603 | 32.84 | 24.12 | 33.54 | 7.57 | 0.13 | 0.71 | . | 0.085 | 0.061 | . | 0.14 | 0.042 | 0.015 | 0.015 | . | . | 70 g |
| NCS HC14602 | 32.82 | 19.21 | 38.09 | 6.52 | 0.14 | 0.85 | . | 0.017 | 0.137 | . | 0.25 | 0.014 | 0.015 | 0.013 | . | . | 70 g |
| NCS HC13602 | 32.55 | 32.01 | 20.59 | 7.41 | 0.27 | 1.17 | . | . | . | 0.85 | 0.197 | . | 0.017 | 0.0096 | . | . | 50 g |
| NCS HC93614 | 31.91 | 33.75 | 27.84 | . | . | . | . | . | . | . | . | . | . | . | . | . | 50 g |
| NCS HC93633 | 29.67 | 28.31 | 37.44 | 0.45 | . | . | . | . | . | . | 0.426 | . | 0.023 | 0.022 | . | . | 50 g |
| NCS HC14604 | 25.44 | 19.21 | 49.14 | 2.64 | 0.24 | 0.44 | . | 0.053 | 0.172 | . | 0.25 | 0.018 | 0.011 | 0.011 | . | . | 70 g |
| NCS HC28635 | 16.63 | 43.60 | 17.53 | 1.64 | 1.00 | 15.18 | . | 0.054 | 0.046 | 0.027 | 0.095 | 0.026 | 0.051 | 0.040 | 0.023 | . | 50 g |
| NCS HC14609 | 14.46 | 33.41 | 35.46 | 7.72 | 0.22 | 5.74 | . | 0.116 | 0.32 | 0.18 | 0.33 | 0.016 | 0.018 | 0.017 | 0.092 | 0.055 | 60 g |
| NCS HC14610 | 13.47 | 40.58 | 23.25 | 10.70 | 0.24 | 8.25 | 0.0032 | 0.032 | 0.29 | 0.12 | 0.23 | 0.012 | 0.021 | 0.025 | 0.094 | 0.124 | 60 g |
| NCS HC14608 | 9.14 | 53.39 | 14.22 | 12.39 | 0.13 | 8.28 | 0.0022 | 0.021 | 0.176 | 0.21 | 0.17 | 0.0061 | 0.022 | 0.021 | 0.132 | 0.084 | 60 g |
| NCS HC28637 | 5.42 | 52.76 | 10.54 | 15.33 | 0.63 | 11.10 | . | 0.031 | 0.017 | 0.019 | 0.15 | 0.0070 | 0.024 | 0.073 | 0.042 | . | 50 g |
| NCS HC28636 | 4.07 | 50.36 | 16.68 | 24.26 | 0.34 | 1.44 | . | 0.083 | 0.032 | 0.032 | 0.11 | 0.021 | 0.016 | 0.038 | 0.095 | . | 50 g |
| NCS HC14611 | 1.47 | 56.74 | 5.77 | 17.00 | 1.56 | 13.61 | 0.0016 | 0.0044 | 0.0097 | 0.045 | 0.065 | 0.0020 | 0.016 | 0.14 | 0.22 | 0.126 | 60 g |

CRM SILICOBARIUM

analysis listed in mass %

50 g units

| Number | Ba | Si | Al | C | Ca | Fe | Mg | Mn | P | S | Sr |
|-------------|-------|-------|-------|------|-------|-------|-------|-------|-------|-------|-------|
| NCS HC93632 | 27.54 | 47.56 | 2.78 | 0.99 | . | 11.75 | . | 0.16 | 0.024 | 0.13 | . |
| NCS HC93634 | 14.14 | 52.62 | 1.82 | 0.64 | 14.08 | 12.97 | 0.051 | 0.104 | 0.022 | 0.204 | 0.063 |
| NCS HC93631 | 10.00 | 37.19 | 13.46 | 0.78 | 5.16 | 27.56 | 0.098 | 0.43 | 0.032 | 0.044 | . |

SILICOCALCIUM

= class, where 1 = CRM and 2 = RM

| # | Number | Ca | Si | Al | Ba | C | Cr | Cu | Fe | Mg | Mn | Mo | Ni | P | S | Ti | Units |
|---|---------------|-------|-------|-------------|-------|-------|--------|-------|-------|-------|-------|----------|---------|-------|--------|--------------|-------|
| 1 | NCS HC93613 | 31.67 | 56.20 | 1.77 | . | 1.30 | . | . | 5.58 | . | . | . | . | 0.018 | 0.088 | . | 50 g |
| 1 | BS 119 | 31.3 | 62.9 | 0.44 | . | 0.30 | . | . | 3.03 | . | . | . | . | 0.034 | 0.012 | 17025 | 100 g |
| 1 | NCS HC37620 | 30.70 | 60.09 | So1.Al:1.09 | . | 0.68 | . | . | . | . | . | . | . | 0.017 | 0.033 | . | 50 g |
| 2 | DH 0406 | 30.48 | 60.79 | 0.333 | . | 0.61 | <0.009 | . | 5.17 | 0.010 | 0.056 | . | . | 0.031 | 0.022 | 0.019 | 50 g |
| 1 | NCS HC11604a | 30.45 | 56.02 | 1.97 | . | 0.94 | . | . | 6.93 | . | 0.037 | . | . | 0.054 | 0.073 | . | 50 g |
| 1 | VS F26/3 | 29.9 | 60.1 | 1.52 | . | . | . | . | 6.19 | . | . | . | . | 0.024 | 0.029 | 0.161 | 100 g |
| 1 | VS F26/2 | 29.9 | 59.5 | 1.52 | . | . | . | . | 6.29 | . | . | . | . | 0.024 | 0.030 | 0.156 | 100 g |
| 2 | DH 0402 | 28.48 | 58.68 | 1.13 | . | . | 0.010 | 0.014 | 6.74 | 0.049 | 0.051 | . | . | 0.014 | . | 0.055 | 50 g |
| 1 | 58A CQ42001 | 28.25 | 55.31 | 1.88 | . | 2.44 | . | . | 6.08 | . | . | . | . | 0.019 | 0.0132 | . | 50 g |
| 1 | NCS HC93627 | 28.02 | 57.43 | 1.76 | . | 1.02 | . | . | 6.94 | . | . | . | . | 0.030 | 0.045 | . | 50 g |
| 1 | NCS HC11619 | 27.15 | 61.11 | 2.15 | . | 0.55 | . | . | 6.61 | . | 0.053 | . | . | 0.048 | 0.029 | . | 50 g |
| 1 | NCS HC37621 | 25.25 | 60.19 | So1.Al:1.55 | . | 0.71 | . | . | . | . | . | . | . | 0.031 | 0.020 | . | 50 g |
| 1 | VS F25/3 | 21.3 | 51.5 | 0.67 | . | . | . | . | 23.06 | . | . | . | . | 0.011 | 0.0056 | . | 100 g |
| 1 | NCS HC11605 | 13.22 | 53.46 | 2.34 | 14.02 | 0.385 | 0.054 | 0.079 | 13.57 | 0.022 | 0.075 | Sr:0.235 | 0.023 | 0.014 | 0.039 | . | 60 g |
| 1 | VS F44 | 12.6 | 49.7 | 1.68 | . | 0.166 | . | . | 19.91 | 9.6 | . | . | V:(2.5) | 0.014 | 0.0066 | . | 100 g |

CRM SILICOCHROMIUM

| Number | Cr | Si | Fe | Al | B | C | Co | Cu | Mn | Ni | P | S | Ti | V | Units |
|-------------|-------|-------|------|-------|--------|-------|-------|-------|-------|------|--------|---------|------|------|-------|
| SRM 689 | 36.4 | 39.5 | 23.2 | 0.049 | 0.0017 | 0.043 | 0.034 | 0.013 | 0.32 | 0.20 | 0.026 | 0.002 | 0.40 | 0.09 | 100 g |
| NCS HC25633 | 33.90 | 44.06 | . | 1.00 | . | 0.045 | . | . | 0.29 | . | 0.013 | (0.002) | . | . | 50 g |
| NCS HC25643 | 32.62 | 49.17 | . | 1.24 | . | 0.018 | . | . | 0.429 | . | 0.0083 | 0.0025 | . | . | 50 g |

SILICOMANGANESE

= class, where 1 = CRM and 2 = RM

| # | Number | Mn | Si | Fe | C | Co | Cr | Cu | Ni | P | S | Ti | V | Units |
|---|--------------|-------|-------|-------|-------|---------|-------|---------|---------|-------|----------|--------|--------|-------|
| 2 | DH 0107 | 77.82 | 17.36 | 2.64 | 1.65 | 0.030 | . | 0.012 | 0.020 | 0.135 | . | 0.122 | 0.015 | 50 g |
| 1 | BS SiMn-1 | 73.2 | 16.0 | 8.2 | 1.80 | (0.051) | 0.019 | (0.042) | (0.083) | 0.278 | 0.016 | (0.19) | (0.04) | 100 g |
| 1 | MHCX04 | 70.0 | 22.7 | 5.31 | 0.80 | 0.111 | 0.103 | 0.104 | 0.122 | 0.149 | 0.012 | 0.28 | 0.085 | 70 g |
| 1 | NCS HC25605b | 69.77 | 14.20 | . | 2.21 | . | . | . | . | 0.153 | 0.0052 | . | . | 50 g |
| 1 | NCS HC25657 | 67.96 | 25.03 | . | 0.58 | . | . | . | . | 0.065 | 0.011 | 0.18 | . | 50 g |
| 1 | VS F23/1 | 67.53 | 21.18 | . | 1.45 | . | . | . | . | 0.235 | 0.0155 | . | . | 100 g |
| 1 | NCS HC26611b | 67.44 | 18.24 | . | 1.24 | . | . | . | . | 0.080 | 0.009 | . | . | 50 g |
| 1 | NCS HC28618 | 67.40 | 19.34 | 11.65 | 1.05 | 0.017 | 0.045 | 0.051 | 0.036 | 0.107 | 0.017 | 0.255 | 0.063 | 50 g |
| 1 | NCS HC25605c | 67.20 | 21.87 | 10.01 | 0.456 | 0.020 | 0.029 | 0.019 | 0.013 | 0.132 | 0.0076 | 0.175 | 0.040 | 50 g |
| 1 | NCS HC37612 | 67.02 | 18.96 | . | 1.10 | . | . | . | . | 0.178 | 0.016 | 0.276 | . | 50 g |
| 1 | NCS HC18603 | 66.70 | 17.21 | . | 1.70 | . | . | . | . | 0.183 | 0.025 | . | . | 50 g |
| 1 | NCS HC93619 | 66.40 | 17.55 | . | 1.65 | . | . | . | . | 0.137 | 0.025 | . | . | 50 g |
| 1 | NCS HC11603b | 66.37 | 17.63 | . | 1.34 | . | . | . | . | 0.065 | 0.008 | . | . | 100 g |
| 1 | NCS HC25605a | 66.30 | 18.28 | . | 1.09 | . | . | . | . | 0.145 | 0.010 | 0.18 | . | 50 g |
| 1 | NCS HC19607 | 66.20 | 18.41 | . | 1.56 | . | . | . | . | 0.126 | 0.022 | . | . | 50 g |
| 1 | NCS HC25640 | 65.85 | 24.74 | . | 0.181 | . | . | . | . | 0.104 | 0.010 | . | . | 50 g |
| 1 | NCS HC93625 | 65.74 | 17.19 | . | 1.66 | . | . | . | . | 0.151 | 0.026 | . | . | 50 g |
| 1 | NCS HC93637 | 65.70 | 17.54 | . | 1.80 | . | . | . | . | 0.023 | 0.023 | . | . | 100 g |
| 1 | NCS HC11603a | 65.67 | 17.49 | . | 1.33 | . | . | . | . | 0.065 | 0.011 | . | . | 100 g |
| 1 | NCS HC37605 | 65.64 | 17.67 | . | 1.55 | . | . | . | . | 0.14 | 0.024 | 0.46 | . | 50 g |
| 1 | YSBC25616-97 | 65.54 | 21.88 | . | 0.435 | . | . | . | . | 0.141 | 0.009 | . | . | 50 g |
| 1 | NCS HC25640a | 65.50 | 24.47 | . | 0.197 | . | . | . | . | 0.117 | 0.0079 | . | . | 50 g |
| 1 | MHCX03 | 65.5 | 29.2 | 4.92 | 0.04 | . | 0.20 | 0.13 | 0.11 | 0.047 | 0.004 | 0.49 | . | 100 g |
| 1 | NCS HC25654 | 65.29 | 19.26 | . | 0.876 | . | . | . | . | 0.109 | 0.0122 | 0.19 | . | 50 g |
| 2 | DH 0106 | 65.24 | 18.38 | 14.60 | 1.21 | 0.013 | 0.011 | 0.017 | 0.042 | 0.080 | 0.010 | 0.121 | 0.015 | 50 g |
| 1 | NCS HC28617 | 64.97 | 17.59 | 15.16 | 1.57 | 0.035 | 0.055 | 0.096 | 0.092 | 0.127 | 0.018 | 0.221 | 0.060 | 50 g |
| 1 | NCS HC93624 | 64.86 | 16.87 | . | 1.79 | . | . | . | . | 0.120 | 0.024 | . | . | 50 g |
| 1 | NCS HC93618 | 63.91 | 19.04 | . | 1.13 | . | . | . | . | 0.140 | 0.022 | . | . | 50 g |
| 1 | NCS HC93626 | 63.80 | 16.42 | . | 1.91 | . | . | . | . | 0.097 | 0.020 | . | . | 50 g |
| 1 | JSS 705-5 | 62.69 | 14.99 | . | 1.941 | . | . | . | . | 0.239 | (0.0087) | . | . | 150 g |
| 1 | NCS HC28616 | 62.53 | 14.33 | 20.00 | 2.28 | 0.048 | 0.060 | 0.080 | 0.167 | 0.205 | 0.020 | 0.222 | 0.095 | 50 g |
| 2 | DH 0302 | 62.17 | 27.95 | 8.69 | 0.092 | 0.054 | 0.016 | 0.031 | 0.032 | 0.089 | . | 0.288 | 0.015 | 50 g |
| 1 | NCS HC26621 | 61.49 | 27.49 | . | 0.039 | . | . | . | . | 0.072 | 0.009 | 0.24 | . | 50 g |
| 2 | DH 0303 | 60.60 | 30.66 | 7.52 | 0.029 | 0.023 | 0.024 | 0.016 | 0.039 | 0.059 | . | 0.444 | 0.015 | 50 g |
| 1 | NCS HC25641 | 60.29 | 27.88 | . | 0.082 | . | . | . | . | 0.078 | 0.0069 | 0.41 | . | 50 g |
| 1 | NCS HC37606b | 60.13 | 13.87 | . | 2.26 | . | . | . | . | 0.42 | 0.040 | 0.25 | . | 50 g |
| 1 | NCS HC25646 | 59.34 | 32.90 | . | 0.018 | . | . | . | . | 0.043 | 0.0034 | 0.24 | . | 50 g |
| 2 | DH 0301 | 59.06 | 30.16 | 9.91 | 0.015 | 0.028 | 0.035 | 0.019 | 0.033 | 0.050 | . | 0.471 | 0.015 | 50 g |
| 1 | 58A CQ43001 | 58.49 | 15.96 | . | 1.84 | . | . | . | . | 0.32 | 0.026 | . | . | 50 g |
| 1 | NCS HC26620 | 54.97 | 19.15 | . | 0.40 | . | . | . | . | 0.060 | 0.011 | 0.24 | . | 50 g |

the below continuation shows only the samples with more data

| Number | As | B | Ca | Pb | Sb | Zr |
|------------------|---------|----------|---------|----------|-------------|----------|
| BS SiMn-1 | (0.010) | (0.014) | (0.05) | (0.0005) | (0.002) | (0.0007) |
| MHCX04 | 0.004 | 0.03 | (0.010) | 0.008 | . | Zn:0.012 |
| NCS HC28618 | 0.0099 | . | . | 0.0007 | 0.0004 | . |
| NCS HC25605c | . | 0.010 | . | . | . | . |
| NCS HC25605a | . | 0.0063 | . | . | . | . |
| MHCX03 | . | (0.009) | . | . | . | . |
| NCS HC25654 | . | 0.022 | . | . | . | . |
| NCS HC28617 | 0.010 | . | . | 0.0005 | 0.0009 | . |
| JSS 705-5 | . | 0.0231 | . | . | . | . |
| NCS HC28616 | 0.015 | . | . | 0.0012 | 0.0025 | . |
| DH 0302 | . | . | . | . | . | 0.008 |
| NCS HC25641 | . | 0.021 | . | . | . | . |
| NCS HC25646 | . | 0.048 | . | . | . | . |
| DH 0301 | . | 0.0048 | . | . | . | . |
| NCS HC26620 | . | Al:0.015 | . | . | Mo:(0.0009) | . |

17025 many more values on certificate
Al: 0.010, Mo: 0.016, Sn: 0.008

SILICOZIRCONIUM

| Number | Zr | Si | Fe | Al | C | Ca | Cr | Cu | Hf | Mn | N | Ni | P | S | Ti | Units |
|--------------|-------|-------|------|-------|-------|-------|-------|------|-------|-------|-------|-------|-------|---------|-------|-------|
| CRM VS F27/2 | 51.5 | 26.1 | (12) | 7.48 | 0.111 | . | . | 1.47 | . | . | . | . | 0.044 | (0.001) | 0.215 | 100 g |
| RM DH 3001 | 36.06 | 51.14 | 8.87 | 0.852 | 0.338 | 0.157 | 0.004 | . | 0.804 | 0.210 | 0.027 | 0.013 | 0.033 | 0.002 | 0.073 | 50 g |

CRM BASIC SLAG

analysis listed in mass % 100 g units

| Number | Al | B | Ca | Cr | F | Fe | K | Mg | Mn | Na | P | S | Si | Ti | V | Zn |
|-------------|-------|--------|-------|--------|-------|-------|-------|------|-------|-------|-------|-------|-------|-------|-------|--------|
| IRSID 802-1 | 8.53 | 0.0245 | 30.62 | 0.0053 | 0.243 | 0.576 | 0.491 | 2.87 | 0.460 | 0.236 | 0.109 | 0.714 | 15.16 | 0.366 | 0.028 | 0.0025 |
| ECRM 804-1 | 0.407 | . | 36.88 | . | . | 11.92 | . | 0.88 | 1.48 | . | 7.67 | 0.127 | 2.59 | 0.152 | 0.460 | . |

CRM FERROALLOY SLAG

100 g units

| Number | Al ₂ O ₃ | BaO | CaO | Cr ₂ O ₃ | Fe ₂ O ₃ | K ₂ O | MgO | MnO | Na ₂ O | S | SiO ₂ | SrO | TiO ₂ |
|-----------|--------------------------------|------|-------|--------------------------------|--------------------------------|------------------|------|-------|-------------------|-------|------------------|-------|------------------|
| AMIS 0536 | 8.17 | 1.72 | 27.6 | 0.064 | 1.69 | 0.608 | 5.27 | 26.0 | . | 0.85 | 26.7 | 0.57 | 0.35 |
| AMIS 0535 | 7.51 | 1.55 | 27.1 | 0.071 | 0.76 | 0.761 | 5.89 | 26.3 | 0.50 | 0.76 | 27.0 | 0.48 | 0.34 |
| AMIS 0533 | 6.99 | 1.63 | 26.73 | 0.119 | 0.908 | 0.649 | 6.31 | 23.69 | 0.417 | 0.827 | 30.09 | 0.504 | 0.338 |

IRON MAKING SLAG

= class, where 1 = CRM and 2 = RM

| # | Number | CaO | SiO ₂ | Al ₂ O ₃ | C | Fe | FeO | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | S | TiO ₂ | Units |
|---|---------------|-------|------------------|--------------------------------|-------|-------|------|------------------|-------|--------|-------------------|-------------------------------|--------|------------------|-------|
| 1 | NH 7-1-009 | 49.6 | 32.8 | 9.2 | . | 0.47 | . | (0.19) | 1.1 | 0.60 | (0.14) | . | 1.17 | 0.38 | 75 g |
| 2 | BS Slag 2 | 44.6 | 36.9 | 10.3 | (0.2) | 0.24 | . | 0.16 | 5.9 | 0.19 | 0.16 | . | 1.16 | 0.204 | 50 g |
| 1 | IRSID 803-1 | 43.28 | 36.38 | 13.19 | . | 0.613 | . | . | 4.05 | 0.713 | . | 0.270 | 0.767 | 0.502 | 100 g |
| 1 | IRSID 802-1 * | 42.84 | 32.43 | 16.12 | . | 0.576 | . | . | 4.76 | 0.593 | . | 0.250 | 0.714 | 0.611 | 100 g |
| 1 | NH 7-1-008 | 42.1 | 39.1 | 8.4 | . | 0.30 | . | (0.52) | 6.1 | 0.73 | (0.33) | . | (0.65) | 0.30 | 75 g |
| 1 | NH 7-1-005 | 38.8 | 35.3 | 10.0 | . | 0.21 | . | (0.19) | 12.0 | 0.47 | (0.13) | . | (0.85) | 0.32 | 75 g |
| 1 | CAN SL-1 | 37.48 | 35.73 | 9.63 | . | . | 0.92 | (0.51) | 12.27 | (0.86) | (0.39) | . | 1.26 | (0.38) | 200 g |
| 2 | BS 100A | 37.0 | 35.3 | 10.10 | (0.2) | 0.29 | . | (0.5) | 12.85 | 0.33 | (0.2) | 0.0034 | 1.77 | 0.48 | 100 g |
| 1 | NH 7-1-010 | 31.2 | 44.0 | 7.94 | . | 5.5 | . | (0.59) | 0.73 | 3.40 | (0.18) | . | 0.14 | 0.91 | 75 g |
| 1 | NH 7-1-007 | 31.2 | 39.0 | 6.2 | . | 0.55 | . | (0.38) | 18.9 | 0.78 | (0.24) | . | (0.57) | 0.39 | 75 g |
| 1 | NH 7-1-014 | 30.1 | 33.6 | 24.0 | . | 1.27 | . | (0.07) | 9.3 | (0.3) | (0.07) | . | (0.02) | (0.07) | 75 g |
| 1 | NH 7-1-011 | 29.4 | 21.9 | 24.0 | . | 1.9 | . | (0.04) | 17.5 | 1.97 | (0.19) | . | (0.03) | (0.09) | 75 g |
| 1 | NH 7-1-013 | 28.7 | 20.3 | 38.6 | . | 1.12 | . | (0.03) | 8.0 | 0.26 | (0.04) | . | (0.03) | 0.78 | 75 g |
| 1 | NH 7-1-015 | 28.0 | (44.7) | 14.5 | . | 1.7 | . | (0.08) | 9.2 | 0.58 | (0.1) | . | (0.02) | (0.08) | 75 g |

* Oxides Calculated, see previous chart "BASIC SLAG" for actual certified values

STEEL MAKING SLAG

= class, where 1 = CRM and 2 = RM

CMSI, GBW, RH: 50 g units NH: 75 g units all others: 100 g units

| # | Number | CaO | T.Ca | CaF ₂ | SiO ₂ | Al ₂ O ₃ | Cr ₂ O ₃ | F | Fe | FeO | K ₂ O | MgO | MnO | Na ₂ O | P ₂ O ₅ | s.P ₂ O ₅ | S | TiO ₂ | V ₂ O ₅ |
|---|-------------|--------|-------|------------------|------------------|--------------------------------|--------------------------------|-------|--------|--------|------------------|---------|--------|-------------------|-------------------------------|---------------------------------|-------|------------------|-------------------------------|
| 2 | RH02 | 64.7 | . | . | (12.9) | 11.0 | (0.03) | . | (0.2) | (0.2) | . | 3.5 | 0.024 | . | P:0.003 | . | (0.9) | 0.07 | . |
| 1 | JK S11 * | . | 60.0 | . | 26.8 | 2.85 | 0.17 | (7.9) | . | (0.2)* | . | 4.7 | 0.12 | . | (<0.005) | . | 0.30 | 0.95 | (<0.01) |
| 2 | BS 101/3 | 54.4 | . | . | 18.8 | 1.42 | . | . | 10.9 | . | 0.005 | 3.0 | 5.0 | 0.027 | 0.74 | . | 0.18 | (0.9) | . |
| 1 | CMSI 1745 | . | 37.64 | 1.41 | 14.91 | 1.78 | . | . | 13.38 | 12.33 | . | 9.28 | 1.86 | . | 1.02 | . | 0.097 | 0.42 | . |
| 2 | BS 101/1 | 52.9 | . | . | 23.3 | 0.70 | . | . | 5.8 | . | 0.008 | 8.7 | 3.47 | 0.013 | 0.76 | . | 0.19 | 0.8 | . |
| 1 | BS 101/4 | 52.5 | . | . | 16.7 | 0.86 | . | . | (13.3) | . | 0.007 | 4.8 | 4.79 | 0.018 | 0.81 | . | 0.15 | 1.16 | . |
| 1 | IRSID 804-1 | 51.60 | . | . | 5.54 | (0.79) | . | . | 11.92 | . | . | 1.46 | 1.91 | . | 17.58 | . | 0.127 | 0.25 | 0.82 |
| 1 | BCS 381 | 49.0 | . | . | 8.78 | 0.67 | 0.33 | . | 13.3 | 3.69 | . | 1.03 | 3.16 | . | 15.7 | 15.2 | 0.19 | 0.35 | 0.94 |
| 1 | IRSID 805-1 | 48.92 | . | . | 6.63 | 0.616 | . | . | 14.87 | . | . | 1.86 | 2.05 | . | 16.20 | . | 0.092 | 0.342 | 0.918 |
| 2 | BS 101/2 | 47.6 | . | . | 16.9 | 0.91 | . | . | 15.1 | . | 0.008 | 7.0 | 4.8 | 0.031 | 0.63 | . | 0.20 | (0.8) | . |
| 1 | IRSID 806-1 | 46.13 | . | . | 11.72 | 0.901 | . | . | 17.89 | . | . | 3.02 | 5.94 | . | 2.25 | . | 0.110 | 0.504 | 0.514 |
| 2 | BS 101/5 | 46.1 | . | . | 15.2 | 0.74 | . | . | 19.4 | . | 0.0044 | 5.0 | 5.7 | (0.04) | 0.71 | . | 0.12 | 1.2 | . |
| 1 | ECRM 879-1 | 43.70 | . | . | 8.82 | 0.803 | 0.477 | 0.368 | 18.97 | . | . | 2.19 | 4.45 | . | 8.46 | 7.59 | 0.102 | 0.535 | 0.738 |
| 1 | NH 143 | 42.90 | . | . | 4.88 | (0.50) | 0.97 | . | 14.53 | 8.62 | . | 5.29 | 2.84 | . | 16.71 | . | 0.083 | 0.15 | . |
| 1 | NH 146 | 40.56 | . | . | 11.38 | 4.29 | 0.69 | . | 20.30 | 18.47 | . | 5.47 | 5.52 | . | 2.11 | . | 0.165 | 0.39 | . |
| 1 | NH 151 | 34.83 | . | . | 15.97 | 2.06 | 0.65 | . | 14.94 | 0.14 | . | 5.05 | 8.44 | . | 7.92 | . | 0.079 | 0.53 | . |
| 1 | NH 156 | 34.66 | . | . | 15.20 | 7.80 | 0.75 | . | 16.35 | 0.14 | . | 4.66 | 3.81 | . | 5.98 | . | 0.111 | 0.36 | . |
| 1 | CMSI 1744 | 26.73 | . | . | 8.91 | 3.92 | . | . | 34.33 | 36.55 | . | 12.15 | 2.01 | . | 0.87 | . | 0.107 | 0.32 | . |
| 1 | VS W4/4 | 25.5 | . | . | 16.7 | 3.62 | . | . | 23.2 | 25.5 | . | 18.3 | 4.17 | . | P: 0.259 | . | 0.037 | 1.02 | . |
| 1 | NH 150 | 21.77 | . | . | 15.69 | 3.23 | 1.74 | . | 24.23 | 27.30 | . | (14.46) | 8.16 | . | 0.62 | . | 0.044 | 0.15 | . |
| 1 | NH 152 | 21.95 | . | . | 15.91 | 2.60 | 28.67 | . | 14.40 | 12.79 | . | 6.17 | 4.85 | . | (0.12) | . | 0.028 | 0.37 | . |
| 1 | NH 145 | 20.85 | . | . | 22.43 | 2.39 | 0.99 | . | 27.97 | 30.46 | . | 2.71 | 9.26 | . | 2.05 | . | 0.089 | 0.56 | . |
| 1 | NH 149 | 9.85 | . | . | 8.42 | 3.36 | 53.81 | . | 14.09 | 8.12 | . | 2.89 | 3.74 | . | (0.03) | . | 0.040 | 0.22 | . |
| 1 | SARM 77 | 3.64 | . | . | 26.8 | 27.5 | 12.5 | . | 5.31T | . | . | 22.99 | . | . | . | . | 0.32T | . | . |
| 1 | NH 154 | (1.15) | . | . | 48.67 | 3.68 | 1.54 | . | 10.65 | 13.36 | . | 2.44 | (28.0) | . | (0.03) | . | 0.074 | 0.27 | . |

* JK S11 lists total Fe as FeO

CONVERTER SLAG

= class, where 1 = CRM and 2 = RM

| # | Number | CaO | Ca | SiO ₂ | Al ₂ O ₃ | Fe | K ₂ O | MgO | Mn | MnO | Nb ₂ O ₅ | P ₂ O ₅ | S | TiO ₂ | V ₂ O ₅ |
|---|-------------|-------|-------|------------------|--------------------------------|-------|------------------|-------|------|------|--------------------------------|-------------------------------|-------|------------------|-------------------------------|
| 2 | DH 3911 | 50.50 | . | 8.58 | 0.933 | 18.51 | . | 1.54 | 4.42 | . | 0.055 | 2.65 | 0.160 | 0.350 | 0.590 |
| 2 | DH 3908 | 47.13 | . | 12.70 | 1.096 | 18.96 | 0.008 | 2.513 | 4.31 | . | 0.072 | 1.488 | 0.110 | 0.558 | 0.273 |
| 1 | NCS HC28810 | . | 33.35 | 14.45 | 1.76 | 16.52 | . | 7.10 | . | 2.78 | . | 1.60 | 0.120 | 1.25 | . |
| 1 | NCS HC28809 | . | 32.65 | 15.40 | 4.38 | 13.50 | . | 7.75 | . | 2.30 | . | 1.67 | 0.195 | 1.02 | . |

| Number | Cr | CuO | SrO | ZnO | Units |
|-------------|-------|-------|-------|-------|-------|
| DH 3911 | 0.154 | 0.007 | . | 0.003 | 100 g |
| DH 3908 | 0.331 | . | 0.028 | . | 100 g |
| NCS HC28810 | . | . | . | . | 80 g |
| NCS HC28809 | . | . | . | . | 80 g |

CRM ELECTRIC FURNACE SLAG

50 g units

| Number | Ca(tot) | Al ₂ O ₃ | F | FeO | T.Fe | MgO | MnO | P ₂ O ₅ | S | SiO ₂ | TiO ₂ |
|-----------|---------|--------------------------------|------|---------|-------|-------|-------|-------------------------------|-------|------------------|------------------|
| CMSI 1757 | 28.87 | 8.72 | 0.82 | 1.89 | 2.26 | 15.66 | 2.39 | 0.030 | 0.25 | 24.80 | 0.25 |
| CMSI 1756 | 16.22 | 4.00 | 0.17 | (15.25) | 13.11 | 21.18 | 13.16 | 0.125 | 0.036 | 21.37 | 0.18 |

CRM FLUORINE SLAG

100 g units

| Number | F | T.CaF ₂ | Ca | CaO | Al ₂ O ₃ | C | FeO | MgO | MnO | P | SiO ₂ | TiO ₂ | V ₂ O ₅ |
|-----------|--------|--------------------|-------|------|--------------------------------|-------|------|--------|------|-------|------------------|------------------|-------------------------------|
| JK S10 | 34.4 | 70.7 | 50.8 | 20.3 | 0.54 | 0.022 | 0.10 | 0.30 | 0.03 | 0.002 | 7.8 | 0.05 | (<0.01) |
| IMZ EZP 1 | 31.62 | . | 36.76 | . | 24.85 | . | . | (0.85) | . | . | 2.61 | . | . |
| IMZ EZP 3 | 15.78 | . | 39.53 | . | 19.13 | . | . | 8.44 | . | . | 1.68 | . | . |
| IMZ EZP 2 | (0.89) | . | 24.03 | . | 41.38 | . | . | 16.89 | . | . | 5.81 | . | . |

MANGANESE SLAG

analysis listed in mass %

DH: RM, 100 g units

VS: CRM, 150 g units

| Number | Mn | Mn ₃ O ₄ | Al ₂ O ₃ | C | CaO | CuO | Fe | Fe ₂ O ₃ | K ₂ O | MgO | P | P ₂ O ₅ | S | SiO ₂ | ZnO |
|-----------|------|--------------------------------|--------------------------------|-------|-------|------|-------|--------------------------------|------------------|-------|---|-------------------------------|-------|------------------|-------|
| VS SH11/1 | 48.0 | . | . | . | . | . | . | . | . | . | . | 0.014 | . | . | . |
| DH 7403 | 4.93 | . | 19.84 | . | 15.95 | . | 0.088 | . | 1.30 | 12.34 | . | 0.002 | 0.818 | 43.23 | . |
| DH 7404 | 2.66 | . | 24.61 | . | 26.16 | . | 0.086 | . | 0.630 | 7.04 | . | 0.003 | 0.959 | 37.39 | . |
| DH 7402 | . | 0.113 | 5.99 | 11.92 | 0.405 | 7.02 | . | 3.96 | 0.164 | 0.118 | . | 14.03 | 0.114 | 11.01 | 45.16 |

| Number | Ba | CO ₂ | Cr ₂ O ₃ | Na ₂ O | SnO ₂ | SrO | TiO ₂ | Y ₂ O ₃ | ZrO ₂ | -H ₂ O@900°C |
|-----------|---------|-----------------|--------------------------------|-------------------|------------------|-------|------------------|-------------------------------|------------------|-------------------------|
| VS SH11/1 | . | . | . | . | . | . | . | . | . | . |
| DH 7403 | (0.475) | 0.032 | 0.007 | 0.433 | . | 0.083 | 0.100 | (0.009) | 0.039 | 0.062 |
| DH 7404 | 0.925 | . | 0.007 | (0.229) | . | 0.109 | 0.164 | 0.014 | 0.035 | . |
| DH 7402 | . | . | 0.086 | 0.133 | 0.386 | . | 0.274 | . | 0.024 | 0.077 |

CRM PHOSPHATE SLAG

| Number | total P ₂ O ₅ | citric acid sol. P ₂ O ₅ | CaO | SiO ₂ | Units |
|-----------|--|---|-------|------------------|-------|
| BAM 826-1 | 14.65 | 10.73 | 46.48 | 8.96 | 100 g |
| BAM 827-1 | 20.70 | 18.79 | 47.38 | 6.21 | 100 g |

CRM SLAG

analysis listed in mass %

| Number | Al ₂ O ₃ | C | Ca | CaO | F | Fe | FeO | K ₂ O | MgO | Mn | MnO | Na ₂ O | P | P ₂ O ₅ | S | SiO ₂ | TiO ₂ | V ₂ O ₅ | Units |
|-------------|--------------------------------|-------|-------|-------|------|-------|-------|------------------|-------|-------|-------|-------------------|--------|-------------------------------|-------|------------------|------------------|-------------------------------|-------|
| NCS HC18809 | 21.94 | . | 35.21 | . | . | 0.30 | . | . | 6.55 | . | 0.18 | . | . | 0.024 | 0.69 | 16.50 | 1.03 | . | 100 g |
| NCS HC28808 | 18.05 | . | . | 35.71 | . | 0.48 | 0.55 | 0.42 | 10.92 | . | 0.542 | 0.36 | . | 0.027 | 0.885 | 29.62 | 0.753 | . | 50 g |
| NCS HC28806 | 16.92 | . | . | 37.53 | . | 0.211 | 0.35 | 0.46 | 10.80 | . | 0.414 | 0.39 | . | 0.013 | 1.15 | 30.36 | 0.762 | . | 50 g |
| NCS HC18807 | 16.48 | . | . | 35.77 | . | 1.10 | . | . | 8.77 | . | 0.74 | . | . | 0.009 | 0.90 | 33.04 | 0.73 | . | 100 g |
| NCS HC18806 | 14.11 | . | . | 38.84 | . | 0.60 | . | . | 8.45 | . | 0.30 | . | . | 0.008 | 1.13 | 32.75 | 2.63 | . | 100 g |
| NCS HC25801 | 4.91 | . | . | 7.79 | . | 1.77 | . | . | 3.99 | 35.31 | . | . | 0.0056 | . | 0.66 | 33.47 | . | . | 50 g |
| NCS HC28807 | 3.67 | . | 32.32 | . | 0.76 | 13.54 | 10.44 | 0.033 | 7.27 | . | 4.06 | 0.057 | . | 1.72 | 0.134 | 14.54 | 1.13 | . | 50 g |
| NCS HC18808 | 1.25 | . | 24.10 | . | . | 25.55 | . | . | 11.66 | . | 3.34 | . | . | 2.00 | 0.13 | 13.44 | 2.22 | . | 100 g |
| NCS HC15804 | . | 0.014 | . | . | . | 0.22 | . | . | . | 44.42 | . | . | 0.0032 | . | 0.32 | 25.16 | . | . | 100 g |

CRM TIN SLAG

| Number | Sn | Al ₂ O ₃ | CaO | FeO | SiO ₂ | Units |
|-------------|-------|--------------------------------|-------|-------|------------------|-------|
| NCS HC35801 | 11.96 | 7.36 | 4.12 | 46.18 | 19.61 | 70 g |
| NCS HC35802 | 2.32 | 9.32 | 19.76 | 22.22 | 37.49 | 70 g |

CRM TITANIUM SLAG

100 g units

| Number | TiO ₂ | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | T.Fe | MgO | MnO | S | SiO ₂ | V ₂ O ₅ |
|----------------|------------------|--------------------------------|-------|--------------------------------|------|------|------|-------|------------------|-------------------------------|
| NCS HC19815 | 94.69 | 2.62 | 0.287 | . | 1.02 | 2.67 | 1.21 | 0.166 | 1.92 | . |
| DSZU 123.23-95 | 85.21 | 3.40 | 0.76 | 1.12 | 3.29 | 0.60 | 0.94 | 0.16 | 2.50 | 0.30 |
| DSZU 123.24-01 | 85.19 | 3.28 | . | 0.76 | 3.69 | . | 0.85 | 0.12 | 2.88 | 0.31 |
| NCS HC19814 | 84.94 | 3.04 | 1.83 | . | 1.08 | 7.27 | 0.74 | 0.247 | 4.13 | . |
| NCS HC19813 | 77.66 | 2.64 | 1.52 | . | 6.43 | 5.28 | 1.08 | 0.118 | 5.50 | . |

RM TUNDISH SLAG

typical analysis listed in mass %

100 g units

| Number | CaO | SiO ₂ | MgO | Al ₂ O ₃ | CO ₂ | Fe ₂ O ₃ | K ₂ O | MnO | Mn ₃ O ₄ | Na ₂ O | P ₂ O ₅ | S | SO ₃ | TiO ₂ | -H ₂ O 900°C |
|---------|-------|------------------|-------|--------------------------------|-----------------|--------------------------------|------------------|-------|--------------------------------|-------------------|-------------------------------|-------|-----------------|------------------|-------------------------|
| DH 5903 | 61.81 | 21.33 | 10.07 | 3.70 | 0.016 | 1.36 | 0.599 | . | 0.046 | 0.174 | 0.151 | 0.194 | . | 0.183 | 0.023 |
| DH 5904 | 39.47 | 30.68 | 14.53 | 12.37 | . | 0.591 | 0.622 | . | 0.076 | 0.530 | 0.019 | . | . | 0.038 | . |
| DH 6604 | 1.609 | 24.75 | 64.45 | 1.884 | 0.35 | 4.62 | 0.089 | . | 0.098 | 0.516 | 0.084 | . | 0.026 | 0.141 | 1.02 |
| DH 6606 | 1.37 | 27.46 | 62.63 | 1.30 | 0.16 | 4.93 | 0.070 | 0.093 | . | . | 0.055 | 0.022 | . | 0.103 | 1.15 |
| DH 6605 | . | . | . | . | 0.40179 | . | . | . | . | 0.347 | . | . | . | . | 1.15451 |

| Number | BaO | C Tot | Cr ₂ O ₃ | F | NiO | SrO |
|---------|-------|-------|--------------------------------|-------|-------|-------|
| DH 5903 | . | . | . | 0.046 | . | . |
| DH 5904 | 0.180 | . | . | . | . | 0.025 |
| DH 6604 | . | 0.471 | 0.255 | . | 0.165 | . |
| DH 6606 | . | . | . | . | . | . |
| DH 6605 | . | . | . | . | . | . |

VACUUM SLAG

NCS: CRM, 80 g units

DH: RM, 100 g units

| Number | CaO | Al ₂ O ₃ | Cr | Cr ₂ O ₃ | Fe | K ₂ O | MgO | Mn | MnO | Nb ₂ O ₅ | P ₂ O ₅ | S | SiO ₂ | SrO | TiO ₂ | V ₂ O ₅ | ZrO ₂ | Other |
|-------------|-------|--------------------------------|-------|--------------------------------|-------|------------------|-------|-------|-------|--------------------------------|-------------------------------|-------|------------------|-------|------------------|-------------------------------|------------------|--------------|
| DH 5119 | 53.43 | 17.04 | . | 0.075 | 2.54 | 0.012 | 10.84 | 2.48 | . | 0.270 | 0.102 | 0.188 | 8.30 | 0.033 | 1.37 | 0.053 | 0.240 | |
| DH 5120 | 52.90 | 20.33 | . | 0.039 | 1.55 | 0.011 | 11.68 | 1.27 | . | 0.202 | 0.039 | 0.281 | 8.13 | 0.032 | 1.28 | 0.016 | 0.230 | |
| DH 5121 | 51.14 | 23.56 | 0.039 | . | 1.27 | 0.011 | 11.98 | 0.769 | . | 0.109 | 0.028 | 0.369 | 7.63 | 0.031 | 0.869 | 0.012 | 0.232 | |
| NCS HC19810 | 2.04 | 1.25 | . | 0.93 | 31.26 | . | 1.9 | . | 10.67 | . | 0.046 | 0.052 | 18.25 | . | 10.02 | 17.2 | . | Met.Fe: 0.22 |

CRM VANADIUM SLAG

| Number | Al ₂ O ₃ | CaO | Cr ₂ O ₃ | Fe | MgO | MnO | P | SiO ₂ | TiO ₂ | V ₂ O ₅ | Units |
|----------|--------------------------------|------|--------------------------------|------|------|------|-------|------------------|------------------|-------------------------------|-------|
| VS SH9/3 | 1.76 | 1.61 | 3.32 | 28.9 | 3.53 | 9.73 | 0.015 | 16.63 | 7.39 | 22.2 | 150 g |

CRM SLUDGE

analysis listed in mg/kg except % which is mass %

SRM 2781: 40 g

SRM 2782: 70 g

all others: 40 g units

| Number | Type | Ag | Al% | As | Ba | Be | Bi | Ca% | Cd | Ce | Cl | Co | Cr | Cu | Fe% | Ga | Hg | In |
|-------------|-----------------|------|------|------|-----|------|------|------|-------|------|----|------|-------|-------|------|----|------|-----|
| SRM 2782 | industrial | 30.6 | 1.37 | 166 | 254 | . | . | 0.67 | 4.17 | 1240 | . | 66.3 | 109 | 2594 | 26.9 | 35 | 1.10 | 238 |
| BCR 146R | industrial | . | . | . | . | . | . | . | 18.8 | . | . | 7.39 | 196 | 838 | . | . | 8.62 | . |
| BCR 145R | mixed | . | . | . | . | . | . | . | 3.50 | . | . | 5.61 | (313) | 696 | . | . | 2.01 | . |
| SRM 2781 | domestic | 98 | 1.6 | 7.82 | . | . | . | 3.9 | 12.78 | . | . | . | 202 | 627.4 | 2.8 | . | 3.64 | . |
| IRNT WT-L * | water treatment | 11.9 | 3.03 | 8.87 | 781 | 3.73 | 3.73 | 8.80 | 1.97 | . | . | 6.77 | 79.0 | 136 | 1.70 | . | 4.25 | . |
| IRNT WT-M * | water treatment | 40.4 | 2.61 | 9.84 | 787 | 72 | . | 5.15 | 11.9 | . | . | 8.20 | 939 | 959 | 1.74 | . | 14.3 | . |
| BCR 143R | amended soil | . | . | . | . | . | . | . | 71.8 | . | . | 12.3 | (577) | 130.6 | . | . | 1.10 | . |

continued SRM 2782 also contains (2.1%) Carbon and trace informational values for Au, Eu, Hf, Rb, Sc, Sm, Ta, Tb, Th, U, Y, and Yb.

| Number | K% | La | Li | Mg% | Mn | Mo | N% | Na% | Ni | P% | Pb | S% | Sb | Se | Se | Sn | Sr | Ti | V | Zn |
|-------------|---------|------|-------|-------|-------|-------|------|---------|-------|-------|-------|-------|-------|-------|--------|------|------|--------|------|------|
| SRM 2782 | 0.32 | 58.1 | (5.0) | 0.26 | (300) | 10.07 | . | 1.30 | 154.1 | 0.50 | 574 | (0.2) | (2.0) | 0.44 | (20.3) | . | . | 0.0880 | 80 | 1254 |
| BCR 146R | . | . | . | . | 324 | . | . | . | 69.7 | . | 609 | . | . | . | . | . | . | . | . | 3061 |
| BCR 145R | . | . | . | . | 156 | . | . | . | 247 | . | 286 | . | . | . | . | . | . | . | . | 2122 |
| SRM 2781 | 0.49 | . | . | 0.59 | 46.7 | 4.78 | 0.21 | 80.2 | 2.42 | 202.1 | . | . | 16.0 | 5.1 | . | . | 0.32 | . | 1273 | |
| IRNT WT-L * | (0.695) | . | . | 0.781 | 390 | . | . | (0.414) | 32.0 | 0.881 | 122 | 1.02 | 17.8 | . | . | . | 170 | . | 41.3 | 1310 |
| IRNT WT-M * | 0.589 | . | . | 0.613 | 942 | . | . | (0.303) | 240 | 1.58 | 841 | 1.03 | 12.7 | . | . | 20.3 | 160 | . | 34.2 | 3080 |
| BCR 143R | . | . | . | . | (904) | . | . | . | 299 | . | 179.7 | . | . | (0.6) | . | . | . | . | . | 1055 |

* IRNT certificates expired, however use and sales continue without problems worldwide

CRM RED SLURRY SET

analysis listed in mass % 50 g units, SET ONLY

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | SiO ₂ | TiO ₂ |
|-----------|--------------------------------|------|--------------------------------|------------------|------------------|
| ShK411-01 | 12.7 | 5.67 | 58.7 | 4.57 | 4.85 |
| ShK412-01 | 13.3 | 7.0 | 57.3 | 4.67 | 4.32 |
| ShK413-01 | 11.2 | 11.8 | 35.6 | 22.3 | 2.09 |
| ShK414-01 | 11.4 | 3.35 | 67.2 | 3.44 | 3.27 |
| ShK415-01 | 17.3 | 4.04 | 44.4 | 10.5 | 7.11 |
| ShK416-01 | 15.1 | 5.13 | 52.1 | 7.35 | 5.92 |

RM SODA ASH

analysis listed in mass % 100 g units

| Number | Na ₂ CO ₃ | NaCl | Fe ₂ O ₃ | Na ₂ SO ₄ |
|---------|---------------------------------|-------|--------------------------------|---------------------------------|
| BCS 526 | 99.74 | 0.126 | 0.0005 | 0.008 |

CRM SURFACE AREAdata listed in m²/g

| Number | Multipoint +/- | Single Point +/- | Units |
|----------|----------------|------------------|----------------------------|
| SRM 2206 | 10.99 0.68 | 10.73 0.68 | 5 g granulated glass |
| SRM 1900 | 2.85 0.09 | 2.79 0.07 | 4 g silicon nitride powder |

CRM TENSILE CREEP

| Number | Creep Rate at 400 h | Time to 2% Strain | Time to 4% Strain | Units |
|---------|---|-------------------|-------------------|-------------------------|
| BCR 425 | 72 x 10 ⁻⁶ h ⁻¹ ± 5 | 278 h ± 16 | 557 h ± 30 | 3 rods 14 mm Ø x 150 mm |

CRM TENSILE STRENGTH and HARDNESS

data shows estimates of (material, measurement) uncertainty

| Number | ksi Tensile Strength | ksi Yield Strength | % Total Elongation | % Reduction | Hardness | Material | Units |
|-----------|----------------------|--------------------|--------------------|------------------|-----------------------|---------------|----------------------|
| BS TRM-2A | 125.6 (0.7, 2.4) | 116.2 (0.9, 2.7) | 20.1 (0.6, 1.8) | 57.7 (0.9, 2.7) | HRB 103 (2, 6) | 600 nickel | rod 25 mm Ø x 158 mm |
| BS TRM-3 | 98.2 (0.6, 5.5) | 44.7 (0.3, 3.1) | 52.0 (1.2, 10.8) | 57.1 (1.9, 17.3) | HRB 86.3 (0.7, 6.3) | 304 steel | sheet 30 cm x 30 cm |
| BS TRM-1 | 93.3 (0.3, 2.1) | 89.3 (0.5, 3.2) | 15.6 (0.2, 1.6) | 55.0 (0.4, 2.7) | . | 1018 steel | rod 25 mm Ø x 158 mm |
| BS TRM-1A | 83.9 (0.3, 1.7) | 70.2 (0.2, 1.5) | 18.8 (0.3, 1.8) | 56.9 (0.5, 3.2) | . | 1018 steel | rod 25 mm Ø x 158 mm |
| BS TRM-4 | 36.0 (0.1, 0.8) | 28.4 (0.1, 0.7) | 11.4 (0.1, 1.1) | (37.0) - - | HR15T 71.9 (0.6, 5.4) | 5056 aluminum | sheet 30 cm x 30 cm |

CRM TENSILE STRENGTH

| Number | 0.2% Proof Stress (MPa) | 0.5% Proof Stress (MPa) | Tensile Strength (MPa) | Elongation Fracture (A in %) | Reduction in Area at Fracture (Z in %) | Units |
|----------|-------------------------|-------------------------|------------------------|------------------------------|--|------------------------|
| BCR 661B | 300 ± 7 | 318 ± 7 | 750 ± 13 | 40.9 ± 0.9 | 60 ± 4 | 1 rod 14 mm Ø x 500 mm |

CRM BORON CARBIDE

analysis listed in mass %

analysis listed in mg/kg

100 g

| Number | Tot.B | Sol.B | B Isotopic Abundance | C | N | O | Al | Ca | Co | Cr | Cu | Fe | Mn | Na | Ni | Si | Ti | Zr |
|-----------|-------|-------|----------------------|-------|-------|------|-----|----|------|-----|-----|-----|------|----|-----|-----|----|------|
| BAM ED102 | 78.47 | 0.116 | 19.907 | 21.01 | 0.209 | 0.10 | 157 | 97 | 0.39 | 5.6 | 2.2 | 686 | 10.4 | 63 | 8.0 | 268 | 96 | 48.9 |

CRM CHROMIUM CARBIDE

analysis listed in mass %

| Number | C | Cr | S | Si | Units |
|-------------|-------|-------|-------|------|-------|
| NCS HC37619 | 12.53 | 83.83 | 0.008 | 0.22 | 50 g |

CRM SILICON CARBIDE

analysis listed in mass %

| Number | SiC | Free C | Si | SiO ₂ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | MgO | Units |
|-------------|-------|--------|------|------------------|--------------------------------|-------|--------------------------------|-------|-------|
| NCS DC93028 | 97.87 | 0.48 | 0.18 | 0.55 | 0.10 | 0.055 | 0.39 | 0.008 | 50 g |
| NCS DC93026 | 84.09 | 1.71 | 1.45 | 6.15 | 1.41 | 0.17 | 0.86 | 0.082 | 50 g |
| NCS DC93027 | 90.86 | 3.48 | 0.24 | 2.00 | 0.77 | 0.47 | 1.12 | 0.039 | 50 g |

CRM TUNGSTEN CARBIDE

analysis listed in mass %

SRM 276b: 75 g units

all others: 100 g units

| Number | Grade | C | Free C | Co | Fe | Mo | Nb | Ni | Ta | Ti |
|------------|-----------------|-------|--------|-------|---------|---------|---------|---------|---------|---------|
| ECRM 783-1 | W94-C6 | 6.188 | (0.04) | . | 0.0022 | . | . | . | . | . |
| SRM 889 | W75-Co9-Ta5-Ti4 | (6.0) | . | 9.50 | (<0.05) | (<0.05) | (<0.05) | (<0.05) | 4.60 | 4.03 |
| SRM 887 | W83-Co10 | (5.5) | . | 10.35 | (<0.05) | (<0.05) | (<0.05) | (<0.01) | (<0.01) | (<0.05) |
| SRM 888 | W64-Co25-Ta-5 | (4.6) | . | 24.7 | (<0.05) | (<0.05) | (<0.05) | (<0.05) | 4.77 | (0.04) |

CRM URBAN AEROSOLS analysis listed in mass %

| Number | Al | Ba | Ca | Cl | Cu | Fe | K | Mg | Mn | Na | P | Pb | S | Si | Sr | Ti | Zn |
|---------|------|--------|------|---------|--------|------|------|------|--------|-------|---------|--------|--------|--------|--------|-------|-------|
| NIES 28 | 5.04 | 0.0874 | 6.69 | (0.807) | 0.0104 | 2.92 | 1.37 | 1.40 | 0.0686 | 0.796 | (0.145) | 0.0403 | (3.91) | (14.9) | 0.0469 | 0.292 | 0.114 |

analysis listed in mg/kg

| Number | As | Be | Cd | Co | Cr | La | Mo | Ni | Rb | Sb | Sc | Se | Sn | Th | U | V | Y | Units |
|---------|------|--------|------|--------|--------|--------|--------|------|--------|--------|--------|--------|--------|--------|------|------|--------|-------|
| NIES 28 | 90.2 | (5.09) | 5.60 | (22.0) | (65.6) | (32.7) | (28.4) | 63.8 | (64.1) | (20.1) | (10.7) | (14.4) | (21.5) | (11.1) | 4.33 | 73.2 | (21.9) | 1.5 g |

CRM URBAN PARTICULATE MATTER analysis listed in mass % Org = organic Elem = Elemental powder 2 g

| Number | Al | C | C.Org | C.Elem | Ca | Cl | Cu | Fe | K | Mg | Mn | Na | Pb | S | Si | Ti | Zn |
|-----------|------|--------|--------|--------|------|--------|--------|------|-------|-------|--------|--------|-------|------|------|--------|--------|
| SRM 1648a | 3.43 | (12.7) | (10.5) | (2.3) | 5.84 | 0.4543 | 0.0610 | 3.92 | 1.056 | 0.813 | 0.0790 | 0.4240 | 0.655 | 5.51 | 12.8 | 0.4021 | 0.4800 |

analysis listed in mg/kg

| Number | Ag | As | B | Br | Cd | Ce | Co | Cr | Cs | Hf | La | Ni | Rb | Sb | Sc | Se | Sm | Sr | Th | V | W |
|-----------|-----|-------|-----|-----|------|------|-------|-----|-----|-------|----|------|------|------|---------|------|-----|-----|---------|-----|-----|
| SRM 1648a | 6.0 | 115.5 | 161 | 502 | 73.7 | 54.6 | 17.93 | 402 | 3.4 | (5.2) | 39 | 81.1 | 51.0 | 45.4 | (6-120) | 28.4 | 4.3 | 215 | (7-107) | 127 | 4.6 |

CRM VANADIUM NITROGEN ALLOY analysis listed in mass %

| Number | V | N | C | O | Al | Mn | P | S | Si | As | Ca | Cr | Fe | Units |
|-------------|-------|-------|------|-------|-------|--------|--------|--------|-------|--------|-------|--------|------|-------|
| NCS HC28641 | 78.04 | 14.13 | 5.71 | (0.6) | 0.26 | 0.0065 | 0.012 | 0.0013 | 0.26 | 0.0014 | 0.064 | 0.082 | 0.65 | 25 g |
| NCS HC28642 | 77.73 | 16.64 | 3.39 | (0.6) | 0.24 | 0.0050 | 0.010 | 0.0016 | 0.23 | 0.0012 | 0.044 | 0.082 | 0.57 | 25 g |
| NCS HC93630 | 77.73 | 14.57 | 3.96 | . | 0.164 | 0.0082 | 0.0075 | 0.0014 | 0.061 | . | . | . | . | 25 g |
| NCS HC28639 | 77.58 | 9.44 | 9.22 | (0.5) | 0.24 | 0.0091 | 0.147 | 0.0025 | 0.40 | 0.0074 | 0.066 | 0.0032 | 1.95 | 25 g |
| NCS HC28640 | 76.73 | 13.31 | 6.01 | (0.7) | 0.28 | 0.0045 | 0.142 | 0.0019 | 0.40 | 0.012 | 0.10 | 0.019 | 1.76 | 25 g |
| Y 19606 | 76.57 | 15.55 | 3.07 | 2.26 | 0.043 | 0.076 | 0.011 | 0.0099 | 0.24 | . | . | . | . | 100 g |
| Y 19607 | 76.56 | 15.50 | 2.93 | 2.60 | 0.050 | 0.050 | 0.011 | 0.018 | 0.23 | . | . | . | . | 100 g |

CRM ZIRCON CONCENTRATE * SARM 62 lists Total Fe as Fe₂O₃ and Ti as TiO₂ DSU: 20 or 50 g others: 100 g

| Number | ZrO ₂ | ZrO ₂ +HfO ₂ | SiO ₂ | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | K ₂ O | HfO ₂ | MgO | Na ₂ O | P ₂ O ₅ | SnO ₂ | TiO ₂ | ThO ₂ | U ₃ O ₈ | LOI |
|----------------|------------------|------------------------------------|------------------|--------------------------------|--------|--------------------------------|------------------|------------------|--------|-------------------|-------------------------------|------------------|------------------|------------------|-------------------------------|------|
| DSZU 123.47-03 | . | 66.1 | . | 0.75 | . | 0.074 | . | . | . | 0.099 | . | . | 0.22 | . | . | . |
| SARM 62 * | 64.2 | . | 32.8 | 0.88 | (0.11) | 0.07 | . | 1.31 | (0.04) | . | 0.12 | . | 0.13 | 0.0158 | 0.0354 | . |
| BCS 204A | . | 53.8 | 37.6 | 0.74 | 0.15 | 0.18 | 0.017 | . | 0.012 | 0.014 | 0.77 | 1.69 | 2.22 | . | . | 0.50 |

CRM ZIRCONIA - Yttrium Stabilized Zirconium Oxide

| Number | Al | Ca | Fe | Hf | Mg | P | Si | Th | Ti | U | Y | monoclinic ZrO ₂ | Units |
|-----------|--------|--------|--------|-------|---------|-----------|--------|--------|--------|--------|------|-----------------------------|-------------|
| ERM-ED105 | 0.0660 | 0.0242 | 0.0095 | 1.535 | 0.00129 | (<0.0075) | 0.0195 | 0.0112 | 0.0497 | 0.0292 | 6.11 | (1.94) | powder 47 g |

CRM ZIRCONIA SET available in SET/4 only 50 g units

| Number | Al ₂ O ₃ | CaO | Fe ₂ O ₃ | HfO ₂ | K ₂ O | MgO | Na ₂ O | Nb ₂ O ₅ | SiO ₂ | TiO ₂ | LOI |
|-----------|--------------------------------|--------|--------------------------------|------------------|------------------|--------|-------------------|--------------------------------|------------------|------------------|------|
| JCRM R051 | . | 0.0017 | 0.0017 | 1.96 | . | 0.0004 | 0.015 | . | (0.005) | (0.0005) | 0.71 |
| JCRM R052 | . | 0.019 | (0.0004) | 1.81 | 0.0013 | 0.0042 | 0.0021 | . | 0.019 | 0.0012 | 0.25 |
| JCRM R053 | . | 0.021 | 0.030 | 1.67 | (0.0007) | 0.0020 | 0.028 | 0.054 | 0.036 | 0.127 | 0.65 |
| JCRM R054 | 0.136 | 0.535 | 0.132 | 1.60 | (0.0003) | 0.208 | 0.0027 | 0.427 | 0.300 | 0.138 | 0.15 |